

Package: scCustomize (via r-universe)

October 25, 2024

Type Package

Title Custom Visualizations & Functions for Streamlined Analyses of Single Cell Sequencing

Description Collection of functions created and/or curated to aid in the visualization and analysis of single-cell data using 'R'. 'scCustomize' aims to provide 1) Customized visualizations for aid in ease of use and to create more aesthetic and functional visuals. 2) Improve speed/reproducibility of common tasks/pieces of code in scRNA-seq analysis with a single or group of functions. For citation please use: Marsh SE (2021) ``Custom Visualizations & Functions for Streamlined Analyses of Single Cell Sequencing" <[doi:10.5281/zenodo.5706430](https://doi.org/10.5281/zenodo.5706430)> RRID:SCR_024675.

Version 2.1.2

Date 2024-02-27

URL <https://github.com/samuel-marsh/scCustomize>,
<https://samuel-marsh.github.io/scCustomize/>,
<https://doi.org/10.5281/zenodo.5706431>

BugReports <https://github.com/samuel-marsh/scCustomize/issues>

Depends R (>= 4.0.0), Seurat (>= 4.3.0.1)

Imports circlize, cli (>= 3.2.0), cowplot, data.table, dplyr, forcats, ggbeeswarm, ggplot2, ggprism, ggrastr, ggrepel, glue, grDevices, grid, janitor, lifecycle, magrittr, Matrix (>= 1.5.0), methods, paletteer, patchwork, pbapply, purrr, rlang (>= 1.1.3), scales, scattermore (>= 1.2), SeuratObject (>= 5.0.0), stats, stringi, stringr, tibble, tidyr

Suggests BiocFileCache, ComplexHeatmap, dittoSeq, DropletUtils, ggpubr, hdf5r, knitr, Nebulosa, remotes, reticulate, rliker, rmarkdown, scuttle, tidysselect, qs, viridis

License GPL (>= 3)

Encoding UTF-8

LazyData true
Roxygen list(markdown = TRUE)
RoxygenNote 7.3.1
Repository https://samuel-marsh.r-universe.dev
RemoteUrl https://github.com/samuel-marsh/sccustomize
RemoteRef HEAD
RemoteSha fc7a282af3bef6e1cb816fac3b9229f536159238

Contents

Add_Alt_Feature_ID	5
Add_CellBender_Diff	6
Add_Cell_Complexity	7
Add_Cell_QC_Metrics	8
Add_Mito_Ribo	10
Add_Pct_Diff	12
Add_Sample_Meta	13
Add_Top_Gene_Pct_Seurat	15
as.anndata	16
as.LIGER	18
as.Seurat.liger	20
Barcode_Plot	22
Blank_Theme	23
Case_Check	23
CellBender_Diff_Plot	24
CellBender_Feature_Diff	26
Cell_Highlight_Plot	27
Change_Delim_All	28
Change_Delim_Prefix	29
Change_Delim_Suffix	30
CheckMatrix_scCustom	31
Clustered_DotPlot	32
Cluster_Highlight_Plot	35
Cluster_Stats_All_Samples	36
ColorBlind_Pal	37
Convert_Assay	38
Copy_From_GCP	38
Copy_To_GCP	39
Create_10X_H5	40
Create_CellBender_Merged_Seurat	41
Create_Cluster_Annotation_File	42
Dark2_Pal	42
DimPlot_All_Samples	43
DimPlot_LIGER	44
DimPlot_scCustom	46
DiscretePalette_scCustomize	48

DotPlot_scCustom	49
ensembl_mito_id	51
ensembl_ribo_id	51
Extract_Modality	52
Extract_Sample_Meta	52
Extract_Top_Markers	53
FeaturePlot_DualAssay	55
FeaturePlot_scCustom	56
FeatureScatter_scCustom	58
Feature_Present	60
Fetch_Meta	61
Gene_Present	62
Hue_Pal	63
ieg_gene_list	64
Iterate_Barcode_Rank_Plot	64
Iterate_Cluster_Highlight_Plot	66
Iterate_DimPlot_bySample	67
Iterate_FeaturePlot_scCustom	68
Iterate_Meta_Highlight_Plot	70
Iterate_PC>Loading_Plots	72
Iterate_Plot_Density_Custom	73
Iterate_Plot_Density_Joint	74
Iterate_VlnPlot_scCustom	76
JCO_Four	77
LIGER_Features	78
Liger_to_Seurat	79
MAD_Stats	80
Median_Stats	81
Merge_Seurat_List	82
Merge_Sparse_Data_All	83
Merge_Sparse_Multimodal_All	84
Meta_Highlight_Plot	85
Meta_Numeric	86
Meta_Present	87
Meta_Remove_Seurat	88
Move_Legend	89
msigdb_qc_gene_list	89
NavyAndOrange	91
PalettePlot	91
PC_Plotting	92
Percent_Expressing	93
plotFactors_scCustom	94
Plot_Cells_per_Sample	96
Plot_Density_Custom	97
Plot_Density_Joint_Only	98
Plot_Median_Genes	99
Plot_Median_Mito	100
Plot_Median_Other	102

Plot_Median_UMIs	103
Pull_Cluster_Annotation	104
Pull_Directory_List	105
QC_Histogram	106
QC_Plots_Combined_Vln	107
QC_Plots_Complexity	108
QC_Plots_Feature	110
QC_Plots_Genes	112
QC_Plots_Mito	113
QC_Plots_UMIs	115
QC_Plot_GenevsFeature	116
QC_Plot_UMIvsFeature	118
QC_Plot_UMIvsGene	119
Read10X_GEO	121
Read10X_h5_GEO	123
Read10X_h5_Multi_Directory	124
Read10X_Multi_Directory	126
Read_CellBender_h5_Mat	127
Read_CellBender_h5_Multi_Directory	128
Read_CellBender_h5_Multi_File	130
Read_GEO_Delim	131
Read_Metrics_10X	132
Reduction_Loading_Present	133
Rename_Clusters	134
Replace_Suffix	135
scCustomize_Palette	136
Seq_QC_Plot_Alignment_Combined	137
Seq_QC_Plot_Antisense	138
Seq_QC_Plot_Basic_Combined	139
Seq_QC_Plot_Exonic	140
Seq_QC_Plot_Genes	141
Seq_QC_Plot_Genome	142
Seq_QC_Plot_Intergenic	143
Seq_QC_Plot_Intronic	144
Seq_QC_Plot_Number_Cells	145
Seq_QC_Plot_Reads_in_Cells	146
Seq_QC_Plot_Reads_per_Cell	147
Seq_QC_Plot_Saturation	148
Seq_QC_Plot_Total_Genes	149
Seq_QC_Plot_Transcriptome	150
Seq_QC_Plot_UMIs	151
Setup_scRNAseq_Project	152
Single_Color_Palette	153
Split_Layers	153
Split_Vector	154
Stacked_VlnPlot	155
Store_Misc_Info_Seurat	157
Store_Palette_Seurat	158

theme_ggprism_mod	159
Top_Genes_Factor	160
UnRotate_X	161
Updated_HGNC_Symbols	161
VariableFeaturePlot_scCustom	162
Variable_Features_ALL_LIGER	164
viridis_plasma_dark_high	165
VlnPlot_scCustom	166

Index	168
--------------	------------

Add_Alt_Feature_ID	<i>Add Alternative Feature IDs</i>
--------------------	------------------------------------

Description

Add alternative feature ids to the assay level meta.data slot in Assay5 compatible object (Seurat V5.0.0 or greater)

Usage

```
Add_Alt_Feature_ID(
  seurat_object,
  features_tsv_file = NULL,
  hdf5_file = NULL,
  assay = NULL
)
```

Arguments

seurat_object	object name.
features_tsv_file	output file from Cell Ranger used for creation of Seurat object. (Either provide this or hdf5_file)
hdf5_file	output file from Cell Ranger used for creation of Seurat object. (Either provide this or features_tsv_file)
assay	name of assay(s) to add the alternative features to. Can specify "all" to add to all assays.

Value

Seurat Object with new entries in the obj@assays\$ASSAY@meta.data slot.

Examples

```
## Not run:
# Using features.tsv.gz file
# Either file from filtered or raw outputs can be used as they are identical.
obj <- Add_Alt_Feature_ID(seurat_object = obj,
  features_tsv = "sample01/outs/filtered_feature_bc_matrix/features.tsv.gz", assay = "RNA")

#' # Using hdf5 file
# Either filtered_feature_bc or raw_feature_bc can be used as the features slot is identical
# Though it is faster to load filtered_feature_bc file due to droplet filtering
obj <- Add_Alt_Feature_ID(seurat_object = obj,
  hdf5_file = "sample01/outs/outs/filtered_feature_bc_matrix.h5", assay = "RNA")

## End(Not run)
```

Add_CellBender_Diff *Calculate and add differences post-cell bender analysis*

Description

Calculate the difference in features and UMIs per cell when both cell bender and raw assays are present.

Usage

```
Add_CellBender_Diff(seurat_object, raw_assay_name, cell_bender_assay_name)
```

Arguments

```
seurat_object    object name.
raw_assay_name   name of the assay containing the raw data.
cell_bender_assay_name
                  name of the assay containing the Cell Bender'ed data.
```

Value

Seurat object with 2 new columns in the meta.data slot.

Examples

```
## Not run:
object <- Add_CellBender_Diff(seurat_object = obj, raw_assay_name = "RAW",
  cell_bender_assay_name = "RNA")

## End(Not run)
```

Add_Cell_Complexity *Add Cell Complexity*

Description

Add measure of cell complexity/novelty (\log_{10} GenesPerUMI) for data QC.

Usage

```
Add_Cell_Complexity(object, ...)  
  
## S3 method for class 'liger'  
Add_Cell_Complexity(  
  object,  
  meta_col_name = "log10GenesPerUMI",  
  overwrite = FALSE,  
  ...  
)  
  
## S3 method for class 'Seurat'  
Add_Cell_Complexity(  
  object,  
  meta_col_name = "log10GenesPerUMI",  
  assay = "RNA",  
  overwrite = FALSE,  
  ...  
)
```

Arguments

object	Seurat or LIGER object
...	Arguments passed to other methods
meta_col_name	name to use for new meta data column. Default is "log10GenesPerUMI".
overwrite	Logical. Whether to overwrite existing an meta.data column. Default is FALSE meaning that function will abort if column with name provided to meta_col_name is present in meta.data slot.
assay	assay to use in calculation. Default is "RNA". <i>Note</i> This should only be changed if storing corrected and uncorrected assays in same object (e.g. outputs of both Cell Ranger and Cell Bender).

Value

An object of the same class as object with columns added to object meta data.

Examples

```
## Not run:
# Liger
liger_object <- Add_Cell_Complexity(object = liger_object)

## End(Not run)

# Seurat
library(Seurat)
pbmc_small <- Add_Cell_Complexity(object = pbmc_small)
```

Add_Cell_QC_Metrics *Add Multiple Cell Quality Control Values with Single Function*

Description

Add Mito/Ribo %, Cell Complexity (log10GenesPerUMI), Top Gene Percent with single function call

Usage

```
Add_Cell_QC_Metrics(
  seurat_object,
  add_mito_ribo = TRUE,
  add_complexity = TRUE,
  add_top_pct = TRUE,
  add_MSigDB = TRUE,
  add_IEG = TRUE,
  add_cell_cycle = TRUE,
  species,
  mito_name = "percent_mito",
  ribo_name = "percent_ribo",
  mito_ribo_name = "percent_mito_ribo",
  complexity_name = "log10GenesPerUMI",
  top_pct_name = NULL,
  oxphos_name = "percent_oxphos",
  apop_name = "percent_apop",
  dna_repair_name = "percent_dna_repair",
  ieg_name = "percent_ieg",
  mito_pattern = NULL,
  ribo_pattern = NULL,
  mito_features = NULL,
  ribo_features = NULL,
  ensembl_ids = FALSE,
  num_top_genes = 50,
  assay = NULL,
```



```

    overwrite = FALSE
  )

```

Arguments

<code>seurat_object</code>	object name.
<code>add_mito_ribo</code>	logical, whether to add percentage of counts belonging to mitochondrial/ribosomal genes to object (Default is TRUE).
<code>add_complexity</code>	logical, whether to add Cell Complexity to object (Default is TRUE).
<code>add_top_pct</code>	logical, whether to add Top Gene Percentages to object (Default is TRUE).
<code>add_MSigDB</code>	logical, whether to add percentages of counts belonging to genes from of mSigDB hallmark gene lists: "HALLMARK_OXIDATIVE_PHOSPHORYLATION", "HALLMARK_APOPTOSIS", and "HALLMARK_DNA_REPAIR" to object (Default is TRUE).
<code>add_IEG</code>	logical, whether to add percentage of counts belonging to IEG genes to object (Default is TRUE).
<code>add_cell_cycle</code>	logical, whether to add cell cycle scores and phase based on CellCycleScoring . Only applicable if <code>species = "human"</code> . (Default is TRUE).
<code>species</code>	Species of origin for given Seurat Object. If mouse, human, marmoset, zebrafish, rat, drosophila, or rhesus macaque (name or abbreviation) are provided the function will automatically generate <code>mito_pattern</code> and <code>ribo_pattern</code> values.
<code>mito_name</code>	name to use for the new meta.data column containing percent mitochondrial counts. Default is "percent_mito".
<code>ribo_name</code>	name to use for the new meta.data column containing percent ribosomal counts. Default is "percent_ribo".
<code>mito_ribo_name</code>	name to use for the new meta.data column containing percent mitochondrial+ribosomal counts. Default is "percent_mito_ribo".
<code>complexity_name</code>	name to use for new meta data column for <code>Add_Cell_Complexity_Seurat</code> . Default is "log10GenesPerUMI".
<code>top_pct_name</code>	name to use for new meta data column for <code>Add_Top_Gene_Pct_Seurat</code> . Default is "percent_topXX", where XX is equal to the value provided to <code>num_top_genes</code> .
<code>oxphos_name</code>	name to use for new meta data column for percentage of MSigDB oxidative phosphorylation counts. Default is "percent_oxphos".
<code>apop_name</code>	name to use for new meta data column for percentage of MSigDB apoptosis counts. Default is "percent_apop".
<code>dna_repair_name</code>	name to use for new meta data column for percentage of MSigDB DNA repair counts. Default is "percent_dna_repair"..
<code>ieg_name</code>	name to use for new meta data column for percentage of IEG counts. Default is "percent_ieg".
<code>mito_pattern</code>	A regex pattern to match features against for mitochondrial genes (will set automatically if species is mouse or human; marmoset features list saved separately).

ribo_pattern	A regex pattern to match features against for ribosomal genes (will set automatically if species is mouse, human, or marmoset).
mito_features	A list of mitochondrial gene names to be used instead of using regex pattern. Will override regex pattern if both are present (including default saved regex patterns).
ribo_features	A list of ribosomal gene names to be used instead of using regex pattern. Will override regex pattern if both are present (including default saved regex patterns).
ensembl_ids	logical, whether feature names in the object are gene names or ensembl IDs (default is FALSE; set TRUE if feature names are ensembl IDs).
num_top_genes	An integer vector specifying the size(s) of the top set of high-abundance genes. Used to compute the percentage of library size occupied by the most highly expressed genes in each cell.
assay	assay to use in calculation. Default is "RNA". <i>Note</i> This should only be changed if storing corrected and uncorrected assays in same object (e.g. outputs of both Cell Ranger and Cell Bender).
overwrite	Logical. Whether to overwrite existing an meta.data column. Default is FALSE meaning that function will abort if column with name provided to meta_col_name is present in meta.data slot.

Value

A Seurat Object

Examples

```
## Not run:
obj <- Add_Cell_QC_Metrics(seurat_object = obj, species = "Human")

## End(Not run)
```

Add_Mito_Ribo

Add Mito and Ribo percentages

Description

Add Mito, Ribo, & Mito+Ribo percentages to meta.data slot of Seurat Object or cell.data slot of Liger object

Usage

```
Add_Mito_Ribo(object, ...)

## S3 method for class 'liger'
Add_Mito_Ribo(
```

```

    object,
    species,
    mito_name = "percent_mito",
    ribo_name = "percent_ribo",
    mito_ribo_name = "percent_mito_ribo",
    mito_pattern = NULL,
    ribo_pattern = NULL,
    mito_features = NULL,
    ribo_features = NULL,
    ensembl_ids = FALSE,
    overwrite = FALSE,
    list_species_names = FALSE,
    ...
)

## S3 method for class 'Seurat'
Add_Mito_Ribo(
  object,
  species,
  mito_name = "percent_mito",
  ribo_name = "percent_ribo",
  mito_ribo_name = "percent_mito_ribo",
  mito_pattern = NULL,
  ribo_pattern = NULL,
  mito_features = NULL,
  ribo_features = NULL,
  ensembl_ids = FALSE,
  assay = NULL,
  overwrite = FALSE,
  list_species_names = FALSE,
  ...
)

```

Arguments

object	Seurat or LIGER object
...	Arguments passed to other methods
species	Species of origin for given Seurat Object. If mouse, human, marmoset, zebrafish, rat, drosophila, or rhesus macaque (name or abbreviation) are provided the function will automatically generate mito_pattern and ribo_pattern values.
mito_name	name to use for the new meta.data column containing percent mitochondrial counts. Default is "percent_mito".
ribo_name	name to use for the new meta.data column containing percent ribosomal counts. Default is "percent_ribo".
mito_ribo_name	name to use for the new meta.data column containing percent mitochondrial+ribosomal counts. Default is "percent_mito_ribo".
mito_pattern	A regex pattern to match features against for mitochondrial genes (will set automatically if species is mouse or human; marmoset features list saved separately).

ribo_pattern	A regex pattern to match features against for ribosomal genes (will set automatically if species is mouse, human, or marmoset).
mito_features	A list of mitochondrial gene names to be used instead of using regex pattern. Will override regex pattern if both are present (including default saved regex patterns).
ribo_features	A list of ribosomal gene names to be used instead of using regex pattern. Will override regex pattern if both are present (including default saved regex patterns).
ensembl_ids	logical, whether feature names in the object are gene names or ensembl IDs (default is FALSE; set TRUE if feature names are ensembl IDs).
overwrite	Logical. Whether to overwrite existing meta.data columns. Default is FALSE meaning that function will abort if columns with any one of the names provided to mito_name ribo_name or mito_ribo_name is present in meta.data slot.
list_species_names	returns list of all accepted values to use for default species names which contain internal regex/feature lists (human, mouse, marmoset, zebrafish, rat, drosophila, and rhesus macaque). Default is FALSE.
assay	Assay to use (default is the current object default assay).

Value

An object of the same class as object with columns added to object meta data.

Examples

```
## Not run:
# Liger
liger_object <- Add_Mito_Ribo(object = liger_object, species = "human")

## End(Not run)

## Not run:
# Seurat
seurat_object <- Add_Mito_Ribo(object = seurat_object, species = "human")

## End(Not run)
```

Add_Pct_Diff

Add percentage difference to DE results

Description

Adds new column labeled "pct_diff" to the data.frame output of [FindMarkers](#), [FindAllMarkers](#), or other DE test data.frames.

Usage

```
Add_Pct_Diff(
  marker_dataframe,
  pct.1_name = "pct.1",
  pct.2_name = "pct.2",
  overwrite = FALSE
)
```

Arguments

marker_dataframe	data.frame containing the results of FindMarkers , FindAllMarkers , or other DE test data.frame.
pct.1_name	the name of data.frame column corresponding to percent expressed in group 1. Default is Seurat default "pct.1".
pct.2_name	the name of data.frame column corresponding to percent expressed in group 2. Default is Seurat default "pct.2".
overwrite	logical. If the marker_dataframe already contains column named "pct_diff" whether to overwrite or return error message. Default is FALSE.

Value

Returns input marker_dataframe with additional "pct_diff" column.

Examples

```
## Not run:
marker_df <- FindAllMarkers(object = obj_name)
marker_df <- Add_Pct_Diff(marker_dataframe = marker_df)
# or piped with function
marker_df <- FindAllMarkers(object = obj_name) %>%
  Add_Pct_Diff()

## End(Not run)
```

Add_Sample_Meta

Add Sample Level Meta Data

Description

Add meta data from ample level data.frame/tibble to cell level seurat @meta.data slot

Usage

```
Add_Sample_Meta(
  seurat_object,
  meta_data,
  join_by_seurat,
  join_by_meta,
  na_ok = FALSE,
  overwrite = FALSE
)
```

Arguments

<code>seurat_object</code>	object name.
<code>meta_data</code>	data.frame/tibble containing meta data or path to file to read. Must be formatted as either data.frame or tibble.
<code>join_by_seurat</code>	name of the column in <code>seurat_object@meta.data</code> that contains matching variables to <code>join_by_meta</code> in <code>meta_data</code> .
<code>join_by_meta</code>	name of the column in <code>meta_data</code> that contains matching variables to <code>join_by_seurat</code> in <code>seurat_object@meta.data</code> .
<code>na_ok</code>	logical, is it ok to add NA values to <code>seurat_object@meta.data</code> . Default is FALSE. Be very careful if setting TRUE because if there is error in join operation it may result in all <code>@meta.data</code> values being replaced with NA.
<code>overwrite</code>	logical, if there are shared columns between <code>seurat_object@meta.data</code> and <code>meta_data</code> should the current <code>seurat_object@meta.data</code> columns be overwritten. Default is FALSE. This parameter excludes values provided to <code>join_by_seurat</code> and <code>join_by_meta</code> .

Value

Seurat object with new `@meta.data` columns

Examples

```
## Not run:
# meta_data present in environment
sample_level_meta <- data.frame(...)
obj <- Add_Sample_Meta(seurat_object = obj, meta_data = sample_level_meta,
  join_by_seurat = "orig.ident", join_by_meta = "sample_ID")

# from meta data file
obj <- Add_Sample_Meta(seurat_object = obj, meta_data = "meta_data/sample_level_meta.csv",
  join_by_seurat = "orig.ident", join_by_meta = "sample_ID")

## End(Not run)
```

`Add_Top_Gene_Pct_Seurat`*Add Percent of High Abundance Genes*

Description

Add the percentage of counts occupied by the top XX most highly expressed genes in each cell.

Usage

```
Add_Top_Gene_Pct_Seurat(  
  seurat_object,  
  num_top_genes = 50,  
  meta_col_name = NULL,  
  assay = "RNA",  
  overwrite = FALSE,  
  verbose = TRUE  
)
```

Arguments

<code>seurat_object</code>	object name.
<code>num_top_genes</code>	An integer vector specifying the size(s) of the top set of high-abundance genes. Used to compute the percentage of library size occupied by the most highly expressed genes in each cell.
<code>meta_col_name</code>	name to use for new meta data column. Default is "percent_topXX", where XX is equal to the value provided to <code>num_top_genes</code> .
<code>assay</code>	assay to use in calculation. Default is "RNA". <i>Note</i> This should only be changed if storing corrected and uncorrected assays in same object (e.g. outputs of both Cell Ranger and Cell Bender).
<code>overwrite</code>	Logical. Whether to overwrite existing an meta.data column. Default is FALSE meaning that function will abort if column with name provided to <code>meta_col_name</code> is present in meta.data slot.
<code>verbose</code>	logical, whether to print messages with status updates, default is TRUE.

Value

A Seurat Object

References

This function uses scuttle package (license: GPL-3) to calculate the percent of expression coming from top XX genes in each cell. Parameter description for `num_top_genes` also from scuttle. If using this function in analysis, in addition to citing scCustomize, please cite scuttle: McCarthy DJ, Campbell KR, Lun ATL, Willis QF (2017). "Scater: pre-processing, quality control, normalisation and visualisation of single-cell RNA-seq data in R." *Bioinformatics*, 33, 1179-1186. doi:10.1093/bioinformatics/btw777.

See Also

<https://bioconductor.org/packages/release/bioc/html/scuttle.html>

Examples

```
## Not run:
library(Seurat)
pbmc_small <- Add_Top_Gene_Pct_Seurat(seurat_object = pbmc_small, num_top_genes = 50)

## End(Not run)
```

as.anndata

Convert objects to anndata objects

Description

Convert objects (Seurat & LIGER) to anndata objects

Usage

```
as.anndata(x, ...)
```

```
## S3 method for class 'Seurat'
as.anndata(
  x,
  file_path,
  file_name,
  assay = "RNA",
  main_layer = "data",
  other_layers = "counts",
  transfer_dimreduc = TRUE,
  verbose = TRUE,
  ...
)
```

```
## S3 method for class 'liger'
as.anndata(
  x,
  file_path,
  file_name,
  transfer_norm.data = FALSE,
  reduction_label = NULL,
  add_barcode_names = FALSE,
  barcode_prefix = TRUE,
  barcode_cell_id_delimiter = "_",
  verbose = TRUE,
```



```
    ...
  )
```

Arguments

x	Seurat or LIGER object
...	Arguments passed to other methods
file_path	directory file path and/or file name prefix. Defaults to current wd.
file_name	file name.
assay	Assay containing data to use, (default is "RNA").
main_layer	the layer of data to become default layer in anndata object (default is "data").
other_layers	other data layers to transfer to anndata object (default is "counts").
transfer_dimreduc	logical, whether to transfer dimensionality reduction coordinates from Seurat to anndata object (default is TRUE).
verbose	logical, whether to print status messages during object conversion (default is TRUE).
transfer_norm.data	logical, whether to transfer the norm.data in addition to raw.data, default is FALSE.
reduction_label	What to label the visualization dimensionality reduction. LIGER does not store name of technique and therefore needs to be set manually.
add_barcode_names	logical, whether to add dataset names to the cell barcodes when merging object data, default is FALSE.
barcode_prefix	logical, if add_barcode_names = TRUE should the names be added as prefix to current cell barcodes/names or a suffix (default is TRUE; prefix).
barcode_cell_id_delimiter	The delimiter to use when adding dataset id to barcode prefix/suffix. Default is "_".

Value

an anndata object generated from x, saved at path provided.

References

Seurat version modified and enhanced version of sceasy: : seurat2anndata (sceasy package: <https://github.com/cellgeni/sceasy>; License: GPL-3. Function has additional checks and supports Seurat V3 and V5 object structure.

LIGER version inspired by sceasy: : seurat2anndata modified and updated to apply to LIGER objects (sceasy package: <https://github.com/cellgeni/sceasy>; License: GPL-3.

Examples

```
## Not run:
as.anndata(x = seurat_object, file_path = "/folder_name", file_name = "anndata_converted.h5ad")

## End(Not run)

## Not run:
as.anndata(x = liger_object, file_path = "/folder_name", file_name = "anndata_converted.h5ad")

## End(Not run)
```

as.LIGER

Convert objects to LIGER objects

Description

Convert objects (Seurat & lists of Seurat Objects) to anndata objects

Usage

```
as.LIGER(x, ...)

## S3 method for class 'Seurat'
as.LIGER(
  x,
  group.by = "orig.ident",
  layers_name = NULL,
  assay = "RNA",
  remove_missing = FALSE,
  renormalize = TRUE,
  use_seurat_var_genes = FALSE,
  use_seurat_dimreduc = FALSE,
  reduction = NULL,
  keep_meta = TRUE,
  verbose = TRUE,
  ...
)

## S3 method for class 'list'
as.LIGER(
  x,
  group.by = "orig.ident",
  dataset_names = NULL,
  assay = "RNA",
  remove_missing = FALSE,
  renormalize = TRUE,
```

```

    use_seurat_var_genes = FALSE,
    var_genes_method = "intersect",
    keep_meta = TRUE,
    verbose = TRUE,
    ...
)

```

Arguments

x	An object to convert to class <code>liger</code>
...	Arguments passed to other methods
group.by	Variable in meta data which contains variable to split data by, (default is "orig.ident").
layers_name	name of meta.data column used to split layers if setting group.by = "layers".
assay	Assay containing raw data to use, (default is "RNA").
remove_missing	logical, whether to remove missing genes with no counts when converting to LIGER object (default is FALSE).
renormalize	logical, whether to perform normalization after LIGER object creation (default is TRUE).
use_seurat_var_genes	logical, whether to transfer variable features from Seurat object to new LIGER object (default is FALSE).
use_seurat_dimreduc	logical, whether to transfer dimensionality reduction coordinates from Seurat to new LIGER object (default is FALSE).
reduction	Name of Seurat reduction to transfer if use_seurat_dimreduc = TRUE.
keep_meta	logical, whether to transfer columns in Seurat meta.data slot to LIGER cell.data slot (default is TRUE).
verbose	logical, whether to print status messages during object conversion (default is TRUE).
dataset_names	optional, vector of names to use for naming datasets.
var_genes_method	how variable genes should be selected from Seurat objects if use_seurat_var_genes = TRUE. Can be either "intersect" or "union", (default is "intersect").

Value

a `liger` object generated from x

References

modified and enhanced version of `rliger::seuratToLiger`.

Examples

```
## Not run:
liger_object <- as.LIGER(x = seurat_object)

## End(Not run)

## Not run:
liger_object <- as.LIGER(x = seurat_object_list)

## End(Not run)
```

as.Seurat.liger *Convert objects to Seurat objects*

Description

Merges raw.data and scale.data of object, and creates Seurat object with these values along with slots containing dimensionality reduction coordinates, iNMF factorization, and cluster assignments. Supports Seurat V3/4 and V4.

Usage

```
## S3 method for class 'liger'
as.Seurat(
  x,
  nms = names(x@H),
  renormalize = TRUE,
  use.liger.genes = TRUE,
  by.dataset = FALSE,
  keep_meta = TRUE,
  reduction_label = "UMAP",
  seurat_assay = "RNA",
  assay_type = NULL,
  add_barcode_names = FALSE,
  barcode_prefix = TRUE,
  barcode_cell_id_delimiter = "_",
  ...
)
```

Arguments

x	liger object.
nms	By default, labels cell names with dataset of origin (this is to account for cells in different datasets which may have same name). Other names can be passed here as vector, must have same length as the number of datasets. (default names(H)).
renormalize	Whether to log-normalize raw data using Seurat defaults (default TRUE).

use.liger.genes	Whether to carry over variable genes (default TRUE).
by.dataset	Include dataset of origin in cluster identity in Seurat object (default FALSE).
keep_meta	logical. Whether to transfer additional metadata (nGene/nUMI/dataset already transferred) to new Seurat Object. Default is TRUE.
reduction_label	Name of dimensionality reduction technique used. Enables accurate transfer or name to Seurat object instead of defaulting to "tSNE".
seurat_assay	Name to set for assay in Seurat Object. Default is "RNA".
assay_type	what type of Seurat assay to create in new object (Assay vs Assay5). Default is NULL which will default to the current user settings. See Convert_Assay parameter <code>convert_to</code> for acceptable values.
add_barcode_names	logical, whether to add dataset names to the cell barcodes when creating Seurat object, default is FALSE.
barcode_prefix	logical, if <code>add_barcode_names = TRUE</code> should the names be added as prefix to current cell barcodes/names or a suffix (default is TRUE; prefix).
barcode_cell_id_delimiter	The delimiter to use when adding dataset id to barcode prefix/suffix. Default is "_".
...	unused.

Details

Stores original dataset identity by default in new object metadata if dataset names are passed in `nms`. iNMF factorization is stored in `dim.reduction` object with key "iNMF".

Value

Seurat object with `raw.data`, `scale.data`, `reduction_label`, `iNMF`, and `ident` slots set.
Seurat object.

References

Original function is part of LIGER package <https://github.com/welch-lab/liger> (Licence: GPL-3). Function was modified for use in `scCustomize` with additional parameters/functionality.

Examples

```
## Not run:
seurat_object <- as.Seurat(x = liger_object)

## End(Not run)
```

`Barcode_Plot`*Create Barcode Rank Plot*

Description

Plot UMI vs. Barcode Rank with inflection and knee. Requires input from DropletUtils package.

Usage

```
Barcode_Plot(  
  br_out,  
  pt.size = 6,  
  plot_title = "Barcode Ranks",  
  raster_dpi = c(1024, 1024),  
  plateau = NULL  
)
```

Arguments

<code>br_out</code>	DFrame output from barcodeRanks .
<code>pt.size</code>	point size for plotting, default is 6.
<code>plot_title</code>	Title for plot, default is "Barcode Ranks".
<code>raster_dpi</code>	Pixel resolution for rasterized plots, passed to <code>geom_scattermore()</code> . Default is <code>c(1024, 1024)</code> .
<code>plateau</code>	numerical value at which to add vertical line designating estimated empty droplet plateau (default is <code>NULL</code>).

Value

A ggplot object

Examples

```
## Not run:  
mat <- Read10X_h5(filename = "raw_feature_bc_matrix.h5")  
  
br_results <- DropletUtils::barcodeRanks(mat)  
  
Barcode_Plot(br_out = br_results)  
  
## End(Not run)
```

 Blank_Theme

Blank Theme

Description

Shortcut for thematic modification to remove all axis labels and grid lines

Usage

```
Blank_Theme(...)
```

Arguments

... extra arguments passed to `ggplot2::theme()`.

Value

Returns a list-like object of class *theme*.

Examples

```
# Generate a plot and customize theme
library(ggplot2)
df <- data.frame(x = rnorm(n = 100, mean = 20, sd = 2), y = rbinom(n = 100, size = 100, prob = 0.2))
p <- ggplot(data = df, mapping = aes(x = x, y = y)) + geom_point(mapping = aes(color = 'red'))
p + Blank_Theme()
```

 Case_Check

Check for alternate case features Checks Seurat object for the presence of features with the same spelling but alternate case.

Description

Check for alternate case features Checks Seurat object for the presence of features with the same spelling but alternate case.

Usage

```
Case_Check(
  seurat_object,
  gene_list,
  case_check_msg = TRUE,
  return_features = TRUE,
  assay = NULL
)
```

Arguments

`seurat_object` Seurat object name.
`gene_list` vector of genes to check.
`case_check_msg` logical. Whether to print message to console if alternate case features are found in addition to inclusion in returned list. Default is TRUE.
`return_features` logical. Whether to return vector of alternate case features. Default is TRUE.
`assay` Name of assay to pull feature names from. If NULL will use the result of `DefaultAssay(seurat_object)`.

Value

If features found returns vector of found alternate case features and prints message depending on parameters specified.

Examples

```

## Not run:
alt_features <- Case_Check(seurat_object = obj_name, gene_list = DEG_list)

## End(Not run)

```

CellBender_Diff_Plot *Plot Number of Cells/Nuclei per Sample*

Description

Plot of total cell or nuclei number per sample grouped by another meta data variable.

Usage

```

CellBender_Diff_Plot(
  feature_diff_df,
  pct_diff_threshold = 25,
  num_features = NULL,
  label = TRUE,
  num_labels = 20,
  min_count_label = 1,
  repel = TRUE,
  custom_labels = NULL,
  plot_line = TRUE,
  plot_title = "Raw Counts vs. Cell Bender Counts",
  x_axis_label = "Raw Data Counts",
  y_axis_label = "Cell Bender Counts",
  xnudge = 0,

```



```

    ynudge = 0,
    max.overlaps = 100,
    label_color = "dodgerblue",
    fontface = "bold",
    label_size = 3.88,
    bg.color = "white",
    bg.r = 0.15,
    ...
)

```

Arguments

<code>feature_diff_df</code>	name of data.frame created using CellBender_Feature_Diff .
<code>pct_diff_threshold</code>	threshold to use for feature plotting. Resulting plot will only contain features which exhibit percent change \geq value. Default is 25.
<code>num_features</code>	Number of features to plot. Will ignore <code>pct_diff_threshold</code> and return plot with specified number of features. Default is NULL.
<code>label</code>	logical, whether or not to label the features that have largest percent difference between raw and CellBender counts (Default is TRUE).
<code>num_labels</code>	Number of features to label if <code>label = TRUE</code> , (default is 20).
<code>min_count_label</code>	Minimum number of raw counts per feature necessary to be included in plot labels (default is 1)
<code>repel</code>	logical, whether to use <code>geom_text_repel</code> to create a nicely-repelled labels; this is slow when a lot of points are being plotted. If using <code>repel</code> , set <code>xnudge</code> and <code>ynudge</code> to 0, (Default is TRUE).
<code>custom_labels</code>	A custom set of features to label instead of the features most different between raw and CellBender counts.
<code>plot_line</code>	logical, whether to plot diagonal line with slope = 1 (Default is TRUE).
<code>plot_title</code>	Plot title.
<code>x_axis_label</code>	Label for x axis.
<code>y_axis_label</code>	Label for y axis.
<code>xnudge</code>	Amount to nudge X and Y coordinates of labels by.
<code>ynudge</code>	Amount to nudge X and Y coordinates of labels by.
<code>max.overlaps</code>	passed to geom_text_repel , exclude text labels that overlap too many things. Defaults to 100.
<code>label_color</code>	Color to use for text labels.
<code>fontface</code>	font face to use for text labels ("plain", "bold", "italic", "bold.italic") (Default is "bold").
<code>label_size</code>	text size for feature labels (passed to geom_text_repel).
<code>bg.color</code>	color to use for shadow/outline of text labels (passed to geom_text_repel) (Default is white).

`bg.r` radius to use for shadow/outline of text labels (passed to `geom_text_repel`) (Default is 0.15).
`...` Extra parameters passed to `geom_text_repel` through `LabelPoints`.

Value

A ggplot object

Examples

```
## Not run:
# get cell bender differences data.frame
cb_stats <- CellBender_Feature_Diff(seurat_object = obj, raw_assay = "RAW",
cell_bender_assay = "RNA")

# plot
CellBender_Diff_Plot(feature_diff_df = cb_stats, pct_diff_threshold = 25)

## End(Not run)
```

CellBender_Feature_Diff

CellBender Feature Differences

Description

Get quick values for raw counts, CellBender counts, count differences, and percent count differences per feature.

Usage

```
CellBender_Feature_Diff(
  seurat_object = NULL,
  raw_assay = NULL,
  cell_bender_assay = NULL,
  raw_mat = NULL,
  cell_bender_mat = NULL
)
```

Arguments

`seurat_object` Seurat object name.
`raw_assay` Name of the assay containing the raw count data.
`cell_bender_assay` Name of the assay containing the CellBender count data.
`raw_mat` Name of raw count matrix in environment if not using Seurat object.
`cell_bender_mat` Name of CellBender count matrix in environment if not using Seurat object.

Value

A data.frame containing summed raw counts, CellBender counts, count difference, and percent difference in counts.

Examples

```
## Not run:
cb_stats <- CellBender_Feature_Diff(seurat_object - obj, raw_assay = "RAW",
cell_bender_assay = "RNA")

## End(Not run)
```

Cell_Highlight_Plot *Meta Highlight Plot*

Description

Create Plot with meta data variable of interest highlighted

Usage

```
Cell_Highlight_Plot(
  seurat_object,
  cells_highlight,
  highlight_color = NULL,
  background_color = "lightgray",
  pt.size = NULL,
  aspect_ratio = NULL,
  figure_plot = FALSE,
  raster = NULL,
  raster.dpi = c(512, 512),
  label = FALSE,
  split.by = NULL,
  split_seurat = FALSE,
  ggplot_default_colors = FALSE,
  ...
)
```

Arguments

`seurat_object` Seurat object name.
`cells_highlight`
 Cell names to highlight in named list.
`highlight_color`
 Color to highlight cells.

background_color	non-highlighted cell colors (default is "lightgray").
pt.size	point size for both highlighted cluster and background.
aspect_ratio	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.
figure_plot	logical. Whether to remove the axes and plot with legend on left of plot denoting axes labels. (Default is FALSE). Requires split_seurat = TRUE.
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
raster.dpi	Pixel resolution for rasterized plots, passed to geom_scattermore(). Default is c(512, 512).
label	Whether to label the highlighted meta data variable(s). Default is FALSE.
split.by	Variable in @meta.data to split the plot by.
split_seurat	logical. Whether or not to display split plots like Seurat (shared y axis) or as individual plots in layout. Default is FALSE.
ggplot_default_colors	logical. If highlight_color = NULL, Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
...	Extra parameters passed to <code>DimPlot</code> .

Value

A ggplot object

Examples

```
library(Seurat)

# Creating example non-overlapping vectors of cells
MS4A1 <- WhichCells(object = pbmc_small, expression = MS4A1 > 4)
GZMB <- WhichCells(object = pbmc_small, expression = GZMB > 4)

# Format as named list
cells <- list("MS4A1" = MS4A1,
             "GZMB" = GZMB)

Cell_Highlight_Plot(seurat_object = pbmc_small, cells_highlight = cells)
```

Change_Delim_All

Change all delimiters in cell name

Description

Change all instances of delimiter in cell names from list of data.frames/matrices or single data.frame/matrix

Usage

```
Change_Delim_All(data, current_delim, new_delim)
```

Arguments

data	Either matrix/data.frame or list of matrices/data.frames with the cell barcodes in the column names.
current_delim	a single value of current delimiter.
new_delim	a single value of new delimiter desired.

Value

matrix or data.frame with new column names.

Examples

```
## Not run:  
dge_matrix <- Change_Delim_All(data = dge_matrix, current_delim = ".", new_delim = "-")  
  
## End(Not run)
```

Change_Delim_Prefix *Change barcode prefix delimiter*

Description

Change barcode prefix delimiter from list of data.frames/matrices or single data.frame/matrix

Usage

```
Change_Delim_Prefix(data, current_delim, new_delim)
```

Arguments

data	Either matrix/data.frame or list of matrices/data.frames with the cell barcodes in the column names.
current_delim	a single value of current delimiter.
new_delim	a single value of new delimiter desired.

Value

matrix or data.frame with new column names.

Examples

```
## Not run:  
dge_matrix <- Change_Delim_Prefix(data = dge_matrix, current_delim = ".", new_delim = "-")  
  
## End(Not run)
```

Change_Delim_Suffix *Change barcode suffix delimiter*

Description

Change barcode suffix delimiter from list of data.frames/matrices or single data.frame/matrix

Usage

```
Change_Delim_Suffix(data, current_delim, new_delim)
```

Arguments

data	Either matrix/data.frame or list of matrices/data.frames with the cell barcodes in the column names.
current_delim	a single value of current delimiter.
new_delim	a single value of new delimiter desired.

Value

matrix or data.frame with new column names.

Examples

```
## Not run:  
dge_matrix <- Change_Delim_Suffix(data = dge_matrix, current_delim = ".", new_delim = "-")  
  
## End(Not run)
```

CheckMatrix_scCustom *Check Matrix Validity*

Description

Native implementation of SeuratObjects CheckMatrix but with modified warning messages.

Usage

```
CheckMatrix_scCustom(  
  object,  
  checks = c("infinite", "logical", "integer", "na")  
)
```

Arguments

object	A matrix
checks	Type of checks to perform, choose one or more from: <ul style="list-style-type: none">• “infinite”: Emit a warning if any value is infinite• “logical”: Emit a warning if any value is a logical• “integer”: Emit a warning if any value is <i>not</i> an integer• “na”: Emit a warning if any value is an NA or NaN

Value

Emits warnings for each test and invisibly returns NULL

References

Re-implementing CheckMatrix only for sparse matrices with modified warning messages. Original function from SeuratObject <https://github.com/satijalab/seurat-object/blob/9c0eda946e162d8595696e5280a6eR/utils.R#L625-L650> (License: MIT).

Examples

```
## Not run:  
mat <- Read10X(...)  
CheckMatrix_scCustom(object = mat)  
  
## End(Not run)
```

Clustered_DotPlot *Clustered DotPlot*

Description

Clustered DotPlots using ComplexHeatmap

Usage

```
Clustered_DotPlot(  
  seurat_object,  
  features,  
  split.by = NULL,  
  colors_use_exp = viridis_plasma_dark_high,  
  exp_color_min = -2,  
  exp_color_middle = NULL,  
  exp_color_max = 2,  
  exp_value_type = "scaled",  
  print_exp_quantiles = FALSE,  
  colors_use_idents = NULL,  
  x_lab_rotate = TRUE,  
  plot_padding = NULL,  
  flip = FALSE,  
  k = 1,  
  feature_km_repeats = 1000,  
  ident_km_repeats = 1000,  
  row_label_size = 8,  
  row_label_fontface = "plain",  
  grid_color = NULL,  
  cluster_feature = TRUE,  
  cluster_ident = TRUE,  
  column_label_size = 8,  
  legend_label_size = 10,  
  legend_title_size = 10,  
  raster = FALSE,  
  plot_km_elbow = TRUE,  
  elbow_kmax = NULL,  
  assay = NULL,  
  group.by = NULL,  
  idents = NULL,  
  show_parent_dend_line = TRUE,  
  ggplot_default_colors = FALSE,  
  color_seed = 123,  
  seed = 123  
)
```


Arguments

<code>seurat_object</code>	Seurat object name.
<code>features</code>	Features to plot.
<code>split.by</code>	Variable in <code>@meta.data</code> to split the identities plotted by.
<code>colors_use_exp</code>	Color palette to use for plotting expression scale. Default is <code>viridis::plasma(n = 20, direction = -1)</code> .
<code>exp_color_min</code>	Minimum scaled average expression threshold (everything smaller will be set to this). Default is -2.
<code>exp_color_middle</code>	What scaled expression value to use for the middle of the provided <code>colors_use_exp</code> . By default will be set to value in middle of <code>exp_color_min</code> and <code>exp_color_max</code> .
<code>exp_color_max</code>	Minimum scaled average expression threshold (everything smaller will be set to this). Default is 2.
<code>exp_value_type</code>	Whether to plot average normalized expression or scaled average normalized expression. Only valid when <code>split.by</code> is provided.
<code>print_exp_quantiles</code>	Whether to print the quantiles of expression data in addition to plots. Default is FALSE. NOTE: These values will be altered by choices of <code>exp_color_min</code> and <code>exp_color_max</code> if there are values below or above those cutoffs, respectively.
<code>colors_use_idents</code>	specify color palette to used for identity labels. By default if number of levels plotted is less than or equal to 36 it will use "polychrome" and if greater than 36 will use "varibow" with <code>shuffle = TRUE</code> both from <code>DiscretePalette_scCustomize</code> .
<code>x_lab_rotate</code>	How to rotate column labels. By default set to TRUE which rotates labels 45 degrees. If set FALSE rotation is set to 0 degrees. Users can also supply custom angle for text rotation.
<code>plot_padding</code>	if plot needs extra white space padding so no plot or labels are cutoff. The parameter accepts TRUE or numeric vector of length 4. If TRUE padding will be set to <code>c(2, 10, 0, 0)</code> (bottom, left, top, right). Can also be customized further with numeric vector of length 4 specifying the amount of padding in millimeters. Default is NULL, no padding.
<code>flip</code>	logical, whether to flip the axes of final plot. Default is FALSE; rows = features and columns = idents.
<code>k</code>	Value to use for k-means clustering on features Sets (km) parameter in <code>ComplexHeatmap::Heatmap()</code> . From <code>ComplexHeatmap::Heatmap()</code> : Apply k-means clustering on rows. If the value is larger than 1, the heatmap will be split by rows according to the k-means clustering. For each row slice, hierarchical clustering is still applied with parameters above.
<code>feature_km_repeats</code>	Number of k-means runs to get a consensus k-means clustering for features. Note if <code>feature_km_repeats</code> is set to value greater than one, the final number of groups might be smaller than <code>row_km</code> , but this might mean the original <code>row_km</code> is not a good choice. Default is 1000.

<code>ident_km_repeats</code>	Number of k-means runs to get a consensus k-means clustering. Similar to <code>feature_km_repeats</code> . Default is 1000.
<code>row_label_size</code>	Size of the feature labels. Provided to <code>row_names_gp</code> in Heatmap call.
<code>row_label_fontface</code>	Fontface to use for row labels. Provided to <code>row_names_gp</code> in Heatmap call.
<code>grid_color</code>	color to use for heatmap grid. Default is NULL which "removes" grid by using NA color.
<code>cluster_feature</code>	logical, whether to cluster and reorder feature axis. Default is TRUE.
<code>cluster_ident</code>	logical, whether to cluster and reorder identity axis. Default is TRUE.
<code>column_label_size</code>	Size of the feature labels. Provided to <code>column_names_gp</code> in Heatmap call.
<code>legend_label_size</code>	Size of the legend text labels. Provided to <code>labels_gp</code> in Heatmap legend call.
<code>legend_title_size</code>	Size of the legend title text labels. Provided to <code>title_gp</code> in Heatmap legend call.
<code>raster</code>	Logical, whether to render in raster format (faster plotting, smaller files). Default is FALSE.
<code>plot_km_elbow</code>	Logical, whether or not to return the Sum Squared Error Elbow Plot for k-means clustering. Estimating elbow of this plot is one way to determine "optimal" value for k. Based on: https://stackoverflow.com/a/15376462/15568251 .
<code>elbow_kmax</code>	The maximum value of k to use for <code>plot_km_elbow</code> . Suggest setting larger value so the true shape of plot can be observed. Value must be 1 less than number of features provided. If NULL parameter will be set dependent on length of feature list up to <code>elbow_kmax = 20</code> .
<code>assay</code>	Name of assay to use, defaults to the active assay.
<code>group.by</code>	Group (color) cells in different ways (for example, <code>orig.ident</code>).
<code>idents</code>	Which classes to include in the plot (default is all).
<code>show_parent_dend_line</code>	Logical, Sets parameter of same name in <code>ComplexHeatmap::Heatmap()</code> . From <code>ComplexHeatmap::Heatmap()</code> : When heatmap is split, whether to add a dashed line to mark parent dendrogram and children dendrograms. Default is TRUE.
<code>ggplot_default_colors</code>	logical. If <code>colors_use = NULL</code> , Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
<code>color_seed</code>	random seed for the "varibow" palette shuffle if <code>colors_use = NULL</code> and number of groups plotted is greater than 36. Default = 123.
<code>seed</code>	Sets seed for reproducible plotting (ComplexHeatmap plot).

Value

A `ComplexHeatmap` or if `plot_km_elbow = TRUE` a list containing ggplot2 object and `ComplexHeatmap`.

Author(s)

Ming Tang (Original Code), Sam Marsh (Wrap single function, added/modified functionality)

References

<https://divingintogeneticsandgenomics.rbind.io/post/clustered-dotplot-for-single-cell-rnaseq/>

See Also

<https://twitter.com/tangming2005>

Examples

```
library(Seurat)
Clustered_DotPlot(seurat_object = pbmc_small, features = c("CD3E", "CD8", "GZMB", "MS4A1"))
```

Cluster_Highlight_Plot

Cluster Highlight Plot

Description

Create Plot with cluster of interest highlighted

Usage

```
Cluster_Highlight_Plot(
  seurat_object,
  cluster_name,
  highlight_color = NULL,
  background_color = "lightgray",
  pt.size = NULL,
  aspect_ratio = NULL,
  figure_plot = FALSE,
  raster = NULL,
  raster.dpi = c(512, 512),
  label = FALSE,
  split.by = NULL,
  split_seurat = FALSE,
  ggplot_default_colors = FALSE,
  ...
)
```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>cluster_name</code>	Name(s) (or number(s)) identity of cluster to be highlighted.
<code>highlight_color</code>	Color(s) to highlight cells. The default is NULL and plot will use <code>scCustomize_Palette()</code> .
<code>background_color</code>	non-highlighted cell colors.
<code>pt.size</code>	point size for both highlighted cluster and background.
<code>aspect_ratio</code>	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.
<code>figure_plot</code>	logical. Whether to remove the axes and plot with legend on left of plot denoting axes labels. (Default is FALSE). Requires <code>split_seurat = TRUE</code> .
<code>raster</code>	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
<code>raster.dpi</code>	Pixel resolution for rasterized plots, passed to <code>geom_scattermore()</code> . Default is <code>c(512, 512)</code> .
<code>label</code>	Whether to label the highlighted cluster(s). Default is FALSE.
<code>split.by</code>	Feature to split plots by (i.e. "orig.ident").
<code>split_seurat</code>	logical. Whether or not to display split plots like Seurat (shared y axis) or as individual plots in layout. Default is FALSE.
<code>ggplot_default_colors</code>	logical. If <code>colors_use = NULL</code> , Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
<code>...</code>	Extra parameters passed to <code>DimPlot</code> .

Value

A ggplot object

Examples

```
Cluster_Highlight_Plot(seurat_object = pbmc_small, cluster_name = "1", highlight_color = "gold",
background_color = "lightgray", pt.size = 2)
```

Cluster_Stats_All_Samples

Calculate Cluster Stats

Description

Calculates both overall and per sample cell number and percentages per cluster based on `orig.ident`

Usage

```
Cluster_Stats_All_Samples(seurat_object, group_by_var = "orig.ident")
```

Arguments

```
seurat_object  Seurat object name.  
group_by_var   meta data column to classify samples (default = "orig.ident").
```

Value

A data.frame

Examples

```
## Not run:  
stats <- Cluster_Stats_All_Samples(seurat_object = object, group_by_var = "orig.ident")  
  
## End(Not run)
```

ColorBlind_Pal

Color Universal Design Short Palette

Description

Shortcut to a modified 8 color palette based on Color Universal Design (CUD) colorblindness friendly palette.

Usage

```
ColorBlind_Pal()
```

Value

modified/reordered color palette (8 colors) based on ditto-seq

References

palette is slightly modified version of the Color Universal Design (CUD) colorblindness friendly palette <https://jfly.uni-koeln.de/color/>.

Examples

```
cols <- ColorBlind_Pal()  
PalettePlot(pal = cols)
```

Convert_Assay *Convert between Seurat Assay types*

Description

Will convert assays within a Seurat object between "Assay" and "Assay5" types.

Usage

```
Convert_Assay(seurat_object, assay = NULL, convert_to)
```

Arguments

seurat_object Seurat object name.

assay name(s) of assays to convert. Default is NULL and will check with users which assays they want to convert.

convert_to value of what assay type to convert current assays to. #'

- Accepted values for V3/4 are: "Assay", "assay", "V3", or "v3".
- Accepted values for V5 are: "Assay5", "assay5", "V5", or "v5".

Examples

```
## Not run:
# Convert to V3/4 assay
obj <- Convert_Assay(seurat_object = obj, convert_to = "V3")

# Convert to 5 assay
obj <- Convert_Assay(seurat_object = obj, convert_to = "V5")

## End(Not run)
```

Copy_From_GCP *Copy folder from GCP bucket from R Console*

Description

Run command from R console without moving to terminal to copy folder from GCP bucket to local storage

Usage

```
Copy_From_GCP(folder_file_path, gcp_bucket_path)
```

Arguments

folder_file_path
 folder to be copied to GCP bucket.

gcp_bucket_path
 GCP bucket path to copy to files.

Value

No return value. Performs system copy from GCP bucket.

Examples

```
## Not run:
Copy_From_GCP(folder_file_path = "plots/", gcp_bucket_path = "gs://bucket_name_and_folder_path")

## End(Not run)
```

Copy_To_GCP

Copy folder to GCP bucket from R Console

Description

Run command from R console without moving to terminal to copy folder to GCP bucket

Usage

```
Copy_To_GCP(folder_file_path, gcp_bucket_path)
```

Arguments

folder_file_path
 folder to be copied to GCP bucket.

gcp_bucket_path
 GCP bucket path to copy to files.

Value

No return value. Performs system copy to GCP bucket.

Examples

```
## Not run:
Copy_To_GCP(folder_file_path = "plots/", gcp_bucket_path = "gs://bucket_name_and_folder_path")

## End(Not run)
```

Create_10X_H5	<i>Create H5 from 10X Outputs</i>
---------------	-----------------------------------

Description

Creates HDF5 formatted output analogous to the outputs created by Cell Ranger and can be read into Seurat, LIGER, or SCE class object. Requires DropletUtils package from Bioconductor.

Usage

```
Create_10X_H5(  
  raw_data_file_path,  
  source_type = "10X",  
  save_file_path,  
  save_name  
)
```

Arguments

`raw_data_file_path` file path to raw data file(s).

`source_type` type of source data (Default is "10X"). Alternatively can provide "Matrix" or "data.frame".

`save_file_path` file path to directory to save file.

`save_name` name prefix for output H5 file.

Value

A HDF5 format file that will be recognized as 10X Cell Ranger formatted file by Seurat or LIGER.

Examples

```
## Not run:  
Create_10X_H5(raw_data_file_path = "file_path", save_file_path = "file_path2", save_name = "NAME")  
  
## End(Not run)
```

`Create_CellBender_Merged_Seurat`*Create Seurat Object with Cell Bender and Raw data*

Description

Enables easy creation of Seurat object which contains both cell bender data and raw count data as separate assays within the object.

Usage

```
Create_CellBender_Merged_Seurat(  
  raw_cell_bender_matrix = NULL,  
  raw_counts_matrix = NULL,  
  raw_assay_name = "RAW",  
  min_cells = 5,  
  min_features = 200,  
  ...  
)
```

Arguments

<code>raw_cell_bender_matrix</code>	matrix file containing the cell bender correct counts.
<code>raw_counts_matrix</code>	matrix file contain the uncorrected Cell Ranger (or other) counts.
<code>raw_assay_name</code>	a key value to use specifying the name of assay. Default is "RAW".
<code>min_cells</code>	value to supply to min.cells parameter of CreateSeuratObject . Default is 5.
<code>min_features</code>	value to supply to min.features parameter of CreateSeuratObject . Default is 200.
<code>...</code>	Extra parameters passed to CreateSeuratObject .

Value

A Seurat Object contain both the Cell Bender corrected counts ("RNA" assay) and uncorrected counts ("RAW" assay; or other name specified to `raw_assay_name`).

Examples

```
## Not run:  
seurat_obj <- Create_CellBender_Merged_Seurat(raw_cell_bender_matrix = cb_matrix,  
  raw_counts_matrix = cr_matrix)  
  
## End(Not run)
```

Create_Cluster_Annotation_File
Create cluster annotation csv file

Description

create annotation file

Usage

```
Create_Cluster_Annotation_File(  
  file_path = NULL,  
  file_name = "cluster_annotation"  
)
```

Arguments

file_path	path to directory to save file. Default is current working directory.
file_name	name to use for annotation file. Function automatically adds file type ".csv" suffix. Default is "cluster_annotation".

Value

No value returned. Creates .csv file.

Examples

```
## Not run:  
Create_Cluster_Annotation_File(file_path = "cluster_annotation_folder_name")  
  
## End(Not run)
```

Dark2_Pal *Dark2 Palette*

Description

Shortcut to Dark2 color palette from RColorBrewer (8 Colors)

Usage

```
Dark2_Pal()
```

Value

"Dark2" color palette (8 colors)

References

Dark2 palette from RColorBrewer being called through paletteer. See RColorBrewer for more info on palettes <https://CRAN.R-project.org/package=RColorBrewer>

Examples

```
cols <- Dark2_Pal()
PalettePlot(pal= cols)
```

DimPlot_All_Samples *DimPlot by Meta Data Column*

Description

Creates DimPlot layout containing all samples within Seurat Object from orig.ident column

Usage

```
DimPlot_All_Samples(  
  seurat_object,  
  meta_data_column = "orig.ident",  
  colors_use = "black",  
  pt.size = NULL,  
  aspect_ratio = NULL,  
  title_size = 15,  
  num_columns = NULL,  
  reduction = NULL,  
  dims = c(1, 2),  
  raster = NULL,  
  raster.dpi = c(512, 512),  
  ...  
)
```

Arguments

seurat_object	Seurat object name.
meta_data_column	Meta data column to split plots by.
colors_use	single color to use for all plots or a vector of colors equal to the number of plots.
pt.size	Adjust point size for plotting.
aspect_ratio	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.
title_size	size for plot title labels.
num_columns	number of columns in final layout plot.

reduction	Dimensionality Reduction to use (if NULL then defaults to Object default).
dims	Which dimensions to plot. Defaults to c(1,2) if not specified.
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
raster.dpi	Pixel resolution for rasterized plots, passed to geom_scattermore(). Default is c(512, 512).
...	Extra parameters passed to DimPlot .

Value

A ggplot object

Examples

```
library(Seurat)

pbmc_small$sample_id <- sample(c("sample1", "sample2"), size = ncol(pbmc_small), replace = TRUE)

DimPlot_All_Samples(seurat_object = pbmc_small, meta_data_column = "sample_id", color = "black",
num_columns = 2)
```

DimPlot_LIGER

DimPlot LIGER Version

Description

Standard and modified version of LIGER's plotByDatasetAndCluster

Usage

```
DimPlot_LIGER(
  liger_object,
  group_by = NULL,
  split_by = NULL,
  colors_use_cluster = NULL,
  colors_use_meta = NULL,
  pt_size = NULL,
  shuffle = TRUE,
  shuffle_seed = 1,
  reduction_label = "UMAP",
  aspect_ratio = NULL,
  label = TRUE,
  label_size = NA,
  label_repel = FALSE,
  label_box = FALSE,
  label_color = "black",
```

```

    combination = FALSE,
    raster = NULL,
    raster.dpi = c(512, 512),
    num_columns = NULL,
    ggplot_default_colors = FALSE,
    color_seed = 123
  )

```

Arguments

<code>liger_object</code>	liger liger_object. Need to perform clustering before calling this function
<code>group_by</code>	Variable to be plotted. If NULL will plot clusters from <code>liger@clusters</code> slot. If <code>combination = TRUE</code> will plot both clusters and meta data variable.
<code>split_by</code>	Variable to split plots by.
<code>colors_use_cluster</code>	colors to use for plotting by clusters. By default if number of levels plotted is less than or equal to 36 will use "polychrome" and if greater than 36 will use "varibow" with <code>shuffle = TRUE</code> both from DiscretePalette_scCustomize .
<code>colors_use_meta</code>	colors to use for plotting by meta data (cell.data) variable. By default if number of levels plotted is less than or equal to 36 it will use "polychrome" and if greater than 36 will use "varibow" with <code>shuffle = TRUE</code> both from DiscretePalette_scCustomize .
<code>pt_size</code>	Adjust point size for plotting.
<code>shuffle</code>	logical. Whether to randomly shuffle the order of points. This can be useful for crowded plots if points of interest are being buried. (Default is TRUE).
<code>shuffle_seed</code>	Sets the seed if randomly shuffling the order of points.
<code>reduction_label</code>	What to label the x and y axes of resulting plots. LIGER does not store name of technique and therefore needs to be set manually. Default is "UMAP".
<code>aspect_ratio</code>	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.
<code>label</code>	logical. Whether or not to label the clusters. ONLY applies to plotting by cluster. Default is TRUE.
<code>label_size</code>	size of cluster labels.
<code>label_repel</code>	logical. Whether to repel cluster labels from each other if plotting by cluster (if <code>group_by = NULL</code> or <code>group_by = "cluster"</code>). Default is FALSE.
<code>label_box</code>	logical. Whether to put a box around the label text (uses <code>geom_text</code> vs <code>geom_label</code>). Default is FALSE.
<code>label_color</code>	Color to use for cluster labels. Default is "black".
<code>combination</code>	logical, whether to return patchwork displaying both plots side by side. (Default is FALSE).
<code>raster</code>	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.

raster.dpi	Pixel resolution for rasterized plots, passed to geom_scattermore(). Default is c(512, 512).
num_columns	Number of columns in plot layout. Only valid if split.by != NULL.
ggplot_default_colors	logical. If colors_use = NULL, Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.

Value

A ggplot/patchwork object

Examples

```
## Not run:
DimPlot_LIGER(liger_object = obj_name, reduction_label = "UMAP")

## End(Not run)
```

DimPlot_scCustom

DimPlot with modified default settings

Description

Creates DimPlot with some of the settings modified from their Seurat defaults (colors_use, shuffle, label).

Usage

```
DimPlot_scCustom(
  seurat_object,
  colors_use = NULL,
  pt.size = NULL,
  reduction = NULL,
  group.by = NULL,
  split.by = NULL,
  split_seurat = FALSE,
  figure_plot = FALSE,
  aspect_ratio = NULL,
  shuffle = TRUE,
  seed = 1,
  label = NULL,
  label.size = 4,
  label.color = "black",
  label.box = FALSE,
```

```

    dims = c(1, 2),
    repel = FALSE,
    raster = NULL,
    raster.dpi = c(512, 512),
    num_columns = NULL,
    ggplot_default_colors = FALSE,
    color_seed = 123,
    ...
)

```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>colors_use</code>	color palette to use for plotting. By default if number of levels plotted is less than or equal to 36 it will use "polychrome" and if greater than 36 will use "varibow" with <code>shuffle = TRUE</code> both from <code>DiscretePalette_scCustomize</code> .
<code>pt.size</code>	Adjust point size for plotting.
<code>reduction</code>	Dimensionality Reduction to use (if NULL then defaults to Object default).
<code>group.by</code>	Name of one or more metadata columns to group (color) cells by (for example, <code>orig.ident</code>); default is the current <code>active.ident</code> of the object.
<code>split.by</code>	Feature to split plots by (i.e. "orig.ident").
<code>split_seurat</code>	logical. Whether or not to display split plots like Seurat (shared y axis) or as individual plots in layout. Default is FALSE.
<code>figure_plot</code>	logical. Whether to remove the axes and plot with legend on left of plot denoting axes labels. (Default is FALSE). Requires <code>split_seurat = TRUE</code> .
<code>aspect_ratio</code>	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.
<code>shuffle</code>	logical. Whether to randomly shuffle the order of points. This can be useful for crowded plots if points of interest are being buried. (Default is TRUE).
<code>seed</code>	Sets the seed if randomly shuffling the order of points.
<code>label</code>	Whether to label the clusters. By default if <code>group.by = NULL</code> <code>label = TRUE</code> , and otherwise it is FALSE.
<code>label.size</code>	Sets size of labels.
<code>label.color</code>	Sets the color of the label text.
<code>label.box</code>	Whether to put a box around the label text (<code>geom_text</code> vs <code>geom_label</code>).
<code>dims</code>	Which dimensions to plot. Defaults to <code>c(1,2)</code> if not specified.
<code>repel</code>	Repel labels.
<code>raster</code>	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
<code>raster.dpi</code>	Pixel resolution for rasterized plots, passed to <code>geom_scattermore()</code> . Default is <code>c(512, 512)</code> .
<code>num_columns</code>	Number of columns in plot layout. Only valid if <code>split.by != NULL</code> .

`ggplot_default_colors` logical. If `colors_use = NULL`, Whether or not to return plot using default `ggplot2` "hue" palette instead of default "polychrome" or "varibow" palettes.

`color_seed` random seed for the "varibow" palette shuffle if `colors_use = NULL` and number of groups plotted is greater than 36. Default = 123.

... Extra parameters passed to `DimPlot`.

Value

A `ggplot` object

References

Many of the param names and descriptions are from Seurat to facilitate ease of use as this is simply a wrapper to alter some of the default parameters <https://github.com/satijalab/seurat/blob/master/R/visualization.R> (License: GPL-3). `figure_plot` parameter/code modified from code by Tim Stuart via twitter: <https://twitter.com/timoast/status/1526237116035891200?s=20&t=foJ0F81aPSjr1t7pk1cUPg>.

Examples

```
library(Seurat)
DimPlot_scCustom(seurat_object = pbmc_small)
```

DiscretePalette_scCustomize
Discrete color palettes

Description

Helper function to return a number of discrete color palettes.

Usage

```
DiscretePalette_scCustomize(
  num_colors,
  palette = NULL,
  shuffle_pal = FALSE,
  seed = 123
)
```


Arguments

num_colors	Number of colors to be generated.
palette	Options are "alphabet", "alphabet2", "glasbey", "polychrome", "stepped", "ditto_seq", "varibow".
shuffle_pal	randomly shuffle the outputted palette. Most useful for varibow palette as that is normally an ordered palette.
seed	random seed for the palette shuffle. Default = 123.

Value

A vector of colors

References

This function uses the paletteer package <https://github.com/EmilHvitfeldt/paletteer> to provide simplified access to color palettes from many different R package sources while minimizing scCustomize current and future dependencies.

The following packages & palettes are called by this function (see individual packages for palette references/citations):

1. pals (via paletteer) <https://CRAN.R-project.org/package=pals>
 - alphabet, alphabet2, glasbey, polychrome, and stepped.
2. dittoSeq <https://bioconductor.org/packages/release/bioc/html/dittoSeq.html>
 - dittoColors.
3. colorway <https://github.com/hypercompetent/colorway>
 - varibow

Function name and implementation modified from Seurat (License: GPL-3). <https://github.com/satijalab/seurat>

Examples

```
pal <- DiscretePalette_scCustomize(num_colors = 36, palette = "varibow")
PalettePlot(pal= pal)
```

DotPlot_scCustom

Customized DotPlot

Description

Code for creating customized DotPlot

Usage

```
DotPlot_scCustom(
  seurat_object,
  features,
  group.by = NULL,
  colors_use = viridis_plasma_dark_high,
  remove_axis_titles = TRUE,
  x_lab_rotate = FALSE,
  y_lab_rotate = FALSE,
  facet_label_rotate = FALSE,
  flip_axes = FALSE,
  ...
)
```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>features</code>	Features to plot.
<code>group.by</code>	Name of one or more metadata columns to group (color) cells by (for example, <code>orig.ident</code>); default is the current active.ident of the object.
<code>colors_use</code>	specify color palette to used. Default is <code>viridis_plasma_dark_high</code> .
<code>remove_axis_titles</code>	logical. Whether to remove the x and y axis titles. Default = TRUE.
<code>x_lab_rotate</code>	Rotate x-axis labels 45 degrees (Default is FALSE).
<code>y_lab_rotate</code>	Rotate x-axis labels 45 degrees (Default is FALSE).
<code>facet_label_rotate</code>	Rotate facet labels on grouped DotPlots by 45 degrees (Default is FALSE).
<code>flip_axes</code>	whether or not to flip and X and Y axes (Default is FALSE).
<code>...</code>	Extra parameters passed to DotPlot .

Value

A ggplot object

Examples

```
library(Seurat)
DotPlot_scCustom(seurat_object = pbmc_small, features = c("CD3E", "CD8", "GZMB", "MS4A1"))
```

ensembl_mito_id	<i>Ensembl Mito IDs</i>
-----------------	-------------------------

Description

A list of ensembl ids for mitochondrial genes (Ensembl version 105)

Usage

ensembl_mito_id

Format

A list of six vectors

Mus_musculus_mito_ensembl Ensembl IDs for mouse mitochondrial genes

Homo_sapiens_mito_ensembl Ensembl IDs for human mitochondrial genes

Danio_rerio_mito_ensembl Ensembl IDs for zebrafish mitochondrial genes

Rattus_norvegicus_mito_ensembl Ensembl IDs for rat mitochondrial genes

Drosophila_melanogaster_mito_ensembl Ensembl IDs for fly mitochondrial genes

Macaca_mulatta_mito_ensembl Ensembl IDs for macaque mitochondrial genes

ensembl_ribo_id	<i>Ensembl Ribo IDs</i>
-----------------	-------------------------

Description

A list of ensembl ids for ribosomal genes (Ensembl version 105)

Usage

ensembl_ribo_id

Format

A list of seven vectors

Mus_musculus_ribo_ensembl Ensembl IDs for mouse ribosomal genes

Homo_sapiens_ribo_ensembl Ensembl IDs for human ribosomal genes

Callithrix_jacchus_ribo_ensembl Ensembl IDs for marmoset ribosomal genes

Danio_rerio_ribo_ensembl Ensembl IDs for zebrafish ribosomal genes

Rattus_norvegicus_ribo_ensembl Ensembl IDs for rat ribosomal genes

Drosophila_melanogaster_ribo_ensembl Ensembl IDs for fly ribosomal genes

Macaca_mulatta_ribo_ensembl Ensembl IDs for macaque ribosomal genes

Extract_Modality *Extract multi-modal data into list by modality*

Description

Reorganize multi-modal data after import with `Read10X()` or `scCustomize` read functions. Organizes sub-lists by data modality instead of by sample.

Usage

```
Extract_Modality(matrix_list)
```

Arguments

`matrix_list` list of matrices to split by modality

Value

list of lists, with one sublist per data modality. Sub-list contain 1 matrix entry per sample

Examples

```
## Not run:
multi_mat <- Read10X(...)
new_multi_mat <- Extract_Modality(matrix_list = multi_mat)

## End(Not run)
```

Extract_Sample_Meta *Extract sample level meta.data*

Description

Returns a by identity meta.data data.frame with one row per sample. Useful for downstream quick view of sample breakdown, meta data table creation, and/or use in pseudobulk analysis

Usage

```
Extract_Sample_Meta(
  object,
  sample_name = "orig.ident",
  variables_include = NULL,
  variables_exclude = NULL,
  include_all = FALSE
)
```

Arguments

object	Seurat object
sample_name	meta.data column to use as sample. Output data.frame will contain one row per level or unique value in this variable.
variables_include	@meta.data columns to keep in final data.frame. All other columns will be discarded. Default is NULL.
variables_exclude	columns to discard in final data.frame. Many cell level columns are irrelevant at the sample level (e.g., nFeature_RNA, percent_mito). <ul style="list-style-type: none"> • Default parameter value is NULL but internally will set to discard nFeature_ASSAY(s), nCount_ASSAY(s), percent_mito, percent_ribo, percent_mito_ribo, and log10GenesPerUMI. • If sample level median values are desired for these type of variables the output of this function can be joined with output of Median_Stats. • Set parameter to include_all = TRUE to prevent any columns from being excluded.
include_all	logical, whether or not to include all object meta data columns in output data.frame. Default is FALSE.

Value

Returns a data.frame with one row per sample_name.

Examples

```
library(Seurat)
pbmc_small[["batch"]] <- sample(c("batch1", "batch2"), size = ncol(pbmc_small), replace = TRUE)

sample_meta <- Extract_Sample_Meta(object = pbmc_small, sample_name = "orig.ident")

# Only return specific columns from meta data (orig.ident and batch)
sample_meta2 <- Extract_Sample_Meta(object = pbmc_small, sample_name = "orig.ident",
variables_include = "batch")

# Return all columns from meta data
sample_meta3 <- Extract_Sample_Meta(object = pbmc_small, sample_name = "orig.ident",
include_all = TRUE)
```

Extract_Top_Markers *Extract Top N Marker Genes*

Description

Extract vector gene list (or named gene vector) from data.frame results of [FindAllMarkers](#) or similar analysis.

Usage

```
Extract_Top_Markers(
  marker_dataframe,
  num_genes = 10,
  group_by = "cluster",
  rank_by = "avg_log2FC",
  gene_column = "gene",
  gene_rownames_to_column = FALSE,
  data_frame = FALSE,
  named_vector = TRUE,
  make_unique = FALSE
)
```

Arguments

marker_dataframe	data.frame output from FindAllMarkers or similar analysis.
num_genes	number of genes per group (e.g., cluster) to include in output list.
group_by	column name of marker_dataframe to group data by. Default is "cluster" based on FindAllMarkers .
rank_by	column name of marker_dataframe to rank data by when selecting num_genes per group_by. Default is "avg_log2FC" based on FindAllMarkers .
gene_column	column name of marker_dataframe that contains the gene IDs. Default is "gene" based on FindAllMarkers .
gene_rownames_to_column	logical. Whether gene IDs are stored in rownames and should be moved to column. Default is FALSE.
data_frame	Logical, whether or not to return filtered data.frame of the original markers_dataframe or to return a vector of gene IDs. Default is FALSE.
named_vector	Logical, whether or not to name the vector of gene names that is returned by the function. If TRUE will name the vector using the column provided to group_by. Default is TRUE.
make_unique	Logical, whether an unnamed vector should return only unique values. Default is FALSE. Not applicable when data_frame = TRUE or named_vector = TRUE.

Value

filtered data.frame, vector, or named vector containing gene IDs.

Examples

```
## Not run:
top10_genes <- Extract_Top_Markers(marker_dataframe = markers_results, num_genes = 10,
group_by = "cluster", rank_by = "avg_log2FC")

## End(Not run)
```

 FeaturePlot_DualAssay *Customize FeaturePlot of two assays*

Description

Create Custom FeaturePlots and preserve scale (no binning) from same features in two assays simultaneously. Intended for plotting same modality present in two assays.

Usage

```
FeaturePlot_DualAssay(
  seurat_object,
  features,
  assay1 = "RAW",
  assay2 = "RNA",
  colors_use = viridis_plasma_dark_high,
  na_color = "lightgray",
  order = TRUE,
  pt.size = NULL,
  aspect_ratio = NULL,
  reduction = NULL,
  na_cutoff = 1e-09,
  raster = NULL,
  raster.dpi = c(512, 512),
  slot = deprecated(),
  layer = "data",
  num_columns = NULL,
  alpha_exp = NULL,
  alpha_na_exp = NULL,
  ...
)
```

Arguments

seurat_object	Seurat object name.
features	Feature(s) to plot.
assay1	name of assay one. Default is "RAW" as featured in Create_CellBender_Merged_Seurat
assay2	name of assay two Default is "RNA" as featured in Create_CellBender_Merged_Seurat
colors_use	list of colors or color palette to use.
na_color	color to use for points below lower limit.
order	whether to move positive cells to the top (default = TRUE).
pt.size	Adjust point size for plotting.
aspect_ratio	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.

reduction	Dimensionality Reduction to use (if NULL then defaults to Object default).
na_cutoff	Value to use as minimum expression cutoff. To set no cutoff set to NA.
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
raster.dpi	Pixel resolution for rasterized plots, passed to geom_scattermore(). Default is c(512, 512).
slot	[Deprecated] soft-deprecated. See layer
layer	Which layer to pull expression data from? Default is "data".
num_columns	Number of columns in plot layout. If number of features > 1 then num_columns dictates the number of columns in overall layout (num_columns = 1 means stacked layout & num_columns = 2 means adjacent layout).
alpha_exp	new alpha level to apply to expressing cell color palette (colors_use). Must be value between 0-1.
alpha_na_exp	new alpha level to apply to non-expressing cell color palette (na_color). Must be value between 0-1.
...	Extra parameters passed to FeaturePlot .

Value

A ggplot object

Examples

```
## Not run:
FeaturePlot_DualAssay(seurat_object = object, features = "Cx3cr1", assay1 = "RAW", assay2 = "RNA",
  colors_use = viridis_plasma_dark_high, na_color = "lightgray")

## End(Not run)
```

FeaturePlot_scCustom *Customize FeaturePlot*

Description

Create Custom FeaturePlots and preserve scale (no binning)

Usage

```
FeaturePlot_scCustom(
  seurat_object,
  features,
  colors_use = viridis_plasma_dark_high,
  na_color = "lightgray",
  order = TRUE,
```



```

    pt.size = NULL,
    reduction = NULL,
    na_cutoff = 1e-09,
    raster = NULL,
    raster.dpi = c(512, 512),
    split.by = NULL,
    split_collect = NULL,
    aspect_ratio = NULL,
    figure_plot = FALSE,
    num_columns = NULL,
    slot = deprecated(),
    layer = "data",
    alpha_exp = NULL,
    alpha_na_exp = NULL,
    label = FALSE,
    label_feature_yaxis = FALSE,
    combine = TRUE,
    ...
)

```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>features</code>	Feature(s) to plot.
<code>colors_use</code>	list of colors or color palette to use.
<code>na_color</code>	color to use for points below lower limit.
<code>order</code>	whether to move positive cells to the top (default = TRUE).
<code>pt.size</code>	Adjust point size for plotting.
<code>reduction</code>	Dimensionality Reduction to use (if NULL then defaults to Object default).
<code>na_cutoff</code>	Value to use as minimum expression cutoff. This will be lowest value plotted use palette provided to <code>colors_use</code> . Leave as default value to plot only positive non-zero values using color scale and zero/negative values as NA. To plot all values using color palette set to NA.
<code>raster</code>	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
<code>raster.dpi</code>	Pixel resolution for rasterized plots, passed to <code>geom_scattermore()</code> . Default is <code>c(512, 512)</code> .
<code>split.by</code>	Variable in <code>@meta.data</code> to split the plot by.
<code>split_collect</code>	logical, whether to collect the legends/guides when plotting with <code>split.by</code> . Default is TRUE if one value is provided to <code>features</code> otherwise is set to FALSE.
<code>aspect_ratio</code>	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.
<code>figure_plot</code>	logical. Whether to remove the axes and plot with legend on left of plot denoting axes labels. (Default is FALSE). Requires <code>split_seurat = TRUE</code> .
<code>num_columns</code>	Number of columns in plot layout.

slot	[Deprecated] soft-deprecated. See layer
layer	Which layer to pull expression data from? Default is "data".
alpha_exp	new alpha level to apply to expressing cell color palette (colors_use). Must be value between 0-1.
alpha_na_exp	new alpha level to apply to non-expressing cell color palette (na_color). Must be value between 0-1.
label	logical, whether to label the clusters. Default is FALSE.
label_feature_yaxis	logical, whether to place feature labels on secondary y-axis as opposed to above legend key. Default is FALSE. When setting label_feature_yaxis = TRUE the number of columns in plot output will automatically be set to the number of levels in split.by'
combine	Combine plots into a single patchwork ed ggplot object. If FALSE, return a list of ggplot objects.
...	Extra parameters passed to FeaturePlot .

Value

A ggplot object

Examples

```
library(Seurat)
FeaturePlot_scCustom(seurat_object = pbmc_small, features = "CD3E",
  colors_use = viridis_plasma_dark_high, na_color = "lightgray")
```

FeatureScatter_scCustom

Modified version of FeatureScatter

Description

Create customized FeatureScatter plots with scCustomize defaults.

Usage

```
FeatureScatter_scCustom(
  seurat_object,
  feature1 = NULL,
  feature2 = NULL,
  colors_use = NULL,
  pt.size = NULL,
  group.by = NULL,
  split.by = NULL,
  split_seurat = FALSE,
```

```

    shuffle = TRUE,
    aspect_ratio = NULL,
    title_size = 15,
    plot.cor = TRUE,
    num_columns = NULL,
    raster = NULL,
    raster.dpi = c(512, 512),
    ggplot_default_colors = FALSE,
    color_seed = 123,
    ...
)

```

Arguments

seurat_object	Seurat object name.
feature1	First feature to plot.
feature2	Second feature to plot.
colors_use	color for the points on plot.
pt.size	Adjust point size for plotting.
group.by	Name of one or more metadata columns to group (color) cells by (for example, orig.ident). Default is active ident.
split.by	Feature to split plots by (i.e. "orig.ident").
split_seurat	logical. Whether or not to display split plots like Seurat (shared y axis) or as individual plots in layout. Default is FALSE.
shuffle	logical, whether to randomly shuffle the order of points. This can be useful for crowded plots if points of interest are being buried. Default is TRUE.
aspect_ratio	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.
title_size	size for plot title labels. Does NOT apply if split_seurat = TRUE.
plot.cor	Display correlation in plot subtitle (or title if split_seurat = TRUE).
num_columns	number of columns in final layout plot.
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
raster.dpi	Pixel resolution for rasterized plots, passed to geom_scattermore(). Default is c(512, 512).
ggplot_default_colors	logical. If colors_use = NULL, Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.
...	Extra parameters passed to FeatureScatter .

Value

A ggplot object

Examples

```
library(Seurat)

pbmc_small$sample_id <- sample(c("sample1", "sample2"), size = ncol(pbmc_small), replace = TRUE)

FeatureScatter_scCustom(seurat_object = pbmc_small, feature1 = "nCount_RNA",
feature2 = "nFeature_RNA", split.by = "sample_id")
```

Feature_Present	<i>Check if genes/features are present</i>
-----------------	--------------------------------------------

Description

Check if genes are present in object and return vector of found genes. Return warning messages for genes not found.

Usage

```
Feature_Present(
  data,
  features,
  case_check = TRUE,
  case_check_msg = TRUE,
  print_msg = TRUE,
  omit_warn = TRUE,
  return_none = FALSE,
  seurat_assay = NULL
)
```

Arguments

data	Name of input data. Currently only data of classes: Seurat, liger, data.frame, dgCMMatrix, dgTMatrix, tibble are accepted. Gene_IDs must be present in row-names of the data.
features	vector of features to check.
case_check	logical. Whether or not to check if features are found if the case is changed from the input list (Sentence case to Upper and vice versa). Default is TRUE.
case_check_msg	logical. Whether to print message to console if alternate case features are found in addition to inclusion in returned list. Default is TRUE.
print_msg	logical. Whether message should be printed if all features are found. Default is TRUE.
omit_warn	logical. Whether to print message about features that are not found in current object. Default is TRUE.

return_none logical. Whether list of found vs. bad features should still be returned if no features are found. Default is FALSE.

seurat_assay Name of assay to pull feature names from if data is Seurat Object. Default is NULL which will check against features from all assays present.

Value

A list of length 3 containing 1) found features, 2) not found features, 3) features found if case was modified.

Examples

```
## Not run:
features <- Feature_Present(data = obj_name, features = DEG_list, print_msg = TRUE,
case_check = TRUE)
found_features <- features[[1]]

## End(Not run)
```

Fetch_Meta

Get meta data from object

Description

Quick function to properly pull meta.data from objects.

Usage

```
Fetch_Meta(object)

## S3 method for class 'Seurat'
Fetch_Meta(object)

## S3 method for class 'liger'
Fetch_Meta(object)
```

Arguments

object Object of class Seurat or liger.

Value

A data.frame containing cell-level meta data

Examples

```
library(Seurat)
meta_data <- Fetch_Meta(object = pbmc_small)
head(meta_data, 5)
```

Gene_Present	<i>Check if genes/features are present</i> [Soft-deprecated]
--------------	---------------------------------------------------------------------

Description

Check if genes are present in object and return vector of found genes. Return warning messages for genes not found.

Usage

```
Gene_Present(
  data,
  gene_list,
  case_check = TRUE,
  case_check_msg = TRUE,
  print_msg = TRUE,
  omit_warn = TRUE,
  return_none = FALSE,
  seurat_assay = NULL
)
```

Arguments

<code>data</code>	Name of input data. Currently only data of classes: Seurat, liger, data.frame, dgCMatrx, dgTMatrix, tibble are accepted. Gene_IDs must be present in row-names of the data.
<code>gene_list</code>	vector of genes to check.
<code>case_check</code>	logical. Whether or not to check if features are found if the case is changed from the input list (Sentence case to Upper and vice versa). Default is TRUE.
<code>case_check_msg</code>	logical. Whether to print message to console if alternate case features are found in addition to inclusion in returned list. Default is TRUE.
<code>print_msg</code>	logical. Whether message should be printed if all features are found. Default is TRUE.
<code>omit_warn</code>	logical. Whether to print message about features that are not found in current object. Default is TRUE.
<code>return_none</code>	logical. Whether list of found vs. bad features should still be returned if no features are found. Default is FALSE.
<code>seurat_assay</code>	Name of assay to pull feature names from if data is Seurat Object. Default is NULL which will check against features from all assays present.

Value

A list of length 3 containing 1) found features, 2) not found features, 3) features found if case was modified.

Examples

```
## Not run:
features <- Gene_Present(data = obj_name, gene_list = DEG_list, print_msg = TRUE, case_check = TRUE)
found_features <- features[[1]]

## End(Not run)
```

Hue_Pal

Hue_Pal

Description

Shortcut to hue_pal to return to ggplot2 defaults if user desires, from scales package.

Usage

```
Hue_Pal(num_colors)
```

Arguments

num_colors number of colors to return in palette.

Value

hue color palette (as many colors as desired)

Examples

```
cols <- Hue_Pal(num_colors = 8)
PalettePlot(pal= cols)
```

ieg_gene_list *Immediate Early Gene (IEG) gene lists*

Description

Gene symbols for immediate early genes

Usage

ieg_gene_list

Format

A list of seven vectors

Mus_musculus_IEGs Gene symbols for IEGs from source publication (see below)

Homo_sapiens_IEGs Human gene symbols for homologous genes from mouse gene list

Source

Mouse gene list is from: SI Table 4 from [doi:10.1016/j.neuron.2017.09.026](https://doi.org/10.1016/j.neuron.2017.09.026). Human gene list was compiled by first creating homologous gene list using biomaRt and then adding some manually curated homologs according to HGNC.

Iterate_Barcode_Rank_Plot
Iterative Barcode Rank Plots

Description

Read data, calculate DropletUtils::barcodeRanks, create barcode rank plots, and outout single PDF output.

Usage

```
Iterate_Barcode_Rank_Plot(
  dir_path_h5,
  multi_directory = TRUE,
  h5_filename = "raw_feature_bc_matrix.h5",
  cellranger_multi = FALSE,
  parallel = FALSE,
  num_cores = NULL,
  file_path = NULL,
  file_name = NULL,
  pt.size = 6,
  raster_dpi = c(1024, 1024),
```



```

    plateau = NULL,
    ...
)

```

Arguments

<code>dir_path_h5</code>	path to parent directory (if <code>multi_directory = TRUE</code>) or directory containing all h5 files (if <code>multi_directory = FALSE</code>).
<code>multi_directory</code>	logical, whether or not all h5 files are in their own subdirectories or in a single directory (default is <code>TRUE</code> ; each in own subdirectory (e.g. output from Cell Ranger)).
<code>h5_filename</code>	Either the file name of h5 file (if <code>multi_directory = TRUE</code>) or the shared suffix (if <code>multi_directory = FALSE</code>)
<code>cellranger_multi</code>	logical, whether the outputs to be read are from Cell Ranger <code>multi</code> as opposed to Cell Ranger <code>count</code> (default is <code>FALSE</code>). Only valid if <code>multi_directory = FALSE</code> .
<code>parallel</code>	logical, should files be read in parallel (default is <code>FALSE</code>).
<code>num_cores</code>	Number of cores to use in parallel if <code>parallel = TRUE</code> .
<code>file_path</code>	file path to use for saving PDF output.
<code>file_name</code>	Name of PDF output file.
<code>pt.size</code>	point size for plotting, default is 6.
<code>raster_dpi</code>	Pixel resolution for rasterized plots, passed to <code>geom_scattermore()</code> . Default is <code>c(1024, 1024)</code> .
<code>plateau</code>	numerical values at which to add vertical line designating estimated empty droplet plateau (default is <code>NULL</code>). Must be vector equal in length to number of samples.
<code>...</code>	Additional parameters passed to <code>Read10X_h5_Multi_Directory</code> or <code>Read10X_h5_GEO</code> .

Value

pdf document

Examples

```

## Not run:
Iterate_Barcode_Rank_Plot(dir_path_h5 = "H5_PATH/", multi_directory = TRUE,
h5_filename = "raw_feature_bc_matrix", parallel = TRUE, num_cores = 12, file_path = "OUTPUT_PATH",
file_name = "Barcode_Rank_Plots")

## End(Not run)

```

 Iterate_Cluster_Highlight_Plot

Iterate Cluster Highlight Plot

Description

Iterate the create plots with cluster of interest highlighted across all cluster (active.idents) in given Seurat Object

Usage

```
Iterate_Cluster_Highlight_Plot(
  seurat_object,
  highlight_color = "dodgerblue",
  background_color = "lightgray",
  pt.size = NULL,
  reduction = NULL,
  file_path = NULL,
  file_name = NULL,
  file_type = NULL,
  single_pdf = FALSE,
  dpi = 600,
  raster = NULL,
  ...
)
```

Arguments

seurat_object	Seurat object name.
highlight_color	Color to highlight cells (default "navy"). Can provide either single color to use for all clusters/plots or a vector of colors equal to the number of clusters to use (in order) for the clusters/plots.
background_color	non-highlighted cell colors.
pt.size	point size for both highlighted cluster and background.
reduction	Dimensionality Reduction to use (if NULL then defaults to Object default).
file_path	directory file path and/or file name prefix. Defaults to current wd.
file_name	name suffix to append after sample name.
file_type	File type to save output as. Must be one of following: ".pdf", ".png", ".tiff", ".jpeg", or ".svg".
single_pdf	saves all plots to single PDF file (default = FALSE). 'file_type' must be .pdf
dpi	dpi for image saving.
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
...	Extra parameters passed to <code>DimPlot</code> .

Value

Saved plots

Examples

```
## Not run:
Iterate_Cluster_Highlight_Plot(seurat_object = object, highlight_color = "navy",
background_color = "lightgray", file_path = "path/", file_name = "name", file_type = "pdf",
single_pdf = TRUE)

## End(Not run)
```

Iterate_DimPlot_bySample

Iterate DimPlot By Sample

Description

Iterate DimPlot by orig.ident column from Seurat object metadata

Usage

```
Iterate_DimPlot_bySample(
  seurat_object,
  sample_column = "orig.ident",
  file_path = NULL,
  file_name = NULL,
  file_type = NULL,
  single_pdf = FALSE,
  dpi = 600,
  color = "black",
  no_legend = TRUE,
  title_prefix = NULL,
  reduction = NULL,
  dims = c(1, 2),
  pt.size = NULL,
  raster = NULL,
  ...
)
```

Arguments

seurat_object	Seurat object name.
sample_column	name of meta.data column containing sample names/ids (default is "orig.ident").
file_path	directory file path and/or file name prefix. Defaults to current wd.
file_name	name suffix to append after sample name.

file_type	File type to save output as. Must be one of following: ".pdf", ".png", ".tiff", ".jpeg", or ".svg".
single_pdf	saves all plots to single PDF file (default = FALSE). "file_type" must be .pdf
dpi	dpi for image saving.
color	color scheme to use.
no_legend	logical, whether or not to include plot legend, default is TRUE.
title_prefix	Value that should be used for plot title prefix if no_legend = TRUE. If NULL the value of meta_data_column will be used. Default is NULL.
reduction	Dimensionality Reduction to use (default is object default).
dims	Dimensions to plot.
pt.size	Adjust point size for plotting.
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
...	Extra parameters passed to DimPlot .

Value

A ggplot object

Examples

```
## Not run:
Iterate_DimPlot_bySample(seurat_object = object, file_path = "plots/", file_name = "tsne",
file_type = ".jpg", dpi = 600, color = "black")

## End(Not run)
```

Iterate_FeaturePlot_scCustom

Iterative Plotting of Gene Lists using Custom FeaturePlots

Description

Create and Save plots for Gene list with Single Command

Usage

```
Iterate_FeaturePlot_scCustom(
  seurat_object,
  features,
  gene_list = deprecated(),
  colors_use = viridis_plasma_dark_high,
  na_color = "lightgray",
  na_cutoff = 1e-09,
```

```

split.by = NULL,
order = TRUE,
return_plots = FALSE,
file_path = NULL,
file_name = NULL,
file_type = NULL,
single_pdf = FALSE,
features_per_page = 1,
num_columns = NULL,
landscape = TRUE,
dpi = 600,
pt.size = NULL,
reduction = NULL,
raster = NULL,
alpha_exp = NULL,
alpha_na_exp = NULL,
...
)

```

Arguments

seurat_object	Seurat object name.
features	vector of features to plot. If a named vector is provided then the names for each gene will be incorporated into plot title if <code>single_pdf = TRUE</code> or into file name if <code>FALSE</code> .
gene_list	[Deprecated] soft-deprecated. See features.
colors_use	color scheme to use.
na_color	color for non-expressed cells.
na_cutoff	Value to use as minimum expression cutoff. To set no cutoff set to NA.
split.by	Variable in <code>@meta.data</code> to split the plot by.
order	whether to move positive cells to the top (default = <code>TRUE</code>).
return_plots	logical. Whether to return plots to list instead of saving them to file(s). Default is <code>FALSE</code> .
file_path	directory file path and/or file name prefix. Defaults to current wd.
file_name	name suffix and file extension.
file_type	File type to save output as. Must be one of following: ".pdf", ".png", ".tiff", ".jpeg", or ".svg".
single_pdf	saves all plots to single PDF file (default = <code>FALSE</code>).
features_per_page	numeric, number of features to plot on single page if <code>single_pdf = TRUE</code> . Default is 1.
num_columns	Number of columns in plot layout (only applicable if <code>single_pdf = TRUE</code> AND <code>features_per_page > 1</code>).
landscape	logical, when plotting multiple features per page in single PDF whether to use landscape or portrait page dimensions (default is <code>TRUE</code>).

dpi	dpi for image saving.
pt.size	Adjust point size for plotting.
reduction	Dimensionality Reduction to use (if NULL then defaults to Object default).
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
alpha_exp	new alpha level to apply to expressing cell color palette (colors_use). Must be value between 0-1.
alpha_na_exp	new alpha level to apply to non-expressing cell color palette (na_color). Must be value between 0-1.
...	Extra parameters passed to FeaturePlot .

Value

Saved plots

Examples

```
## Not run:
Iterate_FeaturePlot_scCustom(seurat_object = object, gene_list = DEG_list,
  colors_use = viridis_plasma_dark_high, na_color = "lightgray", file_path = "plots/",
  file_name = "tsne", file_type = ".jpg", dpi = 600)

## End(Not run)
```

Iterate_Meta_Highlight_Plot
Iterate Meta Highlight Plot

Description

Iterate the create plots with meta data variable of interest highlighted.

Usage

```
Iterate_Meta_Highlight_Plot(
  seurat_object,
  meta_data_column,
  new_meta_order = NULL,
  meta_data_sort = TRUE,
  highlight_color = "navy",
  background_color = "lightgray",
  pt.size = NULL,
  no_legend = FALSE,
  title_prefix = NULL,
  reduction = NULL,
```

```

    file_path = NULL,
    file_name = NULL,
    file_type = NULL,
    single_pdf = FALSE,
    dpi = 600,
    raster = NULL,
    ...
)

```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>meta_data_column</code>	Name of the column in <code>seurat_object@meta.data</code> slot to pull value from for highlighting.
<code>new_meta_order</code>	The order in which to plot each level within <code>meta_data_column</code> if <code>single_PDF</code> is TRUE.
<code>meta_data_sort</code>	logical. Whether or not to sort and relevel the levels in <code>meta_data_column</code> if <code>single_PDF</code> is TRUE. Default is TRUE.
<code>highlight_color</code>	Color to highlight cells (default "navy"). Can provide either single color to use for all clusters/plots or a vector of colors equal to the number of clusters to use (in order) for the clusters/plots.
<code>background_color</code>	non-highlighted cell colors.
<code>pt.size</code>	point size for both highlighted cluster and background.
<code>no_legend</code>	logical, whether or not to remove plot legend and move to plot title. Default is FALSE.
<code>title_prefix</code>	Value that should be used for plot title prefix if <code>no_legend = TRUE</code> . If NULL the value of <code>meta_data_column</code> will be used. Default is NULL.
<code>reduction</code>	Dimensionality Reduction to use (if NULL then defaults to Object default).
<code>file_path</code>	directory file path and/or file name prefix. Defaults to current wd.
<code>file_name</code>	name suffix to append after sample name.
<code>file_type</code>	File type to save output as. Must be one of following: ".pdf", ".png", ".tiff", ".jpeg", or ".svg".
<code>single_pdf</code>	saves all plots to single PDF file (default = FALSE). "file_type" must be .pdf
<code>dpi</code>	dpi for image saving.
<code>raster</code>	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
<code>...</code>	Extra parameters passed to <code>DimPlot</code> .

Value

Saved plots

Examples

```
## Not run:
Iterate_Meta_Highlight_Plot(seurat_object = object, meta_data_column = "sample_id",
highlight_color = "navy", background_color = "lightgray", file_path = "path/",
file_name = "name", file_type = "pdf", single_pdf = TRUE)

## End(Not run)
```

Iterate_PC>Loading_Plots

Iterate PC>Loading_Plots

Description

Plot PC Heatmaps and Dim Loadings for exploratory analysis

Usage

```
Iterate_PC>Loading_Plots(
  seurat_object,
  dims_plot = NULL,
  file_path = NULL,
  name_prefix = NULL,
  file_name = "PC>Loading_Plots",
  return_plots = FALSE
)
```

Arguments

seurat_object	Seurat object name.
dims_plot	number of PCs to plot (integer). Default is all dims present in PCA.
file_path	directory file path to save file.
name_prefix	prefix for file name (optional).
file_name	suffix for file name. Default is "PC>Loading_Plots".
return_plots	Whether to return the plot list (Default is FALSE). Must assign to environment to save plot list.

Value

A list of plots outputted as pdf

See Also

[PCHeatmap](#) and [VizDimLoadings](#)

Examples

```
## Not run:
Iterate_PC>Loading_Plots(seurat_object = seurat, dims_plot = 25, file_path = "plots/")

## End(Not run)
```

Iterate_Plot_Density_Custom

Iterative Plotting of Gene Lists using Custom Density Plots

Description

Create and save plots for gene list with single command. Requires Nebulosa package from Bioconductor.

Usage

```
Iterate_Plot_Density_Custom(
  seurat_object,
  gene_list,
  viridis_palette = "magma",
  custom_palette = NULL,
  pt.size = 1,
  file_path = NULL,
  file_name = NULL,
  file_type = NULL,
  single_pdf = FALSE,
  dpi = 600,
  reduction = NULL,
  combine = TRUE,
  joint = FALSE,
  ...
)
```

Arguments

seurat_object	Seurat object name.
gene_list	vector of genes to plot. If a named vector is provided then the names for each gene will be incorporated into plot title if <code>single_pdf = TRUE</code> or into file name if <code>FALSE</code> .
viridis_palette	color scheme to use.
custom_palette	color for non-expressed cells.
pt.size	Adjust point size for plotting.
file_path	directory file path and/or file name prefix. Defaults to current wd.

file_name	name suffix and file extension.
file_type	File type to save output as. Must be one of following: ".pdf", ".png", ".tiff", ".jpeg", or ".svg".
single_pdf	saves all plots to single PDF file (default = FALSE). 'file_type' must be .pdf.
dpi	dpi for image saving.
reduction	Dimensionality Reduction to use (if NULL then defaults to Object default)
combine	Create a single plot? If FALSE, a list with ggplot objects is returned.
joint	NULL. This function only supports joint = FALSE. Leave as NULL to generate plots. To iterate joint plots see function: Iterate_Plot_Density_Joint.
...	Extra parameters passed to plot_density .

Value

Saved plots

Examples

```
## Not run:
Iterate_Plot_Density_Custom(seurat_object = object, gene_list = DEG_list, viridis_palette = "magma",
file_path = "plots/", file_name = "_density_plots", file_type = ".jpg", dpi = 600)

## End(Not run)
```

Iterate_Plot_Density_Joint

Iterative Plotting of Gene Lists using Custom Joint Density Plots

Description

Create and save plots for gene list with single command. Requires Nebulosa package from Bioconductor.

Usage

```
Iterate_Plot_Density_Joint(
  seurat_object,
  gene_list,
  viridis_palette = "magma",
  custom_palette = NULL,
  pt.size = 1,
  file_path = NULL,
  file_name = NULL,
  file_type = NULL,
  single_pdf = FALSE,
  dpi = 600,
```

```

    reduction = NULL,
    combine = TRUE,
    joint = NULL,
    ...
)

```

Arguments

`seurat_object` Seurat object name.

`gene_list` a list of vectors of genes to plot jointly. Each entry in the list will be plotted for the joint density. All entries in list must be greater than 2 features. If a named list is provided then the names for each list entry will be incorporated into plot title if `single_pdf = TRUE` or into file name if `FALSE`.

`viridis_palette` color scheme to use.

`custom_palette` color for non-expressed cells.

`pt.size` Adjust point size for plotting.

`file_path` directory file path and/or file name prefix. Defaults to current wd.

`file_name` name suffix and file extension.

`file_type` File type to save output as. Must be one of following: ".pdf", ".png", ".tiff", ".jpeg", or ".svg".

`single_pdf` saves all plots to single PDF file (default = FALSE). 'file_type' must be .pdf.

`dpi` dpi for image saving.

`reduction` Dimensionality Reduction to use (if NULL then defaults to Object default)

`combine` Create a single plot? If FALSE, a list with ggplot objects is returned.

`joint` NULL. This function only supports `joint = FALSE`. Leave as NULL to generate plots. To iterate joint plots see function: `Iterate_Plot_Density_Joint`.

`...` Extra parameters passed to `plot_density`.

Value

Saved plots

Examples

```

## Not run:
Iterate_Plot_Density_Joint(seurat_object = object, gene_list = DEG_list, viridis_palette = "magma",
file_path = "plots/", file_name = "joint_plots", file_type = ".jpg", dpi = 600)

## End(Not run)

```

 Iterate_VlnPlot_scCustom

Iterative Plotting of Gene Lists using VlnPlot_scCustom

Description

Create and Save plots for Gene list with Single Command

Usage

```
Iterate_VlnPlot_scCustom(
  seurat_object,
  features,
  gene_list = deprecated(),
  colors_use = NULL,
  pt.size = NULL,
  group.by = NULL,
  split.by = NULL,
  file_path = NULL,
  file_name = NULL,
  file_type = NULL,
  single_pdf = FALSE,
  raster = NULL,
  dpi = 600,
  ggplot_default_colors = FALSE,
  color_seed = 123,
  ...
)
```

Arguments

seurat_object	Seurat object name.
features	vector of features to plot.
gene_list	[Deprecated] soft-deprecated. See features.
colors_use	color palette to use for plotting. By default if number of levels plotted is less than or equal to 36 it will use "polychrome" and if greater than 36 will use "varibow" with shuffle = TRUE both from DiscretePalette_scCustomize.
pt.size	point size for plotting.
group.by	Name of one or more metadata columns to group (color) plot by (for example, orig.ident); default is the current active.ident of the object.
split.by	Feature to split plots by (i.e. "orig.ident").
file_path	directory file path and/or file name prefix. Defaults to current wd.
file_name	name suffix and file extension.

file_type	File type to save output as. Must be one of following: ".pdf", ".png", ".tiff", ".jpeg", or ".svg".
single_pdf	saves all plots to single PDF file (default = FALSE). ‘file_type‘ must be .pdf.
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 100,000 total points plotted (# Cells x # of features).
dpi	dpi for image saving.
ggplot_default_colors	logical. If colors_use = NULL, Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.
...	Extra parameters passed to VlnPlot .

Value

Saved plots

Examples

```
## Not run:
Iterate_VlnPlot_scCustom(seurat_object = object, gene_list = DEG_list, colors = color_list,
file_path = "plots/", file_name = "_vln", file_type = ".jpg", dpi = 600)

## End(Not run)
```

JCO_Four

Four Color Palette (JCO)

Description

Shortcut to a specific JCO 4 color palette from ggsci package.

Usage

```
JCO_Four()
```

Value

4 color palette from the JCO ggsci palette

References

Selection of colors from the JCO palette from ggsci being called through paletteer. See ggsci for more info on palettes <https://CRAN.R-project.org/package=ggsci>

Examples

```
cols <- JCO_Four()
PalettePlot(pal= cols)
```

LIGER_Features	<i>Extract Features from LIGER Object</i>
----------------	-------------------------------------------

Description

Extract all unique features from LIGER object

Usage

```
LIGER_Features(liger_object, by_dataset = FALSE)
```

Arguments

liger_object	LIGER object name.
by_dataset	logical, whether to return list with vector of features for each dataset in LIGER object or to return single vector of unique features across all datasets in object (default is FALSE; return vector of unique features)

Value

vector or list depending on by_dataset parameter

Examples

```
## Not run:
# return single vector of all unique features
all_features <- LIGER_Features(liger_object = object, by_dataset = FALSE)

# return list of vectors containing features from each individual dataset in object
dataset_features <- LIGER_Features(liger_object = object, by_dataset = TRUE)

## End(Not run)
```

Liger_to_Seurat	<i>Create a Seurat object containing the data from a liger object</i> [Soft-deprecated]
-----------------	--------------------------------------------------------------------------------------------------

Description

Merges raw.data and scale.data of object, and creates Seurat object with these values along with tsne.coords, iNMF factorization, and cluster assignments. Supports Seurat V2 and V3.

Usage

```
Liger_to_Seurat(
  liger_object,
  nms = names(liger_object@H),
  renormalize = TRUE,
  use.liger.genes = TRUE,
  by.dataset = FALSE,
  keep_meta = TRUE,
  reduction_label = "UMAP",
  seurat_assay = "RNA",
  assay_type = NULL,
  add_barcode_names = FALSE,
  barcode_prefix = TRUE,
  barcode_cell_id_delimiter = "_"
)
```

Arguments

liger_object	liger object.
nms	By default, labels cell names with dataset of origin (this is to account for cells in different datasets which may have same name). Other names can be passed here as vector, must have same length as the number of datasets. (default names(H)).
renormalize	Whether to log-normalize raw data using Seurat defaults (default TRUE).
use.liger.genes	Whether to carry over variable genes (default TRUE).
by.dataset	Include dataset of origin in cluster identity in Seurat object (default FALSE).
keep_meta	logical. Whether to transfer additional metadata (nGene/nUMI/dataset already transferred) to new Seurat Object. Default is TRUE.
reduction_label	Name of dimensionality reduction technique used. Enables accurate transfer or name to Seurat object instead of defaulting to "tSNE".
seurat_assay	Name to set for assay in Seurat Object. Default is "RNA".
assay_type	what type of Seurat assay to create in new object (Assay vs Assay5). Default is NULL which will default to the current user settings. See Convert_Assay parameter convert_to for acceptable values.

`add_barcode_names` logical, whether to add dataset names to the cell barcodes when creating Seurat object, default is FALSE.

`barcode_prefix` logical, if `add_barcode_names = TRUE` should the names be added as prefix to current cell barcodes/names or a suffix (default is TRUE; prefix).

`barcode_cell_id_delimiter` The delimiter to use when adding dataset id to barcode prefix/suffix. Default is "_".

Details

Stores original dataset identity by default in new object metadata if dataset names are passed in `nms`. iNMF factorization is stored in `dim.reduction` object with key "iNMF".

Value

Seurat object with `raw.data`, `scale.data`, `reduction_label`, `iNMF`, and `ident` slots set.

References

Original function is part of LIGER package <https://github.com/welch-lab/liger> (Licence: GPL-3). Function was slightly modified for use in `scCustomize` with `keep.meta` parameter. Also posted as PR to `liger` GitHub.

Examples

```
## Not run:
seurat_object <- Liger_to_Seurat(liger_object = LIGER_OBJ, reduction_label = "UMAP")

## End(Not run)
```

MAD_Stats

Median Absolute Deviation Statistics

Description

Get quick values for X x median absolute deviation for Genes, UMIs, %mito per cell grouped by `meta.data` variable.

Usage

```
MAD_Stats(
  seurat_object,
  group_by_var = "orig.ident",
  default_var = TRUE,
  mad_var = NULL,
  mad_num = 2
)
```


Arguments

seurat_object	Seurat object name.
group_by_var	Column in meta.data slot to group results by (default = "orig.ident").
default_var	logical. Whether to include the default meta.data variables of: "nCount_RNA", "nFeature_RNA", "percent_mito", "percent_ribo", "percent_mito_ribo", and "log10GenesPerUMI" in addition to variables supplied to mad_var.
mad_var	Column(s) in @meta.data to calculate medians for in addition to defaults. Must be of class() integer or numeric.
mad_num	integer value to multiply the MAD in returned data.frame (default is 2). Often helpful when calculating a outlier range to base of of median + (X*MAD).

Value

A data.frame.

Examples

```
## Not run:
mad_stats <- MAD_Stats(seurat_object = obj, group_by_var = "orig.ident")

## End(Not run)
```

Median_Stats

Median Statistics

Description

Get quick values for median Genes, UMIs, %mito per cell grouped by meta.data variable.

Usage

```
Median_Stats(
  seurat_object,
  group_by_var = "orig.ident",
  default_var = TRUE,
  median_var = NULL
)
```

Arguments

seurat_object	Seurat object name.
group_by_var	Column in meta.data slot to group results by (default = "orig.ident").
default_var	logical. Whether to include the default meta.data variables of: "nCount_RNA", "nFeature_RNA", "percent_mito", "percent_ribo", "percent_mito_ribo", and "log10GenesPerUMI" in addition to variables supplied to median_var.
median_var	Column(s) in @meta.data to calculate medians for in addition to defaults. Must be of class() integer or numeric.

Value

A data.frame.

Examples

```
## Not run:
med_stats <- Median_Stats(seurat_object - obj, group_by_var = "orig.ident")

## End(Not run)
```

<code>Merge_Seurat_List</code>	<i>Merge a list of Seurat Objects</i>
--------------------------------	---------------------------------------

Description

Enables easy merge of a list of Seurat Objects. See [merge](#) for more information,

Usage

```
Merge_Seurat_List(
  list_seurat,
  add.cell.ids = NULL,
  merge.data = TRUE,
  project = "SeuratProject"
)
```

Arguments

<code>list_seurat</code>	list composed of multiple Seurat Objects.
<code>add.cell.ids</code>	A character vector of equal length to the number of objects in <code>list_seurat</code> . Appends the corresponding values to the start of each objects' cell names. See merge .
<code>merge.data</code>	Merge the data slots instead of just merging the counts (which requires renormalization). This is recommended if the same normalization approach was applied to all objects. See merge .
<code>project</code>	Project name for the Seurat object. See merge .

Value

A Seurat Object

Examples

```
## Not run:
object_list <- list(obj1, obj2, obj3, ...)
merged_object <- Merge_Seurat_List(list_seurat = object_list)

## End(Not run)
```

Merge_Sparse_Data_All *Merge a list of Sparse Matrices*

Description

Enables easy merge of a list of sparse matrices

Usage

```
Merge_Sparse_Data_All(
  matrix_list,
  add_cell_ids = NULL,
  prefix = TRUE,
  cell_id_delimiter = "_"
)
```

Arguments

<code>matrix_list</code>	list of matrices to merge.
<code>add_cell_ids</code>	a vector of sample ids to add as prefix to cell barcode during merge.
<code>prefix</code>	logical. Whether <code>add_cell_ids</code> should be added as prefix to current cell barcodes/names or as suffix to current cell barcodes/names. Default is TRUE, add as prefix.
<code>cell_id_delimiter</code>	The delimiter to use when adding cell id prefix/suffix. Default is "_".

Value

A sparse Matrix

References

Original function is part of LIGER package as non-exported function <https://github.com/welch-lab/liger/blob/master/R/utilities.R> (License: GPL-3). Function was modified for use in scCustomize (add progress bar, prefix vs. suffix, and delimiter options).

Examples

```
## Not run:
data_list <- Read10X_GEO(...)
merged <- Merge_Sparse_Data_All(matrix_list = data_list, add_cell_ids = names(data_list),
prefix = TRUE, cell_id_delimiter = "_")

## End(Not run)
```

Merge_Sparse_Multimodal_All

Merge a list of Sparse Matrices contain multi-modal data.

Description

Enables easy merge of a list of sparse matrices for multi-modal data.

Usage

```
Merge_Sparse_Multimodal_All(
  matrix_list,
  add_cell_ids = NULL,
  prefix = TRUE,
  cell_id_delimiter = "_"
)
```

Arguments

<code>matrix_list</code>	list of matrices to merge.
<code>add_cell_ids</code>	a vector of sample ids to add as prefix to cell barcode during merge.
<code>prefix</code>	logical. Whether <code>add_cell_ids</code> should be added as prefix to current cell barcodes/names or as suffix to current cell barcodes/names. Default is TRUE, add as prefix.
<code>cell_id_delimiter</code>	The delimiter to use when adding cell id prefix/suffix. Default is "_".

Value

A list containing one sparse matrix for each modality

Examples

```
## Not run:
data_list <- Read10X_GEO(...)
merged_list <- Merge_Sparse_Multimodal_All(matrix_list = data_list, add_cell_ids = names(data_list),
prefix = TRUE, cell_id_delimiter = "_")

## End(Not run)
```

Meta_Highlight_Plot *Meta Highlight Plot*

Description

Create Plot with meta data variable of interest highlighted

Usage

```
Meta_Highlight_Plot(
  seurat_object,
  meta_data_column,
  meta_data_highlight,
  highlight_color = NULL,
  background_color = "lightgray",
  pt.size = NULL,
  aspect_ratio = NULL,
  figure_plot = FALSE,
  raster = NULL,
  raster.dpi = c(512, 512),
  label = FALSE,
  split.by = NULL,
  split_seurat = FALSE,
  ggplot_default_colors = FALSE,
  ...
)
```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>meta_data_column</code>	Name of the column in <code>seurat_object@meta.data</code> slot to pull value from for highlighting.
<code>meta_data_highlight</code>	Name of variable(s) within <code>meta_data_name</code> to highlight in the plot.
<code>highlight_color</code>	Color to highlight cells (default "navy").
<code>background_color</code>	non-highlighted cell colors.
<code>pt.size</code>	point size for both highlighted cluster and background.
<code>aspect_ratio</code>	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.
<code>figure_plot</code>	logical. Whether to remove the axes and plot with legend on left of plot denoting axes labels. (Default is FALSE). Requires <code>split_seurat = TRUE</code> .

raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
raster.dpi	Pixel resolution for rasterized plots, passed to <code>geom_scattermore()</code> . Default is <code>c(512, 512)</code> .
label	Whether to label the highlighted meta data variable(s). Default is FALSE.
split.by	Variable in <code>@meta.data</code> to split the plot by.
split_seurat	logical. Whether or not to display split plots like Seurat (shared y axis) or as individual plots in layout. Default is FALSE.
ggplot_default_colors	logical. If <code>colors_use = NULL</code> , Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
...	Extra parameters passed to <code>DimPlot</code> .

Value

A ggplot object

Examples

```
library(Seurat)
pbmc_small$sample_id <- sample(c("sample1", "sample2"), size = ncol(pbmc_small), replace = TRUE)

Meta_Highlight_Plot(seurat_object = pbmc_small, meta_data_column = "sample_id",
  meta_data_highlight = "sample1", highlight_color = "gold", background_color = "lightgray",
  pt.size = 2)
```

Meta_Numeric

Check if meta data columns are numeric

Description

Check if any present meta data columns are numeric and returns vector of valid numeric columns. Issues warning message if any columns not in numeric form.

Usage

```
Meta_Numeric(data)
```

Arguments

data a data.frame contain meta.data.

Value

vector of meta data columns that are numeric.

Examples

```
## Not run:  
numeric_meta_columns <- Meta_Numeric(data = meta_data)  
  
## End(Not run)
```

Meta_Present	<i>Check if meta data are present</i>
--------------	---------------------------------------

Description

Check if meta data columns are present in object and return vector of found columns Return warning messages for meta data columns not found.

Usage

```
Meta_Present(  
  object,  
  seurat_object = deprecated(),  
  meta_col_names,  
  print_msg = TRUE,  
  omit_warn = TRUE,  
  return_none = FALSE  
)
```

Arguments

object	Seurat or Liger object name.
seurat_object	[Deprecated] deprecated. Please use object instead.
meta_col_names	vector of column names to check.
print_msg	logical. Whether message should be printed if all features are found. Default is TRUE.
omit_warn	logical. Whether to print message about features that are not found in current object. Default is TRUE.
return_none	logical. Whether list of found vs. bad features should still be returned if no meta_col_names are found. Default is FALSE.

Value

vector of meta data columns that are present

Examples

```
## Not run:
meta_variables <- Meta_Present(object = obj_name, meta_col_names = "percent_mito", print_msg = TRUE)

## End(Not run)
```

Meta_Remove_Seurat *Remove meta data columns containing Seurat Defaults*

Description

Remove any columns from new meta_data data.frame in preparation for adding back to Seurat Object

Usage

```
Meta_Remove_Seurat(
  meta_data,
  seurat_object,
  barcodes_to_rownames = FALSE,
  barcodes_colname = "barcodes"
)
```

Arguments

meta_data data.frame containing meta data.
seurat_object object name.
barcodes_to_rownames
 logical, are barcodes present as column and should they be moved to rownames
 (to be compatible with Seurat::AddMetaData). Default is FALSE.
barcodes_colname
 name of barcodes column in meta_data. Required if barcodes_to_rownames =
 TRUE.

Value

data.frame with only new columns.

Examples

```
## Not run:
new_meta <- Meta_Remove_Seurat(meta_data = meta_data_df, seurat_object = object)
object <- AddMetaData(object = object, metadata = new_meta)

## End(Not run)
```

Move_Legend	<i>Move Legend Position</i>
-------------	-----------------------------

Description

Shortcut for thematic modification to move legend position.

Usage

```
Move_Legend(position = "right", ...)
```

Arguments

position	valid position to move legend. Default is "right".
...	extra arguments passed to <code>ggplot2::theme()</code> .

Value

Returns a list-like object of class *theme*.

Examples

```
# Generate a plot and customize theme
library(ggplot2)
df <- data.frame(x = rnorm(n = 100, mean = 20, sd = 2), y = rbinom(n = 100, size = 100, prob = 0.2))
p <- ggplot(data = df, mapping = aes(x = x, y = y)) + geom_point(mapping = aes(color = 'red'))
p + Move_Legend("left")
```

msigdb_qc_gene_list	<i>QC Gene Lists</i>
---------------------	----------------------

Description

Gene symbols for qc percentages from MSigDB database. The gene sets are from 3 MSigDB lists: "HALLMARK_OXIDATIVE_PHOSPHORYLATION", "HALLMARK_APOPTOSIS", and "HALLMARK_DNA_REPAIR".

Usage

```
msigdb_qc_gene_list
```

Format

A list of 18 vectors

Homo_sapiens_msigdb_oxphos Genes in msigdb "HALLMARK_OXIDATIVE_PHOSPHORYLATION" list for human

Homo_sapiens_msigdb_apop Genes in msigdb "HALLMARK_APOPTOSIS" list for human

Homo_sapiens_msigdb_dna_repair Genes in msigdb "HALLMARK_DNA_REPAIR" list for human

Mus_musculus_msigdb_oxphos Genes in msigdb "HALLMARK_OXIDATIVE_PHOSPHORYLATION" list for mouse

Mus_musculus_msigdb_apop Genes in msigdb "HALLMARK_APOPTOSIS" list for mouse

Mus_musculus_msigdb_dna_repair Genes in msigdb "HALLMARK_DNA_REPAIR" list for mouse

Rattus_norvegicus_msigdb_oxphos Genes in msigdb "HALLMARK_OXIDATIVE_PHOSPHORYLATION" list for rat

Rattus_norvegicus_msigdb_apop Genes in msigdb "HALLMARK_APOPTOSIS" list for rat

Rattus_norvegicus_msigdb_dna_repair Genes in msigdb "HALLMARK_DNA_REPAIR" list for rat

Drosophila_melanogaster_msigdb_oxphos Genes in msigdb "HALLMARK_OXIDATIVE_PHOSPHORYLATION" list for fly

Drosophila_melanogaster_msigdb_apop Genes in msigdb "HALLMARK_APOPTOSIS" list for fly

Drosophila_melanogaster_msigdb_dna_repair Genes in msigdb "HALLMARK_DNA_REPAIR" list for fly

Dario_erio_msigdb_oxphos Genes in msigdb "HALLMARK_OXIDATIVE_PHOSPHORYLATION" list for zebrafish

Dario_erio_msigdb_apop Genes in msigdb "HALLMARK_APOPTOSIS" list for zebrafish

Dario_erio_msigdb_dna_repair Genes in msigdb "HALLMARK_DNA_REPAIR" list for zebrafish

Macaca_mulatta_msigdb_oxphos Genes in msigdb "HALLMARK_OXIDATIVE_PHOSPHORYLATION" list for macaque

Macaca_mulatta_msigdb_apop Genes in msigdb "HALLMARK_APOPTOSIS" list for macaque

Macaca_mulatta_msigdb_dna_repair Genes in msigdb "HALLMARK_DNA_REPAIR" list for macaque

Source

MSigDB gene sets via msigdb package <https://cran.r-project.org/package=msigdb>

NavyAndOrange	<i>Navy and Orange Dual Color Palette</i>
---------------	-------------------------------------------

Description

Shortcut to navy orange color plot

Usage

```
NavyAndOrange(flip_order = FALSE)
```

Arguments

`flip_order` Whether to flip the order of colors.

Value

Navy orange palette

Examples

```
cols <- NavyAndOrange()  
PalettePlot(pal= cols)
```

PalettePlot	<i>Plot color palette in viewer</i>
-------------	-------------------------------------

Description

Plots given color vector/palette in viewer to evaluate palette before plotting on data.

Usage

```
PalettePlot(pal = NULL, label_color_num = NULL)
```

Arguments

`pal` a vector of colors (either named colors or hex codes).

`label_color_num`

logical, whether or not to numerically label the colors in output plot. Default is TRUE is number of colors in `pal` is less than 75 and FALSE is greater than 75.

Value

Plot of all colors in supplied palette/vector

References

Adapted from colorway package build_palette internals (License: GPL-3). <https://github.com/hypercompetent/colorway>.

Examples

```
pal <- DiscretePalette_scCustomize(num_colors = 36, palette = "varibow")
PalettePlot(pal = pal)
```

PC_Plotting

PC Plots

Description

Plot PC Heatmaps and Dim Loadings for exploratory analysis. Plots a single Heatmap and Gene Loading Plot. Used for PC>Loading_Plots function.

Usage

```
PC_Plotting(seurat_object, dim_number)
```

Arguments

seurat_object Seurat Object.
dim_number A single dim to plot (integer).

Value

A plot of PC heatmap and gene loadings for single

See Also

[PCHeatmap](#) and [VizDimLoadings](#)

Examples

```
library(Seurat)
PC_Plotting(seurat_object = pbmc_small, dim_number = 1)
```

Percent_Expressing *Calculate percent of expressing cells*

Description

Calculates the percent of cells that express a given set of features by various grouping factors

Usage

```
Percent_Expressing(  
  seurat_object,  
  features,  
  threshold = 0,  
  group_by = NULL,  
  split_by = NULL,  
  entire_object = FALSE,  
  slot = deprecated(),  
  layer = "data",  
  assay = NULL  
)
```

Arguments

seurat_object	Seurat object name.
features	Feature(s) to plot.
threshold	Expression threshold to use for calculation of percent expressing (default is 0).
group_by	Factor to group the cells by.
split_by	Factor to split the groups by.
entire_object	logical (default = FALSE). Whether to calculate percent of expressing cells across the entire object as opposed to by cluster or by group_by variable.
slot	[Deprecated] soft-deprecated. See layer
layer	Which layer to pull expression data from? Default is "data".
assay	Assay to pull feature data from. Default is active assay.

Value

A data.frame

References

Part of code is modified from Seurat package as used by [DotPlot](#) to generate values to use for plotting. Source code can be found here: <https://github.com/satijalab/seurat/blob/4e868fcde49dc0a3df47f94f5fb54R/visualization.R#L3391> (License: GPL-3).

Examples

```
## Not run:
percent_stats <- Percent_Expressing(seurat_object = object, features = "Cx3cr1", threshold = 0)

## End(Not run)
```

plotFactors_scCustom *Customized version of plotFactors*

Description

Modified and optimized version of plotFactors function from LIGER package.

Usage

```
plotFactors_scCustom(
  liger_object,
  num_genes = 8,
  colors_use_factors = NULL,
  colors_use_dimreduc = c("lemonchiffon", "red"),
  pt.size_factors = 1,
  pt.size_dimreduc = 1,
  reduction_label = "UMAP",
  plot_legend = TRUE,
  raster = TRUE,
  raster.dpi = c(512, 512),
  order = FALSE,
  plot_dimreduc = TRUE,
  save_plots = TRUE,
  file_path = NULL,
  file_name = NULL,
  return_plots = FALSE,
  cells.highlight = NULL,
  reorder_datasets = NULL,
  ggplot_default_colors = FALSE,
  color_seed = 123
)
```

Arguments

liger_object liger liger_object. Need to perform clustering and factorization before calling this function

num_genes Number of genes to display for each factor (Default 8).

colors_use_factors colors to use for plotting factor loadings By default datasets will be plotted using "varibow" with shuffle = TRUE from both from [DiscretePalette_scCustomize](#).

colors_use_dimreduc	colors to use for plotting factor loadings on dimensionality reduction coordinates (tSNE/UMAP). Default is c('lemonchiffon', 'red'),
pt.size_factors	Adjust point size for plotting in the factor plots.
pt.size_dimreduc	Adjust point size for plotting in dimensionality reduction plots.
reduction_label	What to label the x and y axes of resulting plots. LIGER does not store name of technique and therefore needs to be set manually. Default is "UMAP".
plot_legend	logical, whether to plot the legend on factor loading plots, default is TRUE. Helpful if number of datasets is large to avoid crowding the plot with legend.
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 200,000 cells.
raster.dpi	Pixel resolution for rasterized plots, passed to geom_scattermore(). Default is c(512, 512).
order	logical. Whether to plot higher loading cells on top of cells with lower loading values in the dimensionality reduction plots (Default = FALSE).
plot_dimreduc	logical. Whether to plot factor loadings on dimensionality reduction coordinates. Default is TRUE.
save_plots	logical. Whether to save plots. Default is TRUE
file_path	directory file path and/or file name prefix. Defaults to current wd.
file_name	name suffix to append after sample name.
return_plots	logical. Whether or not to return plots to the environment. (Default is FALSE)
cells.highlight	Names of specific cells to highlight in plot (black) (default NULL).
reorder_datasets	New order to plot datasets in for the factor plots if different from current factor level order in cell.data slot.
ggplot_default_colors	logical. If colors_use_factors = NULL, Whether or not to return plot using default ggplot2 "hue" palette instead of default "varibow" palette.
color_seed	random seed for the palette shuffle if colors_use_factors = NULL. Default = 123.

Value

A list of ggplot/patchwork objects and/or PDF file.

Author(s)

Velina Kozareva (Original code for modified function), Sam Marsh (Added/modified functionality)

References

Based on plotFactors functionality from original LIGER package.

Examples

```
## Not run:
plotFactors_scCustom(liger_object = liger_obj, return_plots = FALSE, plot_dimreduc = TRUE,
raster = FALSE, save_plots = TRUE)

## End(Not run)
```

Plot_Cells_per_Sample *Plot Number of Cells/Nuclei per Sample*

Description

Plot of total cell or nuclei number per sample grouped by another meta data variable.

Usage

```
Plot_Cells_per_Sample(
  seurat_object,
  sample_col = "orig.ident",
  group_by = NULL,
  colors_use = NULL,
  dot_size = 1,
  plot_title = "Cells/Nuclei per Sample",
  y_axis_label = "Number of Cells",
  x_axis_label = NULL,
  legend_title = NULL,
  x_lab_rotate = TRUE,
  color_seed = 123
)
```

Arguments

seurat_object	Seurat object name.
sample_col	Specify which column in meta.data specifies sample ID (i.e. orig.ident).
group_by	Column in meta.data slot to group results by (i.e. "Treatment").
colors_use	List of colors or color palette to use.
dot_size	size of the dots plotted if group_by is not NULL. Default is 1.
plot_title	Plot title.
y_axis_label	Label for y axis.
x_axis_label	Label for x axis.
legend_title	Label for plot legend.
x_lab_rotate	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.

Value

A ggplot object

Examples

```
## Not run:  
Plot_Cells_per_Sample(seurat_object = obj, sample_col = "orig.ident", group_by = "Treatment")  
  
## End(Not run)
```

Plot_Density_Custom *Nebulosa Density Plot*

Description

Allow for customization of Nebulosa plot_density. Requires Nebulosa package from Bioconductor.

Usage

```
Plot_Density_Custom(  
  seurat_object,  
  features,  
  joint = FALSE,  
  viridis_palette = "magma",  
  custom_palette = NULL,  
  pt.size = 1,  
  aspect_ratio = NULL,  
  reduction = NULL,  
  combine = TRUE,  
  ...  
)
```

Arguments

seurat_object	Seurat object name.
features	Features to plot.
joint	logical. Whether to return joint density plot. Default is FALSE.
viridis_palette	default viridis palette to use (must be one of: "viridis", "magma", "cividis", "inferno", "plasma"). Default is "magma".
custom_palette	non-default color palette to be used in place of default viridis options.
pt.size	Adjust point size for plotting.
aspect_ratio	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.

reduction	Dimensionality Reduction to use (if NULL then defaults to Object default).
combine	Create a single plot? If FALSE, a list with ggplot objects is returned.
...	Extra parameters passed to <code>plot_density</code> .

Value

A ggplot object

Examples

```
## Not run:
library(Seurat)
Plot_Density_Custom(seurat_object = pbmc_small, features = "CD3E")

## End(Not run)
```

Plot_Density_Joint_Only

Nebulosa Joint Density Plot

Description

Return only the joint density plot from Nebulosa `plot_density` function. Requires Nebulosa package from Bioconductor.

Usage

```
Plot_Density_Joint_Only(
  seurat_object,
  features,
  viridis_palette = "magma",
  custom_palette = NULL,
  pt.size = 1,
  aspect_ratio = NULL,
  reduction = NULL,
  ...
)
```

Arguments

seurat_object	Seurat object name.
features	Features to plot.
viridis_palette	default viridis palette to use (must be one of: "viridis", "magma", "cividis", "inferno", "plasma"). Default is "magma".
custom_palette	non-default color palette to be used in place of default viridis options.

pt.size	Adjust point size for plotting.
aspect_ratio	Control the aspect ratio (y:x axes ratio length). Must be numeric value; Default is NULL.
reduction	Dimensionality Reduction to use (if NULL then defaults to Object default).
...	Extra parameters passed to <code>plot_density</code> .

Value

A ggplot object

Examples

```
## Not run:
library(Seurat)
Plot_Density_Joint_Only(seurat_object = pbmc_small, features = c("CD8A", "CD3E"))

## End(Not run)
```

Plot_Median_Genes	<i>Plot Median Genes per Cell per Sample</i>
-------------------	----------------------------------------------

Description

Plot of median genes per cell per sample grouped by desired meta data variable.

Usage

```
Plot_Median_Genes(
  seurat_object,
  sample_col = "orig.ident",
  group_by = NULL,
  colors_use = NULL,
  dot_size = 1,
  plot_title = "Median Genes/Cell per Sample",
  y_axis_label = "Median Genes",
  x_axis_label = NULL,
  legend_title = NULL,
  x_lab_rotate = TRUE,
  color_seed = 123
)
```

Arguments

seurat_object	Seurat object name.
sample_col	Specify which column in meta.data specifies sample ID (i.e. orig.ident).
group_by	Column in meta.data slot to group results by (i.e. "Treatment").
colors_use	List of colors or color palette to use. Only applicable if group_by is not NULL.
dot_size	size of the dots plotted if group_by is not NULL. Default is 1.
plot_title	Plot title.
y_axis_label	Label for y axis.
x_axis_label	Label for x axis.
legend_title	Label for plot legend.
x_lab_rotate	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.

Value

A ggplot object

Examples

```
library(Seurat)
# Create example groups
pbmc_small$sample_id <- sample(c("sample1", "sample2"), size = ncol(pbmc_small), replace = TRUE)

# Plot
Plot_Median_Genes(seurat_object = pbmc_small, sample_col = "orig.ident", group_by = "sample_id")
```

Plot_Median_Mito

Plot Median Percent Mito per Cell per Sample

Description

Plot of median percent mito per cell per sample grouped by desired meta data variable.

Usage

```
Plot_Median_Mito(
  seurat_object,
  sample_col = "orig.ident",
  group_by = NULL,
  colors_use = NULL,
  dot_size = 1,
  plot_title = "Median % Mito per Sample",
```

```
  y_axis_label = "Percent Mitochondrial Reads",
  x_axis_label = NULL,
  legend_title = NULL,
  x_lab_rotate = TRUE,
  color_seed = 123
)
```

Arguments

seurat_object	Seurat object name.
sample_col	Specify which column in meta.data specifies sample ID (i.e. orig.ident).
group_by	Column in meta.data slot to group results by (i.e. "Treatment").
colors_use	List of colors or color palette to use. Only applicable if group_by is not NULL.
dot_size	size of the dots plotted if group_by is not NULL. Default is 1.
plot_title	Plot title.
y_axis_label	Label for y axis.
x_axis_label	Label for x axis.
legend_title	Label for plot legend.
x_lab_rotate	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.

Value

A ggplot object

Examples

```
## Not run:
# Add mito
obj <- Add_Mito_Ribo_Seurat(seurat_object = obj, species = "human")

# Plot
Plot_Median_Mito(seurat_object = obj, sample_col = "orig.ident", group_by = "sample_id")

## End(Not run)
```

Plot_Median_Other *Plot Median other variable per Cell per Sample*

Description

Plot of median other variable per cell per sample grouped by desired meta data variable.

Usage

```
Plot_Median_Other(
  seurat_object,
  median_var,
  sample_col = "orig.ident",
  group_by = NULL,
  colors_use = NULL,
  dot_size = 1,
  plot_title = NULL,
  y_axis_label = NULL,
  x_axis_label = NULL,
  legend_title = NULL,
  x_lab_rotate = TRUE,
  color_seed = 123
)
```

Arguments

seurat_object	Seurat object name.
median_var	Variable in meta.data slot to calculate and plot median values for.
sample_col	Specify which column in meta.data specifies sample ID (i.e. orig.ident).
group_by	Column in meta.data slot to group results by (i.e. "Treatment").
colors_use	List of colors or color palette to use. Only applicable if group_by is not NULL.
dot_size	size of the dots plotted if group_by is not NULL. Default is 1.
plot_title	Plot title.
y_axis_label	Label for y axis.
x_axis_label	Label for x axis.
legend_title	Label for plot legend.
x_lab_rotate	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.

Value

A ggplot object

Examples

```
## Not run:
library(Seurat)
cd_features <- list(c('CD79B', 'CD79A', 'CD19', 'CD180', 'CD200', 'CD3D', 'CD2', 'CD3E',
'CD7', 'CD8A', 'CD14', 'CD1C', 'CD68', 'CD9', 'CD247'))

pbmc_small <- AddModuleScore(object = pbmc_small, features = cd_features, ctrl = 5,
name = 'CD_Features')

Plot_Median_Other(seurat_object = pbmc_small, median_var = "CD_Features1",
sample_col = "orig.ident", group_by = "Treatment")

## End(Not run)
```

Plot_Median_UMIs

Plot Median UMIs per Cell per Sample

Description

Plot of median UMIs per cell per sample grouped by desired meta data variable.

Usage

```
Plot_Median_UMIs(
  seurat_object,
  sample_col = "orig.ident",
  group_by = NULL,
  colors_use = NULL,
  dot_size = 1,
  plot_title = "Median UMIs/Cell per Sample",
  y_axis_label = "Median UMIs",
  x_axis_label = NULL,
  legend_title = NULL,
  x_lab_rotate = TRUE,
  color_seed = 123
)
```

Arguments

seurat_object	Seurat object name.
sample_col	Specify which column in meta.data specifies sample ID (i.e. orig.ident).
group_by	Column in meta.data slot to group results by (i.e. "Treatment").
colors_use	List of colors or color palette to use. Only applicable if group_by is not NULL.
dot_size	size of the dots plotted if group_by is not NULL. Default is 1.
plot_title	Plot title.

y_axis_label	Label for y axis.
x_axis_label	Label for x axis.
legend_title	Label for plot legend.
x_lab_rotate	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.

Value

A ggplot object

Examples

```
library(Seurat)
# Create example groups
pbmc_small$sample_id <- sample(c("sample1", "sample2"), size = ncol(pbmc_small), replace = TRUE)

# Plot
Plot_Median_UMIs(seurat_object = pbmc_small, sample_col = "orig.ident", group_by = "sample_id")
```

Pull_Cluster_Annotation

Pull cluster information from annotation csv file.

Description

shortcut filter and pull function compatible with annotation files created by Create_Cluster_Annotation_File by default but also any other csv file.

Usage

```
Pull_Cluster_Annotation(
  annotation = NULL,
  cluster_name_col = "cluster",
  cell_type_col = "cell_type"
)
```

Arguments

annotation name of the data.frame/tibble or path to CSV file containing cluster annotation.

cluster_name_col name of column containing cluster names/numbers (default is "cluster").

cell_type_col name of column contain the cell type annotation (default is "cell_type").

Value

a list of named vectors for every cell type in the `cell_type_col` column of the annotation table and vectors new cluster names (for use with `Rename_Clusters` function or manual identity renaming).

Examples

```
## Not run:
# If pulling from a data.frame/tibble
cluster_annotation <- Pull_Cluster_Annotation(annotation = annotation_df,
cluster_name_col = "cluster", cell_type_col = "cell_type")

# If pulling from csv file
cluster_annotation <- Pull_Cluster_Annotation(annotation = "file_path/file_name.csv",
cluster_name_col = "cluster", cell_type_col = "cell_type")

## End(Not run)
```

Pull_Directory_List *Pull Directory List*

Description

Enables easy listing of all sub-directories for use as input library lists in Read10X multi functions.

Usage

```
Pull_Directory_List(base_path)
```

Arguments

`base_path` path to the parent directory which contains all of the subdirectories of interest.

Value

A vector of sub-directories within `base_path`.

Examples

```
## Not run:
data_dir <- 'path/to/data/directory'
library_list <- Pull_Directory_List(base_path = data_dir)

## End(Not run)
```

 QC_Histogram

QC Histogram Plots

Description

Custom histogram for initial QC checks including lines for thresholding

Usage

```
QC_Histogram(
  seurat_object,
  features,
  low_cutoff = NULL,
  high_cutoff = NULL,
  split.by = NULL,
  bins = 250,
  colors_use = "dodgerblue",
  num_columns = NULL,
  plot_title = NULL,
  assay = NULL,
  print_defaults = FALSE
)
```

Arguments

seurat_object	Seurat object name.
features	Feature from meta.data, assay features, or feature name shortcut to plot.
low_cutoff	Plot line a potential low threshold for filtering.
high_cutoff	Plot line a potential high threshold for filtering.
split.by	Feature to split plots by (i.e. "orig.ident").
bins	number of bins to plot default is 250.
colors_use	color to fill histogram bars, default is "dodgerblue".
num_columns	Number of columns in plot layout.
plot_title	optional, vector to use for plot title. Default is the name of the variable being plotted.
assay	assay to pull features from, default is active assay.
print_defaults	return list of accepted default shortcuts to provide to features instead of full name.

Value

A patchwork object

Examples

```
## Not run:
QC_Histogram(seurat_object = object, features = "nFeature_RNA")

## End(Not run)
```

QC_Plots_Combined_Vln *QC Plots Genes, UMIs, & % Mito*

Description

Custom VlnPlot for initial QC checks including lines for thresholding

Usage

```
QC_Plots_Combined_Vln(
  seurat_object,
  group.by = NULL,
  feature_cutoffs = NULL,
  UMI_cutoffs = NULL,
  mito_cutoffs = NULL,
  mito_name = "percent_mito",
  pt.size = NULL,
  plot_median = FALSE,
  median_size = 15,
  plot_boxplot = FALSE,
  colors_use = NULL,
  x_lab_rotate = TRUE,
  y_axis_log = FALSE,
  raster = NULL,
  ggplot_default_colors = FALSE,
  color_seed = 123,
  ...
)
```

Arguments

seurat_object	Seurat object name.
group.by	Name of one or more metadata columns to group (color) cells by (for example, orig.ident); default is the current active.ident of the object.
feature_cutoffs	Numeric vector of length 1 or 2 to plot lines for potential low/high threshold for filtering.
UMI_cutoffs	Numeric vector of length 1 or 2 to plot lines for potential low/high threshold for filtering.

<code>mito_cutoffs</code>	Numeric vector of length 1 or 2 to plot lines for potential low/high threshold for filtering.
<code>mito_name</code>	The column name containing percent mitochondrial counts information. Default value is "percent_mito" which is default value created when using <code>Add_Mito_Ribo_Seurat()</code> .
<code>pt.size</code>	Point size for plotting
<code>plot_median</code>	logical, whether to plot median for each ident on the plot (Default is FALSE).
<code>median_size</code>	Shape size for the median is plotted.
<code>plot_boxplot</code>	logical, whether to plot boxplot inside of violin (Default is FALSE).
<code>colors_use</code>	vector of colors to use for plot.
<code>x_lab_rotate</code>	Rotate x-axis labels 45 degrees (Default is TRUE).
<code>y_axis_log</code>	logical. Whether to change y axis to log10 scale (Default is FALSE).
<code>raster</code>	Convert points to raster format. Default is NULL which will rasterize by default if greater than 100,000 total points plotted (# Cells x # of features).
<code>ggplot_default_colors</code>	logical. If <code>colors_use = NULL</code> , Whether or not to return plot using default <code>ggplot2</code> "hue" palette instead of default "polychrome" or "varibow" palettes.
<code>color_seed</code>	random seed for the "varibow" palette shuffle if <code>colors_use = NULL</code> and number of groups plotted is greater than 36. Default = 123.
<code>...</code>	Extra parameters passed to <code>VlnPlot</code> .

Value

A ggplot object

Examples

```
## Not run:
QC_Plots_Combined_Vln(seurat_object = object)

## End(Not run)
```

`QC_Plots_Complexity` *QC Plots Cell "Complexity"*

Description

Custom `VlnPlot` for initial QC checks including lines for thresholding

Usage

```

QC_Plots_Complexity(
  seurat_object,
  feature = "log10GenesPerUMI",
  group.by = NULL,
  x_axis_label = NULL,
  y_axis_label = "log10(Genes) / log10(UMIs)",
  plot_title = "Cell Complexity",
  low_cutoff = NULL,
  high_cutoff = NULL,
  pt.size = NULL,
  plot_median = FALSE,
  plot_boxplot = FALSE,
  median_size = 15,
  colors_use = NULL,
  x_lab_rotate = TRUE,
  y_axis_log = FALSE,
  raster = NULL,
  ggplot_default_colors = FALSE,
  color_seed = 123,
  ...
)

```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>feature</code>	Feature from Meta Data to plot.
<code>group.by</code>	Name of one or more metadata columns to group (color) cells by (for example, <code>orig.ident</code>); default is the current <code>active.ident</code> of the object.
<code>x_axis_label</code>	Label for x axis.
<code>y_axis_label</code>	Label for y axis.
<code>plot_title</code>	Plot Title.
<code>low_cutoff</code>	Plot line a potential low threshold for filtering.
<code>high_cutoff</code>	Plot line a potential high threshold for filtering.
<code>pt.size</code>	Point size for plotting
<code>plot_median</code>	logical, whether to plot median for each ident on the plot (Default is FALSE).
<code>plot_boxplot</code>	logical, whether to plot boxplot inside of violin (Default is FALSE).
<code>median_size</code>	Shape size for the median is plotted.
<code>colors_use</code>	vector of colors to use for plot.
<code>x_lab_rotate</code>	Rotate x-axis labels 45 degrees (Default is TRUE).
<code>y_axis_log</code>	logical. Whether to change y axis to log10 scale (Default is FALSE).
<code>raster</code>	Convert points to raster format. Default is NULL which will rasterize by default if greater than 100,000 total points plotted (# Cells x # of features).

```

ggplot_default_colors logical. If colors_use = NULL, Whether or not to return plot using default gg-
                        plot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
color_seed             random seed for the "varibow" palette shuffle if colors_use = NULL and number
                        of groups plotted is greater than 36. Default = 123.
...                   Extra parameters passed to VlnPlot.

```

Value

A ggplot object

Examples

```

library(Seurat)
pbmc_small <- Add_Cell_Complexity(pbmc_small)

QC_Plots_Complexity(seurat_object = pbmc_small)

```

QC_Plots_Feature	<i>QC Plots Feature</i>
------------------	-------------------------

Description

Custom VlnPlot for initial QC checks including lines for thresholding

Usage

```

QC_Plots_Feature(
  seurat_object,
  feature,
  group.by = NULL,
  x_axis_label = NULL,
  y_axis_label = NULL,
  plot_title = NULL,
  low_cutoff = NULL,
  high_cutoff = NULL,
  pt.size = NULL,
  plot_median = FALSE,
  median_size = 15,
  plot_boxplot = FALSE,
  colors_use = NULL,
  x_lab_rotate = TRUE,
  y_axis_log = FALSE,
  raster = NULL,
  ggplot_default_colors = FALSE,
  color_seed = 123,
  ...
)

```

Arguments

seurat_object	Seurat object name.
feature	Feature from Meta Data to plot.
group.by	Name of one or more metadata columns to group (color) cells by (for example, orig.ident); default is the current active.ident of the object.
x_axis_label	Label for x axis.
y_axis_label	Label for y axis.
plot_title	Plot Title.
low_cutoff	Plot line a potential low threshold for filtering.
high_cutoff	Plot line a potential high threshold for filtering.
pt.size	Point size for plotting.
plot_median	logical, whether to plot median for each ident on the plot (Default is FALSE).
median_size	Shape size for the median is plotted.
plot_boxplot	logical, whether to plot boxplot inside of violin (Default is FALSE).
colors_use	vector of colors to use for plot.
x_lab_rotate	Rotate x-axis labels 45 degrees (Default is TRUE).
y_axis_log	logical. Whether to change y axis to log10 scale (Default is FALSE).
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 100,000 total points plotted (# Cells x # of features).
ggplot_default_colors	logical. If colors_use = NULL, Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.
...	Extra parameters passed to VlnPlot .

Value

A ggplot object

Examples

```
## Not run:
QC_Plots_Feature(seurat_object = object, feature = "FEATURE_NAME",
y_axis_label = "FEATURE per Cell", plot_title = "FEATURE per Cell", high_cutoff = 10,
low_cutoff = 2)

## End(Not run)
```

 QC_Plots_Genes

QC Plots Genes

Description

Custom VlnPlot for initial QC checks including lines for thresholding

Usage

```
QC_Plots_Genes(
  seurat_object,
  plot_title = "Genes Per Cell/Nucleus",
  group.by = NULL,
  x_axis_label = NULL,
  y_axis_label = "Features",
  low_cutoff = NULL,
  high_cutoff = NULL,
  pt.size = NULL,
  plot_median = FALSE,
  plot_boxplot = FALSE,
  median_size = 15,
  colors_use = NULL,
  x_lab_rotate = TRUE,
  y_axis_log = FALSE,
  raster = NULL,
  ggplot_default_colors = FALSE,
  color_seed = 123,
  ...
)
```

Arguments

seurat_object	Seurat object name.
plot_title	Plot Title.
group.by	Name of one or more metadata columns to group (color) cells by (for example, orig.ident); default is the current active.ident of the object.
x_axis_label	Label for x axis.
y_axis_label	Label for y axis.
low_cutoff	Plot line a potential low threshold for filtering.
high_cutoff	Plot line a potential high threshold for filtering.
pt.size	Point size for plotting.
plot_median	logical, whether to plot median for each ident on the plot (Default is FALSE).
plot_boxplot	logical, whether to plot boxplot inside of violin (Default is FALSE).
median_size	Shape size for the median is plotted.

colors_use	vector of colors to use for plot.
x_lab_rotate	Rotate x-axis labels 45 degrees (Default is TRUE).
y_axis_log	logical. Whether to change y axis to log10 scale (Default is FALSE).
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 100,000 total points plotted (# Cells x # of features).
ggplot_default_colors	logical. If colors_use = NULL, Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.
...	Extra parameters passed to VlnPlot .

Value

A ggplot object

Examples

```
library(Seurat)
QC_Plots_Genes(seurat_object = pbmc_small, plot_title = "Genes per Cell", low_cutoff = 40,
high_cutoff = 85)
```

QC_Plots_Mito

QC Plots Mito

Description

#' Custom VlnPlot for initial QC checks including lines for thresholding

Usage

```
QC_Plots_Mito(
  seurat_object,
  mito_name = "percent_mito",
  plot_title = "Mito Gene % per Cell/Nucleus",
  group.by = NULL,
  x_axis_label = NULL,
  y_axis_label = "% Mitochondrial Gene Counts",
  low_cutoff = NULL,
  high_cutoff = NULL,
  pt.size = NULL,
  plot_median = FALSE,
  median_size = 15,
  plot_boxplot = FALSE,
  colors_use = NULL,
```

```

    x_lab_rotate = TRUE,
    y_axis_log = FALSE,
    raster = NULL,
    ggplot_default_colors = FALSE,
    color_seed = 123,
    ...
)

```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>mito_name</code>	The column name containing percent mitochondrial counts information. Default value is "percent_mito" which is default value created when using <code>Add_Mito_Ribo_Seurat()</code> .
<code>plot_title</code>	Plot Title.
<code>group.by</code>	Name of one or more metadata columns to group (color) cells by (for example, <code>orig.ident</code>); default is the current <code>active.ident</code> of the object.
<code>x_axis_label</code>	Label for x axis.
<code>y_axis_label</code>	Label for y axis.
<code>low_cutoff</code>	Plot line a potential low threshold for filtering.
<code>high_cutoff</code>	Plot line a potential high threshold for filtering.
<code>pt.size</code>	Point size for plotting.
<code>plot_median</code>	logical, whether to plot median for each ident on the plot (Default is FALSE).
<code>median_size</code>	Shape size for the median is plotted.
<code>plot_boxplot</code>	logical, whether to plot boxplot inside of violin (Default is FALSE).
<code>colors_use</code>	vector of colors to use for plot.
<code>x_lab_rotate</code>	Rotate x-axis labels 45 degrees (Default is TRUE).
<code>y_axis_log</code>	logical. Whether to change y axis to log10 scale (Default is FALSE).
<code>raster</code>	Convert points to raster format. Default is NULL which will rasterize by default if greater than 100,000 total points plotted (# Cells x # of features).
<code>ggplot_default_colors</code>	logical. If <code>colors_use = NULL</code> , Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
<code>color_seed</code>	random seed for the "varibow" palette shuffle if <code>colors_use = NULL</code> and number of groups plotted is greater than 36. Default = 123.
<code>...</code>	Extra parameters passed to <code>VlnPlot</code> .

Value

A ggplot object

Examples

```
## Not run:
QC_Plots_Mito(seurat_object = object, plot_title = "Percent Mito per Cell", high_cutoff = 10)

## End(Not run)
```

QC_Plots_UMIs

QC Plots UMIs

Description

#' Custom VlnPlot for initial QC checks including lines for thresholding

Usage

```
QC_Plots_UMIs(
  seurat_object,
  plot_title = "UMIs per Cell/Nucleus",
  group.by = NULL,
  x_axis_label = NULL,
  y_axis_label = "UMIs",
  low_cutoff = NULL,
  high_cutoff = NULL,
  pt.size = NULL,
  plot_median = FALSE,
  median_size = 15,
  plot_boxplot = FALSE,
  colors_use = NULL,
  x_lab_rotate = TRUE,
  y_axis_log = FALSE,
  raster = NULL,
  ggplot_default_colors = FALSE,
  color_seed = 123,
  ...
)
```

Arguments

seurat_object	Seurat object name.
plot_title	Plot Title.
group.by	Name of one or more metadata columns to group (color) cells by (for example, orig.ident); default is the current active.ident of the object.
x_axis_label	Label for x axis.
y_axis_label	Label for y axis.
low_cutoff	Plot line a potential low threshold for filtering.

high_cutoff	Plot line a potential high threshold for filtering.
pt.size	Point size for plotting.
plot_median	logical, whether to plot median for each ident on the plot (Default is FALSE).
median_size	Shape size for the median is plotted.
plot_boxplot	logical, whether to plot boxplot inside of violin (Default is FALSE).
colors_use	vector of colors to use for plot.
x_lab_rotate	Rotate x-axis labels 45 degrees (Default is TRUE).
y_axis_log	logical. Whether to change y axis to log10 scale (Default is FALSE).
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 100,000 total points plotted (# Cells x # of features).
ggplot_default_colors	logical. If colors_use = NULL, Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.
...	Extra parameters passed to VlnPlot .

Value

A ggplot object

Examples

```
library(Seurat)
QC_Plots_UMIs(seurat_object = pbmc_small, plot_title = "UMIs per Cell", low_cutoff = 75,
high_cutoff = 600)
```

QC_Plot_GenevsFeature *QC Plots Genes vs Misc*

Description

Custom FeatureScatter for initial QC checks including lines for thresholding

Usage

```
QC_Plot_GenevsFeature(
  seurat_object,
  feature1,
  x_axis_label = NULL,
  y_axis_label = "Genes per Cell/Nucleus",
  low_cutoff_gene = NULL,
  high_cutoff_gene = NULL,
```

```

    low_cutoff_feature = NULL,
    high_cutoff_feature = NULL,
    colors_use = NULL,
    pt.size = 1,
    group.by = NULL,
    raster = NULL,
    raster.dpi = c(512, 512),
    ggplot_default_colors = FALSE,
    color_seed = 123,
    shuffle_seed = 1,
    ...
)

```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>feature1</code>	First feature to plot.
<code>x_axis_label</code>	Label for x axis.
<code>y_axis_label</code>	Label for y axis.
<code>low_cutoff_gene</code>	Plot line a potential low threshold for filtering genes per cell.
<code>high_cutoff_gene</code>	Plot line a potential high threshold for filtering genes per cell.
<code>low_cutoff_feature</code>	Plot line a potential low threshold for filtering feature1 per cell.
<code>high_cutoff_feature</code>	Plot line a potential high threshold for filtering feature1 per cell.
<code>colors_use</code>	vector of colors to use for plotting by identity.
<code>pt.size</code>	Adjust point size for plotting.
<code>group.by</code>	Name of one or more metadata columns to group (color) cells by (for example, <code>orig.ident</code>). Default is <code>@active.ident</code> .
<code>raster</code>	Convert points to raster format. Default is <code>NULL</code> which will rasterize by default if greater than 100,000 cells.
<code>raster.dpi</code>	Pixel resolution for rasterized plots, passed to <code>geom_scattermore()</code> . Default is <code>c(512, 512)</code> .
<code>ggplot_default_colors</code>	logical. If <code>colors_use = NULL</code> , Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
<code>color_seed</code>	random seed for the "varibow" palette shuffle if <code>colors_use = NULL</code> and number of groups plotted is greater than 36. Default = 123.
<code>shuffle_seed</code>	Sets the seed if randomly shuffling the order of points (Default is 1).
<code>...</code>	Extra parameters passed to FeatureScatter .

Value

A ggplot object

Examples

```
## Not run:
QC_Plot_GenevsFeature(seurat_object = obj, y_axis_label = "Feature per Cell")

## End(Not run)
```

QC_Plot_UMIvsFeature *QC Plots UMI vs Misc*

Description

Custom FeatureScatter for initial QC checks including lines for thresholding

Usage

```
QC_Plot_UMIvsFeature(
  seurat_object,
  feature1,
  x_axis_label = NULL,
  y_axis_label = "UMIs per Cell/Nucleus",
  low_cutoff_UMI = NULL,
  high_cutoff_UMI = NULL,
  low_cutoff_feature = NULL,
  high_cutoff_feature = NULL,
  colors_use = NULL,
  pt.size = 1,
  group.by = NULL,
  raster = NULL,
  raster.dpi = c(512, 512),
  ggplot_default_colors = FALSE,
  color_seed = 123,
  shuffle_seed = 1,
  ...
)
```

Arguments

seurat_object	Seurat object name.
feature1	First feature to plot.
x_axis_label	Label for x axis.
y_axis_label	Label for y axis.
low_cutoff_UMI	Plot line a potential low threshold for filtering UMI per cell.
high_cutoff_UMI	Plot line a potential high threshold for filtering UMI per cell.

low_cutoff_feature	Plot line a potential low threshold for filtering feature1 per cell.
high_cutoff_feature	Plot line a potential high threshold for filtering feature1 per cell.
colors_use	vector of colors to use for plotting by identity.
pt.size	Adjust point size for plotting.
group.by	Name of one or more metadata columns to group (color) cells by (for example, orig.ident). Default is @active.ident.
raster	Convert points to raster format. Default is NULL which will rasterize by default if greater than 100,000 cells.
raster.dpi	Pixel resolution for rasterized plots, passed to geom_scattermore(). Default is c(512, 512).
ggplot_default_colors	logical. If colors_use = NULL, Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.
shuffle_seed	Sets the seed if randomly shuffling the order of points (Default is 1).
...	Extra parameters passed to FeatureScatter .

Value

A ggplot object

Examples

```
## Not run:
QC_Plot_UMIvsFeature(seurat_object = obj, y_axis_label = "Feature per Cell")

## End(Not run)
```

QC_Plot_UMIvsGene *QC Plots Genes vs UMIs*

Description

Custom FeatureScatter for initial QC checks including lines for thresholding

Usage

```

QC_Plot_UMIvsGene(
  seurat_object,
  x_axis_label = "UMIs per Cell/Nucleus",
  y_axis_label = "Genes per Cell/Nucleus",
  low_cutoff_gene = -Inf,
  high_cutoff_gene = Inf,
  low_cutoff_UMI = -Inf,
  high_cutoff_UMI = Inf,
  colors_use = NULL,
  meta_gradient_name = NULL,
  meta_gradient_color = viridis_plasma_dark_high,
  meta_gradient_na_color = "lightgray",
  meta_gradient_low_cutoff = NULL,
  cells = NULL,
  combination = FALSE,
  ident_legend = TRUE,
  pt.size = 1,
  group.by = NULL,
  raster = NULL,
  raster.dpi = c(512, 512),
  ggplot_default_colors = FALSE,
  color_seed = 123,
  shuffle_seed = 1,
  ...
)

```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>x_axis_label</code>	Label for x axis.
<code>y_axis_label</code>	Label for y axis.
<code>low_cutoff_gene</code>	Plot line a potential low threshold for filtering genes per cell.
<code>high_cutoff_gene</code>	Plot line a potential high threshold for filtering genes per cell.
<code>low_cutoff_UMI</code>	Plot line a potential low threshold for filtering UMIs per cell.
<code>high_cutoff_UMI</code>	Plot line a potential high threshold for filtering UMIs per cell.
<code>colors_use</code>	vector of colors to use for plotting by identity.
<code>meta_gradient_name</code>	Name of continuous meta data variable to color points in plot by. (MUST be continuous variable i.e. "percent_mito").
<code>meta_gradient_color</code>	The gradient color palette to use for plotting of meta variable (default is viridis "Plasma" palette with dark colors high).

<code>meta_gradient_na_color</code>	Color to use for plotting values when a <code>meta_gradient_low_cutoff</code> is set (default is "lightgray").
<code>meta_gradient_low_cutoff</code>	Value to use as threshold for plotting. <code>meta_gradient_name</code> values below this value will be plotted using <code>meta_gradient_na_color</code> .
<code>cells</code>	Cells to include on the scatter plot (default is all cells).
<code>combination</code>	logical (default FALSE). Whether or not to return a plot layout with both the plot colored by identity and the meta data gradient plot.
<code>ident_legend</code>	logical, whether to plot the legend containing identities (left plot) when <code>combination = TRUE</code> . Default is TRUE.
<code>pt.size</code>	Passes size of points to both FeatureScatter and <code>geom_point</code> .
<code>group.by</code>	Name of one or more metadata columns to group (color) cells by (for example, <code>orig.ident</code>). Default is <code>@active.ident</code> .
<code>raster</code>	Convert points to raster format. Default is NULL which will rasterize by default if greater than 100,000 cells.
<code>raster.dpi</code>	Pixel resolution for rasterized plots, passed to <code>geom_scattermore()</code> . Default is <code>c(512, 512)</code> .
<code>ggplot_default_colors</code>	logical. If <code>colors_use = NULL</code> , Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
<code>color_seed</code>	Random seed for the "varibow" palette shuffle if <code>colors_use = NULL</code> and number of groups plotted is greater than 36. Default = 123.
<code>shuffle_seed</code>	Sets the seed if randomly shuffling the order of points (Default is 1).
<code>...</code>	Extra parameters passed to FeatureScatter .

Value

A ggplot object

Examples

```
library(Seurat)
QC_Plot_UMIvsGene(seurat_object = pbmc_small, x_axis_label = "UMIs per Cell/Nucleus",
y_axis_label = "Genes per Cell/Nucleus")
```

Read10X_GEO

Load in NCBI GEO data from 10X

Description

Enables easy loading of sparse data matrices provided by 10X genomics. That have file prefixes added to them by NCBI GEO or other repos.

Usage

```

Read10X_GEO(
  data_dir = NULL,
  sample_list = NULL,
  sample_names = NULL,
  gene.column = 2,
  cell.column = 1,
  unique.features = TRUE,
  strip.suffix = FALSE,
  parallel = FALSE,
  num_cores = NULL,
  merge = FALSE
)

```

Arguments

<code>data_dir</code>	Directory containing the matrix.mtx, genes.tsv (or features.tsv), and barcodes.tsv files provided by 10X.
<code>sample_list</code>	A vector of file prefixes/names if specific samples are desired. Default is NULL and will load all samples in given directory.
<code>sample_names</code>	a set of sample names to use for each sample entry in returned list. If NULL will set names to the file name of each sample.
<code>gene.column</code>	Specify which column of genes.tsv or features.tsv to use for gene names; default is 2.
<code>cell.column</code>	Specify which column of barcodes.tsv to use for cell names; default is 1.
<code>unique.features</code>	Make feature names unique (default TRUE).
<code>strip.suffix</code>	Remove trailing "-1" if present in all cell barcodes.
<code>parallel</code>	logical (default FALSE). Whether to use multiple cores when reading in data. Only possible on Linux based systems.
<code>num_cores</code>	if <code>parallel = TRUE</code> indicates the number of cores to use for multicore processing.
<code>merge</code>	logical (default FALSE) whether or not to merge samples into a single matrix or return list of matrices. If TRUE each sample entry in list will have cell barcode prefix added. The prefix will be taken from <code>sample_names</code> .

Value

If features.csv indicates the data has multiple data types, a list containing a sparse matrix of the data from each type will be returned. Otherwise a sparse matrix containing the expression data will be returned.

References

Code used in function has been slightly modified from `Seurat::Read10X` function of Seurat package <https://github.com/satijalab/seurat> (License: GPL-3). Function was modified to sup-

port file prefixes and altered loop by Samuel Marsh for scCustomize (also previously posted as potential PR to Seurat GitHub).

Examples

```
## Not run:
data_dir <- 'path/to/data/directory'
expression_matrices <- Read10X_GEO(data_dir = data_dir)
# To create object from single file
seurat_object = CreateSeuratObject(counts = expression_matrices[[1]])

## End(Not run)
```

Read10X_h5_GEO

Load in NCBI GEO data from 10X in HDF5 file format

Description

Enables easy loading of HDF5 data matrices provided by 10X genomics. That have file prefixes added to them by NCBI GEO or other repos or programs (i.e. Cell Bender)

Usage

```
Read10X_h5_GEO(
  data_dir = NULL,
  sample_list = NULL,
  sample_names = NULL,
  shared_suffix = NULL,
  parallel = FALSE,
  num_cores = NULL,
  merge = FALSE,
  ...
)
```

Arguments

<code>data_dir</code>	Directory containing the .h5 files provided by 10X.
<code>sample_list</code>	A vector of file prefixes/names if specific samples are desired. Default is NULL and will load all samples in given directory.
<code>sample_names</code>	a set of sample names to use for each sample entry in returned list. If NULL will set names to the file name of each sample.
<code>shared_suffix</code>	a suffix and file extension shared by all samples.
<code>parallel</code>	logical (default FALSE). Whether to use multiple cores when reading in data. Only possible on Linux based systems.
<code>num_cores</code>	if <code>parallel = TRUE</code> indicates the number of cores to use for multicore processing.

merge logical (default FALSE) whether or not to merge samples into a single matrix or return list of matrices. If TRUE each sample entry in list will have cell barcode prefix added. The prefix will be taken from `sample_names`.

... Additional arguments passed to [Read10X_h5](#)

Value

If the data has multiple data types, a list containing a sparse matrix of the data from each type will be returned. Otherwise a sparse matrix containing the expression data will be returned.

Examples

```
## Not run:
data_dir <- 'path/to/data/directory'
expression_matrices <- Read10X_h5_GEO(data_dir = data_dir)
# To create object from single file
seurat_object = CreateSeuratObject(counts = expression_matrices[[1]])

## End(Not run)
```

Read10X_h5_Multi_Directory

Load 10X h5 count matrices from multiple directories

Description

Enables easy loading of sparse data matrices provided by 10X genomics that are present in multiple subdirectories. Can function with either default output directory structure of Cell Ranger or custom directory structure.

Usage

```
Read10X_h5_Multi_Directory(
  base_path,
  secondary_path = NULL,
  default_10X_path = TRUE,
  cellranger_multi = FALSE,
  h5_filename = "filtered_feature_bc_matrix.h5",
  cell_bender = deprecated(),
  sample_list = NULL,
  sample_names = NULL,
  replace_suffix = FALSE,
  new_suffix_list = NULL,
  parallel = FALSE,
  num_cores = NULL,
  merge = FALSE,
  ...
)
```

Arguments

<code>base_path</code>	path to the parent directory which contains all of the subdirectories of interest.
<code>secondary_path</code>	path from the parent directory to count matrix files for each sample.
<code>default_10X_path</code>	logical (default TRUE) sets the secondary path variable to the default 10X directory structure.
<code>cellranger_multi</code>	logical, whether samples were processed with Cell Ranger multi, default is FALSE.
<code>h5_filename</code>	name of h5 file (including .h5 suffix). If all h5 files have same name (i.e. Cell Ranger output) then use full file name. By default function uses Cell Ranger name: "filtered_feature_bc_matrix.h5". If h5 files have sample specific prefixes (i.e. from Cell Bender) then use only the shared part of file name (e.g., "_filtered_out.h5").
<code>cell_bender</code>	[Deprecated] CellBender read functions are now independent family of functions. See <code>Read_CellBender_*</code> functions.
<code>sample_list</code>	a vector of sample directory names if only specific samples are desired. If NULL will read in subdirectories in parent directory.
<code>sample_names</code>	a set of sample names to use for each sample entry in returned list. If NULL will set names to the subdirectory name of each sample.
<code>replace_suffix</code>	logical (default FALSE). Whether or not to replace the barcode suffixes of matrices using Replace_Suffix .
<code>new_suffix_list</code>	a vector of new suffixes to replace existing suffixes if <code>replace_suffix = TRUE</code> . See Replace_Suffix for more information. To remove all suffixes set <code>new_suffix_list = ""</code> .
<code>parallel</code>	logical (default FALSE) whether or not to use multi core processing to read in matrices.
<code>num_cores</code>	how many cores to use for parallel processing.
<code>merge</code>	logical (default FALSE) whether or not to merge samples into a single matrix or return list of matrices. If TRUE each sample entry in list will have cell barcode prefix added. The prefix will be taken from <code>sample_names</code> .
<code>...</code>	Extra parameters passed to Read10X_h5 .

Value

a list of sparse matrices (merge = FALSE) or a single sparse matrix (merge = TRUE).

Examples

```
## Not run:
base_path <- 'path/to/data/directory'
expression_matrices <- Read10X_h5_Multi_Directory(base_path = base_path)

## End(Not run)
```

 Read10X_Multi_Directory

Load 10X count matrices from multiple directories

Description

Enables easy loading of sparse data matrices provided by 10X genomics that are present in multiple subdirectories. Can function with either default output directory structure of Cell Ranger or custom directory structure.

Usage

```
Read10X_Multi_Directory(
  base_path,
  secondary_path = NULL,
  default_10X_path = TRUE,
  cellranger_multi = FALSE,
  sample_list = NULL,
  sample_names = NULL,
  parallel = FALSE,
  num_cores = NULL,
  merge = FALSE,
  ...
)
```

Arguments

<code>base_path</code>	path to the parent directory which contains all of the subdirectories of interest.
<code>secondary_path</code>	path from the parent directory to count matrix files for each sample.
<code>default_10X_path</code>	logical (default TRUE) sets the secondary path variable to the default 10X directory structure.
<code>cellranger_multi</code>	logical, whether samples were processed with Cell Ranger multi, default is FALSE.
<code>sample_list</code>	a vector of sample directory names if only specific samples are desired. If NULL will read in subdirectories in parent directory.
<code>sample_names</code>	a set of sample names to use for each sample entry in returned list. If NULL will set names to the subdirectory name of each sample.
<code>parallel</code>	logical (default FALSE) whether or not to use multi core processing to read in matrices.
<code>num_cores</code>	how many cores to use for parallel processing.
<code>merge</code>	logical (default FALSE) whether or not to merge samples into a single matrix or return list of matrices. If TRUE each sample entry in list will have cell barcode prefix added. The prefix will be taken from <code>sample_names</code> .
<code>...</code>	Extra parameters passed to Read10X .

Value

a list of sparse matrices (merge = FALSE) or a single sparse matrix (merge = TRUE).

Examples

```
## Not run:
base_path <- 'path/to/data/directory'
expression_matrices <- Read10X_Multi_Directory(base_path = base_path)

## End(Not run)
```

Read_CellBender_h5_Mat

Load CellBender h5 matrices (corrected)

Description

Extract sparse matrix with corrected counts from CellBender h5 output file.

Usage

```
Read_CellBender_h5_Mat(
  file_name,
  use.names = TRUE,
  unique.features = TRUE,
  h5_group_name = NULL,
  feature_slot_name = "features"
)
```

Arguments

file_name	Path to h5 file.
use.names	Label row names with feature names rather than ID numbers (default TRUE).
unique.features	Make feature names unique (default TRUE).
h5_group_name	Name of the group within H5 file that contains count data. This is only required if H5 file contains multiple subgroups and non-default names. Default is NULL.
feature_slot_name	Name of the slot contain feature names/ids. Must be one of: "features"(Cell Ranger v3+) or "genes" (Cell Ranger v1/v2 or STARsolo). Default is "features".

Value

sparse matrix

References

Code used in function has been modified from `Seurat::Read10X_h5` function of Seurat package <https://github.com/satijalab/seurat> (License: GPL-3).

Examples

```
## Not run:
mat <- Read_CellBender_h5_Mat(file_name = "/SampleA_out_filtered.h5")

## End(Not run)
```

Read_CellBender_h5_Multi_Directory

Load CellBender h5 matrices (corrected) from multiple directories

Description

Extract sparse matrix with corrected counts from CellBender h5 output file across multiple sample subdirectories.

Usage

```
Read_CellBender_h5_Multi_Directory(
  base_path,
  secondary_path = NULL,
  filtered_h5 = TRUE,
  custom_name = NULL,
  sample_list = NULL,
  sample_names = NULL,
  h5_group_name = NULL,
  feature_slot_name = "features",
  replace_suffix = FALSE,
  new_suffix_list = NULL,
  parallel = FALSE,
  num_cores = NULL,
  merge = FALSE,
  ...
)
```

Arguments

`base_path` path to the parent directory which contains all of the subdirectories of interest.

`secondary_path` path from the parent directory to count matrix files for each sample.

`filtered_h5` logical (default TRUE). Will set the shared file name suffix `custom_name` is NULL.

custom_name	if file name was customized in CellBender then this parameter should contain the portion of file name that is shared across all samples. Must include the ".h5" extension as well.
sample_list	a vector of sample directory names if only specific samples are desired. If NULL will read in subdirectories in parent directory.
sample_names	a set of sample names to use for each sample entry in returned list. If NULL will set names to the subdirectory name of each sample.
h5_group_name	Name of the group within H5 file that contains count data. This is only required if H5 file contains multiple subgroups and non-default names. Default is NULL.
feature_slot_name	Name of the slot contain feature names/ids. Must be one of: "features"(Cell Ranger v3+) or "genes" (Cell Ranger v1/v2 or STARsolo). Default is "features".
replace_suffix	logical (default FALSE). Whether or not to replace the barcode suffixes of matrices using Replace_Suffix .
new_suffix_list	a vector of new suffixes to replace existing suffixes if <code>replace_suffix = TRUE</code> . See Replace_Suffix for more information. To remove all suffixes set <code>new_suffix_list = ""</code> .
parallel	logical (default FALSE) whether or not to use multi core processing to read in matrices.
num_cores	how many cores to use for parallel processing.
merge	logical (default FALSE) whether or not to merge samples into a single matrix or return list of matrices. If TRUE each sample entry in list will have cell barcode prefix added. The prefix will be taken from <code>sample_names</code> .
...	Extra parameters passed to Read_CellBender_h5_Mat .

Value

list of sparse matrices

Examples

```
## Not run:
base_path <- 'path/to/data/directory'
mat_list <- Read_CellBender_h5_Multi_Directory(base_path = base_path)

## End(Not run)
```

 Read_CellBender_h5_Multi_File

Load CellBender h5 matrices (corrected) from multiple files

Description

Extract sparse matrix with corrected counts from CellBender h5 output file across multiple samples within the same directory.

Usage

```
Read_CellBender_h5_Multi_File(
  data_dir = NULL,
  filtered_h5 = TRUE,
  custom_name = NULL,
  sample_list = NULL,
  sample_names = NULL,
  h5_group_name = NULL,
  feature_slot_name = "features",
  parallel = FALSE,
  num_cores = NULL,
  merge = FALSE,
  ...
)
```

Arguments

<code>data_dir</code>	Directory containing the .h5 files output by CellBender.
<code>filtered_h5</code>	logical (default TRUE). Will set the shared file name suffix if <code>custom_name</code> is NULL.
<code>custom_name</code>	if file name was customized in CellBender then this parameter should contain the portion of file name that is shared across all samples. Must include the ".h5" extension as well.
<code>sample_list</code>	a vector of sample names if only specific samples are desired. If NULL will read in all files within <code>data_dir</code> directory.
<code>sample_names</code>	a set of sample names to use for each sample entry in returned list. If NULL will set names to the subdirectory name of each sample.
<code>h5_group_name</code>	Name of the group within H5 file that contains count data. This is only required if H5 file contains multiple subgroups and non-default names. Default is NULL.
<code>feature_slot_name</code>	Name of the slot contain feature names/ids. Must be one of: "features"(Cell Ranger v3+) or "genes" (Cell Ranger v1/v2 or STARsolo). Default is "features".
<code>parallel</code>	logical (default FALSE) whether or not to use multi core processing to read in matrices
<code>num_cores</code>	how many cores to use for parallel processing.

merge logical (default FALSE) whether or not to merge samples into a single matrix or return list of matrices. If TRUE each sample entry in list will have cell barcode prefix added. The prefix will be taken from `sample_names`.

... Extra parameters passed to [Read_CellBender_h5_Mat](#).

Value

list of sparse matrices

Examples

```
## Not run:
base_path <- 'path/to/data/directory'
mat_list <- Read_CellBender_h5_Multi_File(data_dir = base_path)

## End(Not run)
```

Read_GEO_Delim	<i>Load in NCBI GEO data formatted as single file per sample</i>
----------------	------------------------------------------------------------------

Description

Can read delimited file types (i.e. csv, tsv, txt)

Usage

```
Read_GEO_Delim(
  data_dir,
  file_suffix,
  move_genes_rownames = TRUE,
  sample_list = NULL,
  full_names = FALSE,
  sample_names = NULL,
  barcode_suffix_period = FALSE,
  parallel = FALSE,
  num_cores = NULL,
  merge = FALSE
)
```

Arguments

`data_dir` Directory containing the files.

`file_suffix` The file suffix of the individual files. Must be the same across all files being imported. This is used to detect files to import and their GEO IDs.

<code>move_genes_rownames</code>	logical. Whether gene IDs are present in first column or in row names of delimited file. If TRUE will move the first column to row names before creating final matrix. Default is TRUE.
<code>sample_list</code>	a vector of samples within directory to read in (can be either with or without <code>file_suffix</code> see <code>full_names</code>). If NULL will read in all subdirectories.
<code>full_names</code>	logical (default FALSE). Whether or not the <code>sample_list</code> vector includes the file suffix. If FALSE the function will add suffix based on <code>file_suffix</code> parameter.
<code>sample_names</code>	a set of sample names to use for each sample entry in returned list. If NULL will set names to the directory name of each sample.
<code>barcode_suffix_period</code>	Is the barcode suffix a period and should it be changed to "-". Default (FALSE; barcodes will be left identical to their format in input files.). If TRUE "." in barcode suffix will be changed to "-".
<code>parallel</code>	logical (default FALSE). Whether to use multiple cores when reading in data. Only possible on Linux based systems.
<code>num_cores</code>	if <code>parallel = TRUE</code> indicates the number of cores to use for multicore processing.
<code>merge</code>	logical (default FALSE) whether or not to merge samples into a single matrix or return list of matrices. If TRUE each sample entry in list will have cell barcode prefix added. The prefix will be taken from <code>sample_names</code> .

Value

List of gene x cell matrices in list format named by sample name.

Examples

```
## Not run:
data_dir <- 'path/to/data/directory'
expression_matrices <- Read_GEO_Delim(data_dir = data_dir)

## End(Not run)
```

Read_Metrics_10X

Read Overall Statistics from 10X Cell Ranger Count

Description

Get data.frame with all metrics from the Cell Ranger count analysis (present in `web_summary.html`)

Usage

```
Read_Metrics_10X(  
  base_path,  
  secondary_path = NULL,  
  default_10X = TRUE,  
  cellranger_multi = FALSE,  
  lib_list = NULL,  
  lib_names = NULL  
)
```

Arguments

<code>base_path</code>	path to the parent directory which contains all of the subdirectories of interest.
<code>secondary_path</code>	path from the parent directory to count "outs/" folder which contains the "metrics_summary.csv" file.
<code>default_10X</code>	logical (default TRUE) sets the secondary path variable to the default 10X directory structure.
<code>cellranger_multi</code>	logical, whether or not metrics come from Cell Ranger count or from Cell Ranger multi. Default is FALSE.
<code>lib_list</code>	a list of sample names (matching directory names) to import. If NULL will read in all samples in parent directory.
<code>lib_names</code>	a set of sample names to use for each sample. If NULL will set names to the directory name of each sample.

Value

A data frame with sample metrics from cell ranger.

Examples

```
## Not run:  
metrics <- Read_Metrics_10X(base_path = "/path/to/directories", default_10X = TRUE)  
  
## End(Not run)
```

Reduction>Loading_Present

Check if reduction loadings are present

Description

Check if reduction loadings are present in object and return vector of found loading names. Return warning messages for genes not found.

Usage

```
Reduction>Loading_Present(
  seurat_object,
  reduction_names,
  print_msg = TRUE,
  omit_warn = TRUE,
  return_none = FALSE
)
```

Arguments

seurat_object	object name.
reduction_names	vector of genes to check.
print_msg	logical. Whether message should be printed if all features are found. Default is TRUE.
omit_warn	logical. Whether to print message about features that are not found in current object. Default is TRUE.
return_none	logical. Whether list of found vs. bad features should still be returned if no features are found. Default is FALSE.

Value

A list of length 3 containing 1) found features, 2) not found features.

Examples

```
## Not run:
reductions <- Reduction>Loading_Present(seurat_object = obj_name, reduction_name = "PC_1")
found_features <- features[[1]]

## End(Not run)
```

Rename_Clusters

Rename Cluster Seurat

Description

Wrapper function to rename active identities in Seurat Object with new ids.

Usage

```
Rename_Clusters(seurat_object, new_ids, meta_col_name = NULL, ...)
```

Arguments

`seurat_object` object name.
`new_idents` vector of new cluster names. Must be equal to the length of current `active.ident` in Seurat Object. Will accept named vector (with old idents as names) or will name the `new_idents` vector internally.
`meta_col_name` (Optional). Whether or not to create new named column in `Object@meta.data` to store the old identities.
 ... Extra parameters passed to [RenameIdents](#).

Value

Seurat Object with new identities placed in `active.ident` slot.

Examples

```

## Not run:
obj <- Rename_Clusters(seurat_object = obj_name, new_idents = new_idents_vec,
meta_col_name = "Round01_Res0.6_Idents")

## End(Not run)

```

Replace_Suffix	<i>Replace barcode suffixes</i>
----------------	---------------------------------

Description

Replace barcode suffixes in matrix, data.frame, or list of matrices/data.frames

Usage

```
Replace_Suffix(data, current_suffix, new_suffix)
```

Arguments

`data` Either matrix/data.frame or list of matrices/data.frames with the cell barcodes in the column names.
`current_suffix` a single value or vector of values representing current barcode suffix. If suffix is the same for all matrices/data.frames in list only single value is required.
`new_suffix` a single value or vector of values representing new barcode suffix to be added. If desired suffix is the same for all matrices/data.frames in list only single value is required. If no suffix is desired set `new_suffix = ""`.

Value

matrix or data.frame with new column names.

Examples

```
## Not run:
dge_matrix <- Replace_Suffix(data = dge_matrix, current_suffix = "-1", new_suffix = "-2")

## End(Not run)
```

scCustomize_Palette *Color Palette Selection for scCustomize*

Description

Function to return package default discrete palettes depending on number of groups plotted.

Usage

```
scCustomize_Palette(
  num_groups,
  ggplot_default_colors = FALSE,
  color_seed = 123
)
```

Arguments

num_groups	number of groups to be plotted. If ggplot_default_colors = FALSE then by default: <ul style="list-style-type: none"> • If number of levels plotted equal to 2 then colors will be NavyAndOrange(). • If number of levels plotted greater than 2 but less than or equal to 36 it will use "polychrome" from DiscretePalette_scCustomize(). • If greater than 36 will use "varibow" with shuffle = TRUE from DiscretePalette_scCustomize.
ggplot_default_colors	logical. Whether to use default ggplot hue palette or not.
color_seed	random seed to use for shuffling the "varibow" palette.

Value

vector of colors to use for plotting.

Examples

```
cols <- scCustomize_Palette(num_groups = 24, ggplot_default_colors = FALSE)
PalettePlot(pal= cols)
```

Seq_QC_Plot_Alignment_Combined

QC Plots Sequencing metrics (Alignment) (Layout)

Description

Plot a combined plot of the Alignment QC metrics from sequencing output.

Usage

```
Seq_QC_Plot_Alignment_Combined(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  patchwork_title = "Sequencing QC Plots: Read Alignment Metrics",
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>patchwork_title</code>	Title to use for the patchworked plot output.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Alignment_Combined(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Antisense *QC Plots Sequencing metrics (Alignment)*

Description

Plot the fraction of reads mapped Antisense to Gene

Usage

```
Seq_QC_Plot_Antisense(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Antisense(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Basic_Combined
QC Plots Sequencing metrics (Layout)

Description

Plot a combined plot of the basic QC metrics from sequencing output.

Usage

```
Seq_QC_Plot_Basic_Combined(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  patchwork_title = "Sequencing QC Plots: Basic Cell Metrics",
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>patchwork_title</code>	Title to use for the patchworked plot output.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Basic_Combined(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Exonic *QC Plots Sequencing metrics (Alignment)*

Description

Plot the fraction of reads confidently mapped to Exonic regions

Usage

```
Seq_QC_Plot_Exonic(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Exonic(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Genes *QC Plots Sequencing metrics*

Description

Plot the median genes per cell per sample

Usage

```
Seq_QC_Plot_Genes(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Genes(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Genome *QC Plots Sequencing metrics (Alignment)*

Description

Plot the fraction of reads confidently mapped to genome

Usage

```
Seq_QC_Plot_Genome(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Genome(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Intergenic

QC Plots Sequencing metrics (Alignment)

Description

Plot the fraction of reads confidently mapped to intergenic regions

Usage

```
Seq_QC_Plot_Intergenic(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Intergenic(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Intronic *QC Plots Sequencing metrics (Alignment)*

Description

Plot the fraction of reads confidently mapped to intronic regions

Usage

```
Seq_QC_Plot_Intronic(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Intronic(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Number_Cells

QC Plots Sequencing metrics

Description

Plot the number of cells per sample

Usage

```
Seq_QC_Plot_Number_Cells(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

metrics_dataframe	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
plot_by	Grouping factor for the plot. Default is to plot as single group with single point per sample.
colors_use	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
dot_size	size of the dots plotted if plot_by is not sample_id Default is 1.
x_lab_rotate	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
significance	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
...	Other variables to pass to ggpubr::stat_compare_means when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Number_Cells(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Reads_in_Cells

QC Plots Sequencing metrics

Description

Plot the fraction of reads in cells per sample

Usage

```
Seq_QC_Plot_Reads_in_Cells(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Reads_in_Cells(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Reads_per_Cell

QC Plots Sequencing metrics

Description

Plot the mean number of reads per cell

Usage

```
Seq_QC_Plot_Reads_per_Cell(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Reads_per_Cell(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Saturation

QC Plots Sequencing metrics

Description

Plot the sequencing saturation percentage per sample

Usage

```
Seq_QC_Plot_Saturation(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Saturation(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Total_Genes

QC Plots Sequencing metrics

Description

Plot the total genes detected per sample

Usage

```
Seq_QC_Plot_Total_Genes(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Total_Genes(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_Transcriptome

QC Plots Sequencing metrics (Alignment)

Description

Plot the fraction of reads confidently mapped to transcriptome

Usage

```
Seq_QC_Plot_Transcriptome(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

<code>metrics_dataframe</code>	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
<code>plot_by</code>	Grouping factor for the plot. Default is to plot as single group with single point per sample.
<code>colors_use</code>	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
<code>dot_size</code>	size of the dots plotted if <code>plot_by</code> is not <code>sample_id</code> Default is 1.
<code>x_lab_rotate</code>	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
<code>significance</code>	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
<code>...</code>	Other variables to pass to <code>ggpubr::stat_compare_means</code> when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:
Seq_QC_Plot_Transcriptome(metrics_dataframe = metrics)

## End(Not run)
```

Seq_QC_Plot_UMIs	<i>QC Plots Sequencing metrics</i>
------------------	------------------------------------

Description

Plot the median UMIs per cell per sample

Usage

```
Seq_QC_Plot_UMIs(
  metrics_dataframe,
  plot_by = "sample_id",
  colors_use = NULL,
  dot_size = 1,
  x_lab_rotate = FALSE,
  significance = FALSE,
  ...
)
```

Arguments

metrics_dataframe	data.frame contain Cell Ranger QC Metrics (see Read_Metrics_10X).
plot_by	Grouping factor for the plot. Default is to plot as single group with single point per sample.
colors_use	colors to use for plot if plotting by group. Defaults to RColorBrewer Dark2 palette if less than 8 groups and DiscretePalette_scCustomize(palette = "polychrome") if more than 8.
dot_size	size of the dots plotted if plot_by is not sample_id Default is 1.
x_lab_rotate	logical. Whether to rotate the axes labels on the x-axis. Default is FALSE.
significance	logical. Whether to calculate and plot p-value comparisons when plotting by grouping factor. Default is FALSE.
...	Other variables to pass to ggpubr::stat_compare_means when doing significance testing.

Value

A ggplot object

Examples

```
## Not run:  
Seq_QC_Plot_UMIs(metrics_dataframe = metrics)  
  
## End(Not run)
```

Setup_scRNAseq_Project

Setup project directory structure

Description

Create reproducible project directory organization when initiating a new analysis.

Usage

```
Setup_scRNAseq_Project(  
  custom_dir_file = NULL,  
  cluster_annotation_path = NULL,  
  cluster_annotation_file_name = "cluster_annotation.csv"  
)
```

Arguments

`custom_dir_file`
file to file containing desired directory structure. Default is NULL and will provide generic built-in directory structure.

`cluster_annotation_path`
path to place cluster annotation file using [Create_Cluster_Annotation_File](#).

`cluster_annotation_file_name`
name to use for annotation file if created (optional).

Value

no return value. Creates system folders.

Examples

```
## Not run:  
# If using built-in directory structure.  
Setup_scRNAseq_Project()  
  
## End(Not run)
```

Single_Color_Palette *Single Color Palettes for Plotting*

Description

Selects colors from modified versions of RColorBrewer single colors palettes

Usage

```
Single_Color_Palette(pal_color, num_colors = NULL, seed_use = 123)
```

Arguments

pal_color	color palette to select (Options are: 'reds', 'blues', 'greens', 'purples', 'oranges', 'grays').
num_colors	set number of colors (max = 7).
seed_use	set seed for reproducibility (default: 123).

Value

A vector of colors

References

See RColorBrewer for more info on palettes <https://CRAN.R-project.org/package=RColorBrewer>

Examples

```
pal <- Single_Color_Palette(pal_color = "reds", num_colors = 7)
PalettePlot(pal= pal)
```

Split_Layers *Split Seurat object into layers*

Description

Split Assay5 of Seurat object into layers by variable in meta.data

Usage

```
Split_Layers(seurat_object, assay = "RNA", split.by)
```

Arguments

seurat_object Seurat object name.
assay name(s) of assays to convert. Defaults to current active assay.
split.by Variable in meta.data to use for splitting layers.

Examples

```
## Not run:  
# Split object by "treatment"  
obj <- Split_Layers(object = obj, assay = "RNA", split.by = "treatment")  
  
## End(Not run)
```

Split_Vector

Split vector into list

Description

Splits vector into chunks of x sizes

Usage

```
Split_Vector(x, chunk_size = 100, verbose = FALSE)
```

Arguments

x vector to split
chunk_size size of chunks for vector to be split into, default is 100.
verbose logical, print details of vector and split, default is FALSE.

Value

list with vector of X length

References

Base code from stackoverflow post: <https://stackoverflow.com/a/3321659/15568251>

Examples

```
vector <- c("gene1", "gene2", "gene3", "gene4", "gene5", "gene6")  
  
vector_list <- Split_Vector(x = vector, chunk_size = 3)
```

Stacked_VlnPlot *Stacked Violin Plot*

Description

Code for creating stacked violin plot gene expression.

Usage

```
Stacked_VlnPlot(
  seurat_object,
  features,
  group.by = NULL,
  split.by = NULL,
  idents = NULL,
  x_lab_rotate = FALSE,
  plot_legend = FALSE,
  colors_use = NULL,
  color_seed = 123,
  ggplot_default_colors = FALSE,
  plot_spacing = 0.15,
  spacing_unit = "cm",
  vln_linewidth = NULL,
  pt.size = NULL,
  raster = NULL,
  add.noise = TRUE,
  ...
)
```

Arguments

seurat_object	Seurat object name.
features	Features to plot.
group.by	Group (color) cells in different ways (for example, orig.ident).
split.by	A variable to split the violin plots by,
idents	Which classes to include in the plot (default is all).
x_lab_rotate	logical or numeric. If logical whether to rotate x-axis labels 45 degrees (Default is FALSE). If numeric must be either 45 or 90. Setting 45 is equivalent to setting TRUE.
plot_legend	logical. Adds plot legend containing idents to the returned plot.
colors_use	specify color palette to used in VlnPlot . By default if number of levels plotted is less than or equal to 36 it will use "polychrome" and if greater than 36 will use "varibow" with shuffle = TRUE both from <code>DiscretePalette_scCustomize</code> .
color_seed	random seed for the "varibow" palette shuffle if colors_use = NULL and number of groups plotted is greater than 36. Default = 123.

<code>ggplot_default_colors</code>	logical. If <code>colors_use = NULL</code> , Whether or not to return plot using default <code>ggplot2</code> "hue" palette instead of default "polychrome" or "varibow" palettes.
<code>plot_spacing</code>	Numerical value specifying the vertical spacing between each plot in the stack. Default is 0.15 ("cm"). Spacing dependent on unit provided to <code>spacing_unit</code> .
<code>spacing_unit</code>	Unit to use in specifying vertical spacing between plots. Default is "cm".
<code>vln_linewidth</code>	Adjust the linewidth of violin outline. Must be numeric.
<code>pt.size</code>	Adjust point size for plotting. Default for <code>Stacked_VlnPlot</code> is 0 to avoid issues with rendering so many points in vector form. Alternatively, see <code>raster</code> parameter.
<code>raster</code>	Convert points to raster format. Default is <code>NULL</code> which will rasterize by default if greater than 100,000 total points plotted (# Cells x # of features).
<code>add.noise</code>	logical, determine if adding a small noise for plotting (Default is <code>TRUE</code>).
<code>...</code>	Extra parameters passed to <code>VlnPlot</code> .

Value

A `ggplot` object

Author(s)

Ming Tang (Original Code), Sam Marsh (Wrap single function, added/modified functionality)

References

<https://divingintogeneticsandgenomics.rbind.io/post/stacked-violin-plot-for-visualizing-single-cell/>

See Also

<https://twitter.com/tangming2005>

Examples

```
library(Seurat)
Stacked_VlnPlot(seurat_object = pbmc_small, features = c("CD3E", "CD8", "GZMB", "MS4A1"),
x_lab_rotate = TRUE)
```

`Store_Misc_Info_Seurat`*Store misc data in Seurat object*

Description

Wrapper function save variety of data types to the `object@misc` slot of Seurat object.

Usage

```
Store_Misc_Info_Seurat(  
  seurat_object,  
  data_to_store,  
  data_name,  
  list_as_list = FALSE,  
  overwrite = FALSE  
)
```

Arguments

<code>seurat_object</code>	object name.
<code>data_to_store</code>	data to be stored in <code>@misc</code> slot. Can be single piece of data or list. If list of data see <code>list_as_list</code> parameter for control over data storage.
<code>data_name</code>	name to give the entry in <code>@misc</code> slot. Must be of equal length of the number of data items being stored.
<code>list_as_list</code>	logical. If <code>data_to_store</code> is a list, this dictates whether to store in <code>@misc</code> slot as list (TRUE) or whether to store each entry in the list separately (FALSE). Default is FALSE.
<code>overwrite</code>	Logical. Whether to overwrite existing items with the same name. Default is FALSE, meaning that function will abort if item with <code>data_name</code> is present in <code>misc</code> slot.

Value

Seurat Object with new entries in the `@misc` slot.

Examples

```
library(Seurat)  
clu_pal <- c("red", "green", "blue")  
  
pbmc_small <- Store_Misc_Info_Seurat(seurat_object = pbmc_small, data_to_store = clu_pal,  
  data_name = "rd1_colors")
```

Store_Palette_Seurat *Store color palette in Seurat object*

Description

Wrapper function around Store_Misc_Info_Seurat to store color palettes.

Usage

```
Store_Palette_Seurat(  
  seurat_object,  
  palette,  
  palette_name,  
  list_as_list = FALSE,  
  overwrite = FALSE  
)
```

Arguments

seurat_object	object name.
palette	vector or list of vectors containing color palettes to store. If list of palettes see list_as_list parameter for control over data storage.
palette_name	name to give the palette(s) in @misc slot. Must be of equal length to the number of data items being stored.
list_as_list	logical. If data_to_store is a list, this dictates whether to store in @misc slot as list (TRUE) or whether to store each entry in the list separately (FALSE). Default is FALSE.
overwrite	Logical. Whether to overwrite existing items with the same name. Default is FALSE, meaning that function will abort if item with data_name is present in misc slot.

Value

Seurat Object with new entries in the @misc slot.

Examples

```
library(Seurat)  
clu_pal <- c("red", "green", "blue")  
  
pbmc_small <- Store_Misc_Info_Seurat(seurat_object = pbmc_small, data_to_store = clu_pal,  
data_name = "rd1_colors")
```

theme_ggprism_mod	<i>Modified ggprism theme</i>
-------------------	-------------------------------

Description

Modified ggprism theme which restores the legend title.

Usage

```
theme_ggprism_mod(
  palette = "black_and_white",
  base_size = 14,
  base_family = "sans",
  base_fontface = "bold",
  base_line_size = base_size/20,
  base_rect_size = base_size/20,
  axis_text_angle = 0,
  border = FALSE
)
```

Arguments

palette	string. Palette name, use <code>names(ggprism_data\$themes)</code> to show all valid palette names.
base_size	numeric. Base font size, given in "pt".
base_family	string. Base font family, default is "sans".
base_fontface	string. Base font face, default is "bold".
base_line_size	numeric. Base linewidth for line elements
base_rect_size	numeric. Base linewidth for rect elements
axis_text_angle	integer. Angle of axis text in degrees. One of: 0, 45, 90, 270.
border	logical. Should a border be drawn around the plot? Clipping will occur unless e.g. <code>coord_cartesian(clip = "off")</code> is used.

Value

Returns a list-like object of class *theme*.

References

theme is a modified version of `theme_prism` from ggprism package <https://github.com/csdaw/ggprism> (License: GPL-3). Param text is from ggprism: `theme_prism()` documentation [theme_prism](#). Theme adaptation based on ggprism vignette <https://csdaw.github.io/ggprism/articles/themes.html#make-your-own-ggprism-theme-1>.

Examples

```
# Generate a plot and customize theme
library(ggplot2)
df <- data.frame(x = rnorm(n = 100, mean = 20, sd = 2), y = rbinom(n = 100, size = 100, prob = 0.2))
p <- ggplot(data = df, mapping = aes(x = x, y = y)) + geom_point(mapping = aes(color = 'red'))
p + theme_ggprism_mod()
```

Top_Genes_Factor	<i>Extract top loading genes for LIGER factor</i>
------------------	---------------------------------------------------

Description

Extract vector to the top loading genes for specified LIGER iNMF factor

Usage

```
Top_Genes_Factor(liger_object, liger_factor, num_genes = 10)
```

Arguments

liger_object	LIGER object name.
liger_factor	LIGER factor number to pull genes from.
num_genes	number of top loading genes to return as vector.

Value

A LIGER Object

Examples

```
## Not run:
top_genes_factor10 <- Top_Genes_Factor(liger_object = object, num_genes = 10)

## End(Not run)
```

`UnRotate_X`*Unrotate x axis on VlnPlot*

Description

Shortcut for thematic modification to unrotate the x axis (e.g., for Seurat VlnPlot is rotated by default).

Usage

```
UnRotate_X(...)
```

Arguments

```
...          extra arguments passed to ggplot2::theme().
```

Value

Returns a list-like object of class *theme*.

Examples

```
library(Seurat)
p <- VlnPlot(object = pbmc_small, features = "CD3E")
p + UnRotate_X()
```

`Updated_HGNC_Symbols`*Update HGNC Gene Symbols*

Description

Update human gene symbols using data from HGNC. This function will store cached data in package directory using (BiocFileCache). Use of this function requires internet connection on first use (or if setting `update_symbol_data = TRUE`). Subsequent use does not require connection and will pull from cached data.

Usage

```
Updated_HGNC_Symbols(  
  input_data,  
  update_symbol_data = NULL,  
  case_check_as_warn = FALSE,  
  verbose = TRUE  
)
```

Arguments

<code>input_data</code>	Data source containing gene names. Accepted formats are: <ul style="list-style-type: none"> • character vector • Seurat Objects • <code>data.frame</code>: genes as rownames • <code>dgCMatrx/dgTMatrix</code>: genes as rownames • <code>tibble</code>: genes in first column
<code>update_symbol_data</code>	logical, whether to update cached HGNC data, default is <code>NULL</code> . If <code>NULL</code> <code>BiocFileCache</code> will check and prompt for update if cache is stale. If <code>FALSE</code> the <code>BiocFileCache</code> stale check will be skipped and current cache will be used. If <code>TRUE</code> the <code>BiocFileCache</code> stale check will be skipped and HGNC data will be downloaded.
<code>case_check_as_warn</code>	logical, whether case checking of features should cause abort or only warn, default is <code>FALSE</code> (abort). Set to <code>TRUE</code> if atypical names (i.e. old LOC naming) are present in <code>input_data</code> .
<code>verbose</code>	logical, whether to print results detailing numbers of symbols, found, updated, and not found; default is <code>TRUE</code> .

Value

`data.frame` containing columns: `input_features`, `Approved_Symbol` (already approved; output unchanged), `Not_Found_Symbol` (symbol not in HGNC; output unchanged), `Updated_Symbol` (new symbol from HGNC; output updated).

Examples

```
## Not run:
new_names <- Updated_HGNC_Symbols(input_data = Seurat_Object)

## End(Not run)
```

VariableFeaturePlot_scCustom
Custom Labeled Variable Features Plot

Description

Creates variable features plot with N number of features already labeled by default.

Usage

```
VariableFeaturePlot_scCustom(
  seurat_object,
  num_features = 10,
  custom_features = NULL,
  label = TRUE,
  pt.size = 1,
  colors_use = c("black", "red"),
  repel = TRUE,
  y_axis_log = FALSE,
  assay = NULL,
  selection.method = NULL,
  ...
)
```

Arguments

<code>seurat_object</code>	Seurat object name.
<code>num_features</code>	Number of top variable features to highlight by color/label.
<code>custom_features</code>	A vector of custom feature names to label on plot instead of labeling top variable genes.
<code>label</code>	logical. Whether to label the top features. Default is TRUE.
<code>pt.size</code>	Adjust point size for plotting.
<code>colors_use</code>	colors to use for plotting. Default is "black" and "red".
<code>repel</code>	logical (default TRUE). Whether or not to repel the feature labels on plot.
<code>y_axis_log</code>	logical. Whether to change y axis to log10 scale (Default is FALSE).
<code>assay</code>	Assay to pull variable features from.
<code>selection.method</code>	If more than one method use to calculate variable features specify which method to use for plotting. See <code>selection.method</code> parameter in VariableFeaturePlot for list of options.
<code>...</code>	Extra parameters passed to VariableFeaturePlot .

Value

A ggplot object

Examples

```
library(Seurat)
VariableFeaturePlot_scCustom(seurat_object = pbmc_small, num_features = 10)
```

 Variable_Features_ALL_LIGER

Perform variable gene selection over whole dataset

Description

Performs variable gene selection for LIGER object across the entire object instead of by dataset and then taking union.

Usage

```
Variable_Features_ALL_LIGER(
  liger_object,
  num_genes = NULL,
  var.thresh = 0.3,
  alpha.thresh = 0.99,
  tol = 1e-04,
  do.plot = FALSE,
  pt.size = 0.3,
  chunk = 1000
)
```

Arguments

<code>liger_object</code>	LIGER object name.
<code>num_genes</code>	Number of genes to find. Optimizes the value of <code>var.thresh</code> to get this number of genes, (Default is NULL).
<code>var.thresh</code>	Variance threshold. Main threshold used to identify variable genes. Genes with expression variance greater than threshold (relative to mean) are selected. (higher threshold -> fewer selected genes).
<code>alpha.thresh</code>	Alpha threshold. Controls upper bound for expected mean gene expression (lower threshold -> higher upper bound). (default 0.99)
<code>tol</code>	Tolerance to use for optimization if <code>num_genes</code> values passed in (default 0.0001).
<code>do.plot</code>	Display log plot of gene variance vs. gene expression. Selected genes are plotted in green. (Default FALSE)
<code>pt.size</code>	Point size for plot.
<code>chunk</code>	size of chunks in hdf5 file. (Default 1000)

Value

A LIGER Object with variable genes in correct slot.

References

Matching function parameter text descriptions are taken from `rliger::selectGenes` which is called by this function after creating new temporary object/dataset. <https://github.com/welch-lab/liger>. (License: GPL-3).

Value

A color palette for plotting

Examples

```
## Not run:
FeaturePlot_scCustom(object = seurat_object, features = "Cx3cr1",
  colors_use = viridis_plasma_dark_high, na_color = "lightgray")

## End(Not run)
```

VlnPlot_scCustom

VlnPlot with modified default settings

Description

Creates DimPlot with some of the settings modified from their Seurat defaults (colors_use, shuffle, label).

Usage

```
VlnPlot_scCustom(
  seurat_object,
  features,
  colors_use = NULL,
  pt.size = NULL,
  group.by = NULL,
  split.by = NULL,
  plot_median = FALSE,
  plot_boxplot = FALSE,
  median_size = 15,
  idents = NULL,
  num_columns = NULL,
  raster = NULL,
  add.noise = TRUE,
  ggplot_default_colors = FALSE,
  color_seed = 123,
  ...
)
```

Arguments

seurat_object Seurat object name.
 features Feature(s) to plot.

<code>colors_use</code>	color palette to use for plotting. By default if number of levels plotted is less than or equal to 36 it will use "polychrome" and if greater than 36 will use "varibow" with <code>shuffle = TRUE</code> both from <code>DiscretePalette_scCustomize</code> .
<code>pt.size</code>	Adjust point size for plotting.
<code>group.by</code>	Name of one or more metadata columns to group (color) cells by (for example, <code>orig.ident</code>); default is the current <code>active.ident</code> of the object.
<code>split.by</code>	Feature to split plots by (i.e. "orig.ident").
<code>plot_median</code>	logical, whether to plot median for each ident on the plot (Default is FALSE).
<code>plot_boxplot</code>	logical, whether to plot boxplot inside of violin (Default is FALSE).
<code>median_size</code>	Shape size for the median is plotted.
<code>idents</code>	Which classes to include in the plot (default is all).
<code>num_columns</code>	Number of columns in plot layout. Only valid if <code>split.by</code> != NULL.
<code>raster</code>	Convert points to raster format. Default is NULL which will rasterize by default if greater than 100,000 total points plotted (# Cells x # of features).
<code>add.noise</code>	logical, determine if adding a small noise for plotting (Default is TRUE).
<code>ggplot_default_colors</code>	logical. If <code>colors_use = NULL</code> , Whether or not to return plot using default ggplot2 "hue" palette instead of default "polychrome" or "varibow" palettes.
<code>color_seed</code>	random seed for the "varibow" palette shuffle if <code>colors_use = NULL</code> and number of groups plotted is greater than 36. Default = 123.
<code>...</code>	Extra parameters passed to <code>VlnPlot</code> .

Value

A ggplot object

References

Many of the param names and descriptions are from Seurat to facilitate ease of use as this is simply a wrapper to alter some of the default parameters <https://github.com/satijalab/seurat/blob/master/R/visualization.R> (License: GPL-3).

Examples

```
library(Seurat)
VlnPlot_scCustom(seurat_object = pbmc_small, features = "CD3E")
```

Index

- * **barcode_util**
 - Change_Delim_All, 28
 - Change_Delim_Prefix, 29
 - Change_Delim_Suffix, 30
 - Replace_Suffix, 135
- * **check_util**
 - Case_Check, 23
 - CheckMatrix_scCustom, 31
 - Feature_Present, 60
 - Gene_Present, 62
 - Meta_Numeric, 86
 - Meta_Present, 87
 - Reduction_Loading_Present, 133
- * **datasets**
 - ensembl_mito_id, 51
 - ensembl_ribo_id, 51
 - ieg_gene_list, 64
 - msigdb_qc_gene_list, 89
 - viridis_plasma_dark_high, 165
- * **data**
 - ensembl_mito_id, 51
 - ensembl_ribo_id, 51
 - ieg_gene_list, 64
 - msigdb_qc_gene_list, 89
- * **get_set_util**
 - Add_Alt_Feature_ID, 5
 - Add_Sample_Meta, 13
 - Extract_Sample_Meta, 52
 - Fetch_Meta, 61
 - Meta_Remove_Seurat, 88
 - Store_Misc_Info_Seurat, 157
 - Store_Palette_Seurat, 158
- * **iterative_plotting**
 - Iterate_Cluster_Highlight_Plot, 66
 - Iterate_DimPlot_bySample, 67
 - Iterate_FeaturePlot_scCustom, 68
 - Iterate_Meta_Highlight_Plot, 70
 - Iterate_PC_Loading_Plots, 72
 - Iterate_Plot_Density_Custom, 73
 - Iterate_Plot_Density_Joint, 74
 - Iterate_VlnPlot_scCustom, 76
- * **liger_object_util**
 - Add_Cell_Complexity, 7
 - Add_Mito_Ribo, 10
 - Fetch_Meta, 61
 - LIGER_Features, 78
 - Top_Genes_Factor, 160
 - Variable_Features_ALL_LIGER, 164
- * **liger_plotting**
 - DimPlot_LIGER, 44
 - plotFactors_scCustom, 94
- * **marker_annotation_util**
 - Add_Pct_Diff, 12
 - Create_Cluster_Annotation_File, 42
 - Extract_Top_Markers, 53
 - Pull_Cluster_Annotation, 104
 - Rename_Clusters, 134
- * **misc_util**
 - Merge_Seurat_List, 82
 - Split_Vector, 154
 - Updated_HGNC_Symbols, 161
- * **object_conversion**
 - as.anndata, 16
 - as.LIGER, 18
 - as.Seurat.liger, 20
 - Convert_Assay, 38
 - Liger_to_Seurat, 79
 - Split_Layers, 153
- * **object_qc_plotting**
 - QC_Histogram, 106
 - QC_Plot_GenevsFeature, 116
 - QC_Plot_UMIvsFeature, 118
 - QC_Plot_UMIvsGene, 119
 - QC_Plots_Combined_Vln, 107
 - QC_Plots_Complexity, 108
 - QC_Plots_Feature, 110
 - QC_Plots_Genes, 112
 - QC_Plots_Mito, 113

- QC_Plots_UMIs, 115
- * **organization_util**
 - Copy_From_GCP, 38
 - Copy_To_GCP, 39
 - Setup_scRNAseq_Project, 152
- * **other_seurat_plotting**
 - Plot_Density_Custom, 97
 - Plot_Density_Joint_Only, 98
- * **palettes**
 - ColorBlind_Pal, 37
 - Dark2_Pal, 42
 - DiscretePalette_scCustomize, 48
 - Hue_Pal, 63
 - JCO_Four, 77
 - NavyAndOrange, 91
 - PalettePlot, 91
 - scCustomize_Palette, 136
 - Single_Color_Palette, 153
 - viridis_plasma_dark_high, 165
- * **qc_util**
 - Add_Cell_Complexity, 7
 - Add_Cell_QC_Metrics, 8
 - Add_CellBender_Diff, 6
 - Add_Mito_Ribo, 10
 - Add_Top_Gene_Pct_Seurat, 15
- * **read_&_write**
 - Create_10X_H5, 40
 - Create_CellBender_Merged_Seurat, 41
 - Pull_Directory_List, 105
 - Read10X_GEO, 121
 - Read10X_h5_GEO, 123
 - Read10X_h5_Multi_Directory, 124
 - Read10X_Multi_Directory, 126
 - Read_CellBender_h5_Mat, 127
 - Read_CellBender_h5_Multi_Directory, 128
 - Read_CellBender_h5_Multi_File, 130
 - Read_GEO_Delim, 131
 - Read_Metrics_10X, 132
- * **read_merge_util**
 - Extract_Modality, 52
 - Merge_Sparse_Data_All, 83
 - Merge_Sparse_Multimodal_All, 84
- * **seq_qc_plotting_alignment**
 - Seq_QC_Plot_Antisense, 138
 - Seq_QC_Plot_Exonic, 140
 - Seq_QC_Plot_Genome, 142
 - Seq_QC_Plot_Intergenic, 143
 - Seq_QC_Plot_Intronic, 144
 - Seq_QC_Plot_Transcriptome, 150
- * **seq_qc_plotting_basic**
 - Barcode_Plot, 22
 - Iterate_Barcode_Rank_Plot, 64
 - Seq_QC_Plot_Genes, 141
 - Seq_QC_Plot_Number_Cells, 145
 - Seq_QC_Plot_Reads_in_Cells, 146
 - Seq_QC_Plot_Reads_per_Cell, 147
 - Seq_QC_Plot_Saturation, 148
 - Seq_QC_Plot_Total_Genes, 149
 - Seq_QC_Plot_UMIs, 151
- * **seq_qc_plotting_layout**
 - Seq_QC_Plot_Alignment_Combined, 137
 - Seq_QC_Plot_Basic_Combined, 139
- * **seurat_plotting**
 - Cell_Highlight_Plot, 27
 - Cluster_Highlight_Plot, 35
 - Clustered_DotPlot, 32
 - DimPlot_All_Samples, 43
 - DimPlot_scCustom, 46
 - DotPlot_scCustom, 49
 - FeaturePlot_DualAssay, 55
 - FeaturePlot_scCustom, 56
 - FeatureScatter_scCustom, 58
 - Meta_Highlight_Plot, 85
 - PC_Plotting, 92
 - Stacked_VlnPlot, 155
 - VariableFeaturePlot_scCustom, 162
 - VlnPlot_scCustom, 166
- * **stats_plotting**
 - CellBender_Diff_Plot, 24
 - Plot_Cells_per_Sample, 96
 - Plot_Median_Genes, 99
 - Plot_Median_Mito, 100
 - Plot_Median_Other, 102
 - Plot_Median_UMIs, 103
- * **stats**
 - CellBender_Feature_Diff, 26
 - Cluster_Stats_All_Samples, 36
 - MAD_Stats, 80
 - Median_Stats, 81
 - Percent_Expressing, 93
- * **themes**
 - Blank_Theme, 23
 - Move_Legend, 89

- theme_ggprism_mod, 159
- UnRotate_X, 161
- Add_Alt_Feature_ID, 5
- Add_Cell_Complexity, 7
- Add_Cell_QC_Metrics, 8
- Add_CellBender_Diff, 6
- Add_Mito_Ribo, 10
- Add_Pct_Diff, 12
- Add_Sample_Meta, 13
- Add_Top_Gene_Pct_Seurat, 15
- as.anndata, 16
- as.LIGER, 18
- as.Seurat.liger, 20
- Barcode_Plot, 22
- barcodeRanks, 22
- Blank_Theme, 23
- Case_Check, 23
- Cell_Highlight_Plot, 27
- CellBender_Diff_Plot, 24
- CellBender_Feature_Diff, 25, 26
- CellCycleScoring, 9
- Change_Delim_All, 28
- Change_Delim_Prefix, 29
- Change_Delim_Suffix, 30
- CheckMatrix_scCustom, 31
- Cluster_Highlight_Plot, 35
- Cluster_Stats_All_Samples, 36
- Clustered_DotPlot, 32
- ColorBlind_Pal, 37
- Convert_Assay, 21, 38, 79
- Copy_From_GCP, 38
- Copy_To_GCP, 39
- Create_10X_H5, 40
- Create_CellBender_Merged_Seurat, 41, 55
- Create_Cluster_Annotation_File, 42, 152
- CreateSeuratObject, 41
- Dark2_Pal, 42
- DimPlot, 28, 36, 44, 48, 66, 68, 71, 86
- DimPlot_All_Samples, 43
- DimPlot_LIGER, 44
- DimPlot_scCustom, 46
- DiscretePalette_scCustomize, 45, 48, 94
- DotPlot, 50, 93
- DotPlot_scCustom, 49
- ensembl_mito_id, 51
- ensembl_ribo_id, 51
- Extract_Modality, 52
- Extract_Sample_Meta, 52
- Extract_Top_Markers, 53
- Feature_Present, 60
- FeaturePlot, 56, 58, 70
- FeaturePlot_DualAssay, 55
- FeaturePlot_scCustom, 56
- FeatureScatter, 59, 117, 119, 121
- FeatureScatter_scCustom, 58
- Fetch_Meta, 61
- FindAllMarkers, 12, 13, 53, 54
- FindMarkers, 12, 13
- Gene_Present, 62
- geom_text_repel, 25, 26
- Hue_Pal, 63
- ieg_gene_list, 64
- Iterate_Barcode_Rank_Plot, 64
- Iterate_Cluster_Highlight_Plot, 66
- Iterate_DimPlot_bySample, 67
- Iterate_FeaturePlot_scCustom, 68
- Iterate_Meta_Highlight_Plot, 70
- Iterate_PC_Loading_Plots, 72
- Iterate_Plot_Density_Custom, 73
- Iterate_Plot_Density_Joint, 74
- Iterate_VlnPlot_scCustom, 76
- JCO_Four, 77
- LabelPoints, 26
- LIGER_Features, 78
- Liger_to_Seurat, 79
- MAD_Stats, 80
- Median_Stats, 53, 81
- merge, 82
- Merge_Seurat_List, 82
- Merge_Sparse_Data_All, 83
- Merge_Sparse_Multimodal_All, 84
- Meta_Highlight_Plot, 85
- Meta_Numeric, 86
- Meta_Present, 87
- Meta_Remove_Seurat, 88
- Move_Legend, 89
- msigdb_qc_gene_list, 89

- NavyAndOrange, [91](#)
- PalettePlot, [91](#)
- patchwork, [58](#)
- PC_Plotting, [92](#)
- PCHeatmap, [72, 92](#)
- Percent_Expressing, [93](#)
- Plot_Cells_per_Sample, [96](#)
- plot_density, [74, 75, 98, 99](#)
- Plot_Density_Custom, [97](#)
- Plot_Density_Joint_Only, [98](#)
- Plot_Median_Genes, [99](#)
- Plot_Median_Mito, [100](#)
- Plot_Median_Other, [102](#)
- Plot_Median_UMIs, [103](#)
- plotFactors_scCustom, [94](#)
- Pull_Cluster_Annotation, [104](#)
- Pull_Directory_List, [105](#)
- QC_Histogram, [106](#)
- QC_Plot_GenevsFeature, [116](#)
- QC_Plot_UMIvsFeature, [118](#)
- QC_Plot_UMIvsGene, [119](#)
- QC_Plots_Combined_Vln, [107](#)
- QC_Plots_Complexity, [108](#)
- QC_Plots_Feature, [110](#)
- QC_Plots_Genes, [112](#)
- QC_Plots_Mito, [113](#)
- QC_Plots_UMIs, [115](#)
- Read10X, [126](#)
- Read10X_GEO, [121](#)
- Read10X_h5, [124, 125](#)
- Read10X_h5_GEO, [123](#)
- Read10X_h5_Multi_Directory, [124](#)
- Read10X_Multi_Directory, [126](#)
- Read_CellBender_h5_Mat, [127, 129, 131](#)
- Read_CellBender_h5_Multi_Directory, [128](#)
- Read_CellBender_h5_Multi_File, [130](#)
- Read_GEO_Delim, [131](#)
- Read_Metrics_10X, [132, 137–151](#)
- Reduction_Loading_Present, [133](#)
- Rename_Clusters, [134](#)
- RenameIdsents, [135](#)
- Replace_Suffix, [125, 129, 135](#)
- scCustomize_Palette, [136](#)
- Seq_QC_Plot_Alignment_Combined, [137](#)
- Seq_QC_Plot_Antisense, [138](#)
- Seq_QC_Plot_Basic_Combined, [139](#)
- Seq_QC_Plot_Exonic, [140](#)
- Seq_QC_Plot_Genes, [141](#)
- Seq_QC_Plot_Genome, [142](#)
- Seq_QC_Plot_Intergenic, [143](#)
- Seq_QC_Plot_Intronic, [144](#)
- Seq_QC_Plot_Number_Cells, [145](#)
- Seq_QC_Plot_Reads_in_Cells, [146](#)
- Seq_QC_Plot_Reads_per_Cell, [147](#)
- Seq_QC_Plot_Saturation, [148](#)
- Seq_QC_Plot_Total_Genes, [149](#)
- Seq_QC_Plot_Transcriptome, [150](#)
- Seq_QC_Plot_UMIs, [151](#)
- Setup_scRNAseq_Project, [152](#)
- Single_Color_Palette, [153](#)
- Split_Layers, [153](#)
- Split_Vector, [154](#)
- Stacked_VlnPlot, [155](#)
- Store_Misc_Info_Seurat, [157](#)
- Store_Palette_Seurat, [158](#)
- theme_ggprism_mod, [159](#)
- theme_prism, [159](#)
- Top_Genes_Factor, [160](#)
- UnRotate_X, [161](#)
- Updated_HGNC_Symbols, [161](#)
- Variable_Features_ALL_LIGER, [164](#)
- VariableFeaturePlot, [163](#)
- VariableFeaturePlot_scCustom, [162](#)
- viridis_dark_high
 - (viridis_plasma_dark_high), [165](#)
- viridis_inferno_dark_high
 - (viridis_plasma_dark_high), [165](#)
- viridis_inferno_light_high
 - (viridis_plasma_dark_high), [165](#)
- viridis_light_high
 - (viridis_plasma_dark_high), [165](#)
- viridis_magma_dark_high
 - (viridis_plasma_dark_high), [165](#)
- viridis_magma_light_high
 - (viridis_plasma_dark_high), [165](#)
- viridis_plasma_dark_high, [165](#)
- viridis_plasma_light_high
 - (viridis_plasma_dark_high), [165](#)
- VizDimLoadings, [72, 92](#)

VInPlot, [77](#), [108](#), [110](#), [111](#), [113](#), [114](#), [116](#),
[155](#), [156](#), [167](#)
VInPlot_scCustom, [166](#)